

Gelest

Silane Coupling Agents

Connecting Across Boundaries

Metal Primers

Bind Biomaterials

Provide Crosslinking

Immobilize Catalysts

*Improve Polymer and
Particle Dispersion*

Enhance Adhesive Bonding

Increase Electrical Properties

Maximize Composite Strength

Increase Mechanical Properties


Version 3.0:

Water-borne Silanes

New Coupling Agents:

Cyclic Aza-Silanes, Azido-Silanes,
Dipodal Silanes

Oligomeric Hydrolysates



Gelest, Inc.

Telephone: General 215-547-1015
Order Entry 888-734-8344
Technical Service: 215-547-1016
FAX: 215-547-2484
Internet: www.gelest.com
e-mail: sales@gelest.com
Correspondence: 11 East Steel Rd.
Morrisville, PA 19067

For further information consult our web site at: www.gelest.com

In Mainland China:

For commercial and research quantities contact:

**A Meryer Chemical Technology
Shanghai Company**

No. 3636, Jiangcheng Road
Shanghai, China 200245
Tel: +86-(0)-21-61259170
Fax: +86-(0)-21-61259169
Email: pur02@meryer.com

In Japan:

For commercial and research quantities contact:

AZmax Co. Ltd. Tokyo Office

Matsuda Yaesudori, Bldg F8
1-10-7 Hatchoubori, Chou-Ku
Tokyo 104-0032
Tel: 81-3-5543-1630
Fax: 81-3-5543-0312
Email: sales@azmax.co.jp
On-line catalog: www.azmax.co.jp

In India:

For commercial and research quantities contact:

Gautavik International

301, A Wing Chandan Co-op Hsg Soc.
Dadabhai Cross Road North
Vile Parle West, Mumbai 400056
India
Tel: 91-22-26703175
Fax: 91-96-19190510
Email: vasantiyadav@yahoo.co.in

In South-East Asia:

For commercial and research quantities contact:

Gulf Chemical

39 Jalan Pemimpin
Tai Lee Industrial Building #04-03
Singapore 577182
Tel: 65-6358-3185
Fax: 65-6353-2542
Email: support@gulfchem.com.sg

In Taiwan:

For commercial and bulk quantities contact:

Kelly Chemical Corporation

9F, No.155, Sec.1, Keelung Rd, Taipei
Taiwan
Tel : +886-2-27621985
Fax: +886-2-27532400
Website: <http://www.kellychemical.com>
Email: chem@kellychemical.com



Gelest  **Silane Coupling Agents**
Connecting Across Boundaries

TABLE OF CONTENTS

What is a Silane Coupling Agent	2
How Does a Silane Coupling Agent Work	3
Selecting a Silane Coupling Agent - Inorganic Substrate Perspective.....	4
Selecting a Silane Coupling Agent - Interphase Considerations	5
Partition, Orientation and Self-Assembly in Bonded Phases.....	6
Selecting a Silane Coupling Agent - Polymer Applications	7
Special Topics:	
Linker Length	11
Dipodal Silanes.....	12
Cyclic Azasilanes	14
Thermal Stability of Silanes.....	16
Aqueous Systems & Water-Borne Silanes.....	17
Masked Silanes - Latent Functionality	18
Coupling Agents for Metal Substrates.....	19
Difficult Substrates.....	20
Applying a Silane Coupling Agent.....	21
Silane Coupling Agents for Polymers - Selection Chart.....	24
Silane Coupling Agents for Biomaterials - Selection Chart	27
Silane Coupling Agents - Properties	28
UV Active and Fluorescent Silanes	65
Chiral Silanes and Biomolecular Probes.....	67
Silyl Hydrides and Trihydridosilanes.....	69
Dipodal Silanes - Non-Functional.....	70
Organosilane Modified Silica Nanoparticles	71
Further Information - Other Resources	72




Silane Coupling Agents

Silane coupling agents have the ability to form a durable bond between organic and inorganic materials. Encounters between dissimilar materials often involve at least one member that's siliceous or has surface chemistry with siliceous properties; silicates, aluminates, borates, etc., are the principal components of the earth's crust. Interfaces involving such materials have become a dynamic area of chemistry in which surfaces have been modified in order to generate desired heterogeneous environments or to incorporate the bulk properties of different phases into a uniform composite structure.



The general formula for a silane coupling agent typically shows the two classes of functionality. X is a hydrolyzable group typically alkoxy, acyloxy, halogen or amine. Following hydrolysis, a reactive silanol group is formed, which can condense with other silanol groups, for example, those on the surface of siliceous fillers, to form siloxane linkages. Stable condensation products are also formed with other oxides such as those of aluminum, zirconium, tin, titanium, and nickel. Less stable bonds are formed with oxides of boron, iron, and carbon. Alkali metal oxides and carbonates do not form stable bonds with Si-O-. The R group is a nonhydrolyzable organic radical that may possess a functionality that imparts desired characteristics.

The final result of reacting an organosilane with a substrate ranges from altering the wetting or adhesion characteristics of the substrate, utilizing the substrate to catalyze chemical transformations at the heterogeneous interface, ordering the interfacial region, and modifying its partition characteristics. Significantly, it includes the ability to effect a covalent bond between organic and inorganic materials.



How Does a Silane Modify a Surface?

Most of the widely used organosilanes have one organic substituent and three hydrolyzable substituents. In the vast majority of surface treatment applications, the alkoxy groups of the trialkoxysilanes are hydrolyzed to form silanol-containing species. Reaction of these silanes involves four steps. Initially, hydrolysis of the three labile groups occurs. Condensation to oligomers follows. The oligomers then hydrogen bond with OH groups of the substrate. Finally, during drying or curing, a covalent linkage is formed with the substrate with concomitant loss of water. Although described sequentially, these reactions can occur simultaneously after the initial hydrolysis step. At the interface, there is usually only one bond from each silicon of the organosilane to the substrate surface. The two remaining silanol groups are present either in condensed or free form. The R group remains available for covalent reaction or physical interaction with other phases.

Silanes can modify surfaces under anhydrous conditions consistent with monolayer and vapor phase deposition requirements. Extended reaction times (4-12 hours) at elevated temperatures (50°-120°C) are typical. Of the alkoxysilanes, only methoxysilanes are effective without catalysis for vapor deposition. The most effective silanes for vapor phase deposition are cyclic azasilanes.


Hydrolysis Considerations

Water for hydrolysis may come from several sources. It may be added, it may be present on the substrate surface, or it may come from the atmosphere. The degree of polymerization of the silane is determined by the amount of water available and the organic substituent. If the silane is added to water and has low solubility, a high degree of polymerization is favored. Multiple organic substitution, particularly if phenyl or tertiary butyl groups are involved, favors formation of stable monomeric silanols.

The thickness of a polysiloxane layer is also determined by the concentration of the siloxane solution. Although a monolayer is generally desired, multilayer adsorption results from solutions customarily used. It has been calculated that deposition from a 0.25% silane solution onto glass could result in three to eight molecular layers. These multilayers could be either inter-connected through a loose network structure, or intermixed, or both, and are, in fact, formed by most deposition techniques. The orientation of functional groups is generally horizontal, but not necessarily planar, on the surface of the substrate.


The formation of covalent bonds to the surface proceeds with a certain amount of reversibility. As water is removed, generally by heating to 120°C for 30 to 90 minutes or evacuation for 2 to 6 hours, bonds may form, break, and reform to relieve internal stress. The same mechanism can permit a positional displacement of interface components.

Hydrolytic Deposition of Silanes



B. Arkles, CHEMTECH, 7, 766, 1977

Anhydrous Deposition of Silanes



Selecting a Silane for Surface Modification Inorganic Substrate Perspective

Factors influencing silane surface modification selection include:

Concentration of surface hydroxyl groups

Type of surface hydroxyl groups

Hydrolytic Stability of the bond formed



Physical dimensions of the substrate or substrate features

Surface modification is maximized when silanes react with the substrate surface and present the maximum number of accessible sites with appropriate surface energies. An additional consideration is the physical and chemical properties of the interphase region. The interphase can promote or detract from total system properties depending on its physical properties such as modulus or chemical properties such as water/hydroxyl content.

Hydroxyl-containing substrates vary widely in concentration and type of hydroxyl groups present. Freshly fused substrates stored under neutral conditions have a minimum number of hydroxyls. Hydrolytically derived oxides aged in moist air have significant amounts of physically adsorbed water which can interfere with coupling. Hydrogen bonded vicinal silanols react more readily with silane coupling agents, while isolated or free hydroxyls react reluctantly.

Silanes with three alkoxy groups are the usual starting point for substrate modification. These materials tend to deposit as polymeric films, effecting total coverage and maximizing the introduction of organic functionality. They are the primary materials utilized in composites, adhesives, sealants, and coatings. Limitations intrinsic in the utilization of a polylayer deposition are significant for nano-particles or nano-composites where the interphase dimensions generated by polylayer deposition may approach those of the substrate. Residual (non-condensed) hydroxyl groups from alkoxy silanes can also interfere in activity. Monoalkoxy-silanes provide a frequently used alternative for nano-featured substrates since deposition is limited to a monolayer.

If the hydrolytic stability of the oxane bond between the silane and the substrate is poor or the application is in an aggressive aqueous environment, dipodal silanes often exhibit substantial performance improvements. These materials form tighter networks and may offer up to 10⁵x greater hydrolysis resistance making them particularly appropriate for primer applications.



Water droplets on a (heptadecafluoro-1,1,2,2-tetrahydrodecyl)trimethoxysilane-treated silicon wafer exhibit high contact angles, indicative of the low surface energy. Surfaces are both hydrophobic and resist wetting by hydrocarbon oils. (water droplets contain dye for photographic purposes).

Silane Effectiveness on Inorganics

	SUBSTRATES
EXCELLENT	Silica
	Quartz
	Glass
	Aluminum (AlO(OH))
	Alumino-silicates (e.g. clays)
	Silicon
	Copper
	Tin (SnO)
	Talc
	Inorganic Oxides (e.g. Fe ₂ O ₃ , TiO ₂ , Cr ₂ O ₃)
GOOD	Steel, Iron
	Asbestos
	Nickel
	Zinc
	Lead
SLIGHT	Marble, Chalk (CaCO ₃)
	Gypsum (CaSO ₄)
	Barytes (BaSO ₄)
	Graphite
	Carbon Black
POOR	

Estimates for Silane Loading on Siliceous Fillers	
Average Particle Size	Amount of Silane (minimum of monolayer coverage)
<1 micron	1.5% or greater
1-10 microns	1.0%
10-20 microns	0.75%
>100 microns	0.1% or less

Critical Surface Tension and Adhesion

While the contact angle of water on a substrate is a good indicator of the relative hydrophobicity or hydrophilicity of a substrate, it is not a good indicator for the wettability of the substrate by other liquids. The contact angle is given by Young's equation:

$$\gamma_{sv} - \gamma_{sl} = \gamma_{lv} \cdot \cos\theta_c$$

where γ_{sl} = interfacial surface tension, γ_{lv} = surface tension of liquid.

Critical surface tension is associated with the wettability or release properties of a solid. It serves as a better predictor of the behavior of a solid with a range of liquids.

Liquids with a surface tension below the critical surface tension (γ_c) of a substrate will wet the surface, i.e., show a contact angle of 0 ($\cos\theta_c = 1$). The critical surface tension is unique for any solid and is determined by plotting the cosine of the contact angles of liquids of different surface tensions and extrapolating to 1.

Hydrophilic behavior is generally observed by surfaces with critical surface tensions greater than 45 dynes/cm. As the critical surface tension increases, the expected decrease in contact angle is accompanied with stronger adsorptive behavior and with increased exotherms.

Hydrophobic behavior is generally observed by surfaces with critical surface tensions less than 35 dynes/cm. At first, the decrease in critical surface tension is associated with oleophilic behavior, i.e. the wetting of the surfaces by hydrocarbon oils. As the critical surface tensions decrease below 20 dynes/cm, the surfaces resist wetting by hydrocarbon oils and are considered oleophobic as well as hydrophobic.

In the reinforcement of thermosets and thermoplastics with glass fibers, one approach for optimizing reinforcement is to match the critical surface tension of the silylated glass surface to the surface tension of the polymer in its melt or uncured condition. This has been most helpful in resins with no obvious functionality such as polyethylene and polystyrene. Silane treatment has allowed control of thixotropic activity of silica and clays in paint and coating applications. Immobilization of cellular organelles, including mitochondria, chloroplasts, and microsomes, has been effected by treating silica with alkylsilanes of C8 or greater substitution.

Critical surface tensions

	γ_c
	mN/m
heneicosafuorododecyltrichlorosilane	6-7
heptadecafluorodecyltrichlorosilane	12.0
poly(tetrafluoroethylene)	18.5
octadecyltrichlorosilane	20-24
methyltrimethoxysilane	22.5
nonafluorohexyltrimethoxysilane	23.0
vinyltriethoxysilane	25
paraffin wax	25.5
ethyltrimethoxysilane	27.0
propyltrimethoxysilane	28.5
glass, soda-lime (wet)	30.0
poly(chlorotrifluoroethylene)	31.0
poly(propylene)	31.0
poly(propylene oxide)	32
polyethylene	33.0
trifluoropropyltrimethoxysilane	33.5
3-(2-aminoethyl)-aminopropyltrimethoxysilane	33.5
poly(styrene)	34
p-tolyltrimethoxysilane	34
cyanoethyltrimethoxysilane	34
aminopropyltriethoxysilane	35
acetoxypolypropyltrimethoxysilane	37.5
polymethylmethacrylate	39
polyvinylchloride	39
phenyltrimethoxysilane	40.0
chloropropyltrimethoxysilane	40.5
mercaptopropyltrimethoxysilane	41
glycidoxypropyltrimethoxysilane	42.5
poly(ethyleneterephthalate)	43
poly(ethylene oxide)	43-45
copper (dry)	44
aluminum (dry)	45
iron (dry)	46
nylon 6/6	45-6
glass, soda-lime (dry)	47
silica, fused	78
titanium dioxide (anatase)	91
ferric oxide	107
tin oxide	111

Note: Critical surface tensions for silanes refer to smooth treated surfaces.


Selecting a Silane Coupling Agent: Polymer Applications

Coupling agents find their largest application in the area of polymers. Since any silane that enhances the adhesion of a polymer is often termed a coupling agent, regardless of whether or not a covalent bond is formed, the definition becomes vague. In this discussion, the parochial outlook will be adopted, and only silanes that form covalent bonds directly to the polymer will be considered. The covalent bond may be formed by reaction with the finished polymer or copolymerized with the monomer. Thermoplastic bonding is achieved through both routes, although principally the former. Thermosets are almost entirely limited to the latter. The mechanism and performance of silane coupling agents is best discussed with reference to specific systems. The most important substrate is E-type fiberglass, which has 6-15 silanol groups per μm^2 .


Thermosets

Acrylates, methacrylates and Unsaturated Polyesters, owing to their facility for undergoing free-radical polymerization, can be modified by copolymerization with silanes that have unsaturated organic substitution. The usual coupling agents for thermoset polyesters undergo radical copolymerization in such systems. These resins, usually of loosely defined structure, often have had their viscosity reduced by addition of a second monomer, typically styrene. In general, better reinforcement is obtained when the silane monomer matches the reactivity of the styrene rather than the maleate portion of the polyester.


Methacrylyl and styryl functional silanes undergo addition much more readily than vinylsilanes. A direct approach to selecting the optimum silane uses the e and Q parameters of the Alfrey-Price treatment of polymerization. Here e indicates the polarity of the monomer radical that forms at the end of a growing chain, while Q represents the resonance stabilization of a radical by adjacent groups. Optimum random copolymerization is obtained from monomers with similar orders of reactivity. Vinyl functional silanes mismatch the reactionary parameters of most unsaturated polyesters. However, they can be used in direct high pressure polymerization with olefins such as ethylene, propylene and dienes.



Acrylate Coupling Reaction



Unsaturated Polyester (Styrene) Coupling Reaction



Polyethylene Graft Coupling Reaction

Urethanes

Thermoset urethane can be effectively coupled with two types of silanes. The first type, including isocyanate functional silanes, may be used to treat the filler directly or integrally blended with the diisocyanate (TDI, MDI, etc.) prior to cure. Amine and alkanolamine functional silanes, on the other hand, are blended with the polyol rather than the diisocyanate. Isocyanate functional silanes couple with the polyol. Alkanolamine functional silanes react with the isocyanate to form urethane linkages, while amine silanes react with the isocyanates to yield urea linkages. A typical application for coupled urethane system is improving bond strength with sand in abrasion-resistant, sand-filled flooring resins.

Moisture-Cureable Urethanes


Secondary aminosilanes have the general ability to convert isocyanate functional urethane prepolymers to systems that crosslink in the presence of water and a tin catalyst. The preferred aminosilanes are secondary containing methyl, ethyl or butyl substitutions on nitrogen.

Epoxyes


Epoxy cyclohexyl and glycidoxy functional silanes are used to pretreat the filler or to blend with the glycidylbisphenol-A ether. Amine functional silanes can likewise be used to pretreat the filler or to blend with the hardener portion. Treatment of fillers in epoxy adhesives improves their dispersibility and increases the mechanical properties of the cured resin. A large application area is glass cloth-reinforced epoxy laminates and prepregs in aerospace and electrical printed circuit board applications.

Phenolics


Phenolic resins are divided into base catalyzed single-step resins called resols or better known acid catalyzed two-step systems called novolaks. Although foundry and molds are formulated with resols such as aminopropyl-methyldialkoxysilanes, the commercial utilization of silanes in phenolic resins is largely limited to novolak/glass fabric laminates and molding compounds. The phenolic hydroxyl group of the resins readily react with the oxirane ring of epoxy silanes to form phenyl ether linkages. When phenolic resins are compounded with rubbers, as in the case with nitrile/phenolic or vinyl butyral/phenolic adhesives, or impact-resistant molding compounds, additional silanes, particularly mercapto-functional silanes, have been found to impart greater bond strength than silanes that couple to the phenolic portion.




Polyurethane Coupling Reactions



Moisture-Cure Silicone Polyurethane (SPUR)



Epoxy Coupling Reaction



Phenolic Coupling Reaction

Thermoplastics

Thermoplastics provide a greater challenge in promoting adhesion through silane coupling agents than thermosets. The silanes must react with the polymer and not the monomeric precursors, which not only limits avenues for coupling, but also presents additional problems in rheology and thermal properties during composite formulation. Moreover mechanical requirements here are stringently determined. Polymers that contain regular sites for covalent reactivity either in the backbone or in a pendant group include polydienes, polyvinylchloride, polyphenylene sulfide, acrylic homopolymers, maleic anhydride, acrylic, vinyl acetate, diene-containing copolymers, and halogen or chlorosulfonyl-modified homopolymers. A surprising number of these are coupled by aminoalkylsilanes. Chlorinated polymers readily form quaternary compounds while the carboxylate and sulfonate groups form amides and sulfonamides under process conditions. At elevated temperatures, the amines add across many double bonds although mercaptoalkylsilanes are the preferred coupling agents.


The most widely used coupling agents, the aminoalkylsilanes are the most economical, but are not necessarily the best. Epoxysilanes, for example, are successfully used with acrylic acid and maleic acid copolymers.

Thermoplastic Condensation Polymers

The group of polymers that most closely approaches theoretical limits of composite strength does not appear to contain regular opportunities for covalent bond formation to substrate. Most of the condensation polymers including polyamides, polyesters, polycarbonates, and polysulfones are in this group. Adhesion is promoted by introducing high energy groups and hydrogen bond potential in the interphase area or by taking advantage of the relatively low molecular weight of these polymers, which results in a significant opportunity for end-group reactions. Aminoalkylsilanes, chloroalkylsilanes, and isocyanatosilanes are the usual candidates for coupling to these resins. This group has the greatest mechanical strength of the thermoplastics, allowing them to replace the cast metals in such typical uses as gears, connectors and bobbins.




Scanning electron micrograph at a broken gear tooth from a non-coupled glass fiber/acetal composite. Note that cleavage occurred between fibers.



Scanning electron micrograph at a broken gear tooth from an aminosilane-coupled glass fiber/nylon 6/6 composite. Note how fibers have broken as well as matrix.



Chopped fiberglass strand sized with aminosilanes is a commonly used reinforcement for high temperature thermoplastics.



Thermoplastic Polyester Coupling Reaction

Polyolefins

The polyolefins and polyethers present no direct opportunity for covalent coupling. Until recently, the principal approach for composite formulation was to match the surface energy of the filler surface, by treating it with an alkylsubstituted silane, with that of the polymer. For optimum reinforcement, preferred resins should be of high molecular weight, linear, and have low melt viscosity. Approaches to improved composite strength have been through compatibility with long-chain alkylsilanes or aminosilanes. Far more effective is coupling with vinyl or methacryloxy groups, particularly if additional coupling sites are created in the resin by addition of peroxides. Dicumyl peroxide and bis(t-butylperoxy) compounds at levels of 0.15% to 0.25% have been introduced into polyethylene compounded with vinylsilane-treated glass fibers for structural composites or vinylsilane-treated clay for wire insulation. Increases of 50% in tensile and flexural properties have been observed in both cases when compared to the same silane systems without peroxides.


Another approach for coupling polypropylene and polyethylene is through silylsulfonylazides. Unlike azide bound to silicon, sulfonyl azides decompose above 150°C to form a molecule of nitrogen and a reactive nitrene that is capable of insertion into carbon-hydrogen bonds, forming sulfonamides, into carbon-carbon double bonds, forming triazoles, and into aromatic bonds, forming sulfonamides. Fillers are treated first with the silane and then the treated filler is fluxed rapidly with polymer melt.

One of the more innovative ways of modifying the surfaces of polyolefins is to apply a multipodal oligomeric coupling agent, such as SSP-055, SSP-056, SSP-058 and SSP-255. Such oligomers provide better adhesion to polyolefins and still have the linking effect of silane dipodal chemistry in attaching to surfaces. The olefin based backbones provide great compatibility with all hydrophobic olefins as well as elastomers of various types.



Finally, an oxygen plasma treatment prior to applying silane coupling agents produces hydroxyl radicals on a polyolefin surface. These hydroxyl radicals provide good linkage sites for any silane coupling agent to link onto the polyolefin surface, and this opens a larger range of applicable silanes.




Vinylsilanes are used in PE and EPDM insulated wire and cable



Polypropylene Coupling Reaction



Dipodal Silanes



Dipodal silanes are a new series of adhesion promoters that have intrinsic hydrolytic stabilities up to $\sim 10,000$ times greater than conventional silanes. These products have a significant impact on substrate bonding and mechanical strength of many composite systems to include epoxy, urethane, epoxy/urethane hybrids, polysulfide, cyanoacrylate and silicone and may be utilized in water-borne, high solids and photo-active chemistries. Dipodal silanes are promising materials that have already achieved commercial success in applications as diverse as plastic optics, multilayer printed circuit boards and as adhesive primers for ferrous and nonferrous metals.

Due to the nature of the silicon molecules the silane coupling agent is a material used to resist deterioration by the intrusion of water between the polymer and the substrate. Through the modification of the interface, silane coupling agents not only provide water resistance, they are responsible for other important changes associated with composite systems. The interface region may exhibit increased strength because of the modification which forms interpenetrating polymer networks of resin and silane.

In silane surface treatment or 'in situ' applications, it has been the practice to hydrolyze the alkoxy groups to form silanol containing species, which are highly reactive and are responsible for hydrogen bonding with the substrate. However, it would be ideal to supply silanes with enhanced hydrolytic stability.


The problem with conventional silanes is that silanols self condense to form siloxanes resulting in phase separation or gelation. Through the addition of dipodal silanes, the enhanced hydrolytic stability will have significant impact on shelf life, substrate bonding and improved mechanical strength of many composite systems.

Functional dipodal silanes and combinations of non-functional dipodal silanes with functional conventional silanes have significant impact on substrate bonding and possess enabling activity in many adhesive systems, particularly primer and aqueous immersion applications. The fundamental step by which silanes provide adhesion is forming a $-\text{Si}-\text{O}-\text{X}$ bond with the substrate. If the substrate is siliceous, the bond durability is dictated by bond dissociation of $\text{Si}-\text{O}-\text{Si}$. According to the equation $\equiv\text{Si}-\text{O}-\text{Si}\equiv + \text{H}_2\text{O} \rightleftharpoons \equiv\text{Si}-\text{OH} + \equiv\text{Si}-\text{OH}$ the equilibrium for bond dissociation is $\sim 10^{-4}$. Recognizing that substrate hydroxyls are not subject to diffusion, the factor is closer to 10^{-2} . By increasing the number of bonds by three, the equilibrium for dissociation is increased to $\sim 10^{-6}$. Theoretically this means that dissociative bond line failure that typically occurs in 1 month is increased to $\sim 10,000$ months. Practically other factors influence the failure, but dipodal silanes clearly have the potential to eliminate failure of adhesive bonds during lifetime requirements of many devices. The effect is thought to be a result of both the increased crosslink density of the interphase and the resistance to hydrolysis of dipodal silanes, which is estimated at $\sim 10,000$ times greater than conventional coupling agents. Dipodal silanes have the ability to form six bonds to a substrate compared to conventional silanes with the ability to form only three bonds to a substrate.


Different substrates, different conditions, varying silane combinations and finally the different applications all have an effect on dipodal silane selection. The key factors determining silane-dipodal silane mixtures are:

1. Improved wet adhesion
2. Improved chemical resistance
3. Improved processing
4. Improved coating performance
(such as improved corrosion protection)

Functional Dipodal Silanes



Non-functional Dipodal Silanes



Many conventional coupling agents are frequently used in combination with 10-40% of a non-functional dipodal silane, where the conventional coupling agent provides the appropriate functionality for the application, and the non-functional dipodal silane provides increased durability. In a typical application a dipodal material such as bis(triethoxysilyl)ethane (SIB1817.0) is combined at a 1:5 to 1:10 ratio with a traditional coupling agent. It is then processed in the same way as the traditional silane coupling agent. With the addition of the non-functional dipodal silane the durability of coatings was extended when compared to the conventional silane alone.

Dipodal silanes are now commonly used in a wide variety of ways and in many diverse applications. Adding such dipodal silanes enhances hydrolytic stability, which impacts on increased product shelf life, ensures better substrate bonding and also leads to improved mechanical properties in coatings as well as composite applications.


1. Zazyczny et al in Adhesives & Sealants Industry, November 2008.

Effect of dipodal -SiCH₂CH₂Si- on the bond strength of a crosslinkable ethylene-vinyl acetate primer formulation

Primer on metal 10% in i-PrOH	Wet adhesion to metals (N/cm)	
	Titanium	Cold-rolled Steel
No silane	Nil	Nil
Methacryloxypropylsilane	0.35	7.0
Methacryloxypropylsilane + 10% dipodal	10.75	28.0 (cohesive failure)
90° peel strength after 2 h in 80°C water		

P. Pape et al, in Silanes and Other Coupling Agents, ed. K. Mittal, 1992, VSP, p105

Dipodal Silane Hydrolytic Stability compared to conventional silane



B. Arkles, et al. Chemistry - A European Journal, 2014, 20, 9442.

Cyclic Azasilanes

Efficient and high yielding, economical reactions are desired throughout chemistry. The rapid reaction of cyclic azasilanes with any and all surface hydroxyl groups is therefore of unique interest for surface modification. Volatile cyclic azasilanes afford high functional density monolayers on inorganic surfaces such as nanoparticles and other nano-featured substrates without a hydrolysis step. Furthermore, byproducts such as alcohol, HCl, and cage-like condensation products typical with the use of conventional silane coupling agents are eliminated by surface modification using cyclic azasilanes. This recently new class of silane coupling agents affords a smooth monolayer and reduces the overall waste stream as there are no longer any hazardous byproducts.

Cyclic azasilanes exploit the Si-N and Si-O bond energy differences affording a thermodynamically favorable ring-opening reaction with surface hydroxyls at ambient temperature, shown below.

Sometimes referred to as “click-chemistry on surfaces,” the ring opening occurs through the cleavage of the inherent Si-N bond in these structures, and promotes a strong covalent attachment to surface hydroxyl groups. This affords an organofunctional amine for further reactivity, depicted below, to link the inorganic surface to an organic moiety. This reaction proceeds to completion in less than a minute,

much faster compared to any conventional silane coupling agent. Importantly, the ring opening reaction depicted below also shows that the Si-OMe groups associated with traditional coupling agents remain unreacted thereby remaining available for hydrolysis and condensation reactions with other conventional silane coupling agents, should this be desired.

Cyclic azasilane coupling agents react with a wide variety of hydroxyl rich surfaces generating a range of organofunctional groups for further surface modification. Common examples of cyclic azasilanes are depicted in Figure 1. Physical properties are tabulated on page 45-46.





Figure 1: Examples of Cyclic azasilanes



Scheme 1: Reaction of one equivalent of a cyclic azasilane, SIM6501.4 (a), and a moisture cross-linking cyclic azasilane, SIB1932.4 (b), with a hydroxyl rich surface.

Cyclic azasilanes react with hydroxyl surfaces to afford a monolayer with amine functionality. The monolayers ranging from 2 to 5 nm, as measured by ellipsometry, and have an average roughness of 0.3 nm as measured by atomic force microscopy. Amine-modified surfaces are traditionally hydrophilic, and the characteristics of cyclic azasilane treated surfaces are consistently hydrophilic. The extent of the reaction is superior to the conventional counterparts as presented in Figure 2.

The rate of reaction with fumed silica can be monitored by diffuse reflectance FTIR as shown in Figure 3.

Consumption of the terminal hydroxyls (3745 cm^{-1}) occurs within 58 seconds of addition of the cyclic azasilane solution, while the C-H stretching vibrations of the $\text{Si}(\text{OMe})_2$ remain at 2864 cm^{-1} , indicating the hydrolysis of these groups, typical of conventional silane coupling agents, remain unaffected in this case, where the initial reaction is solely the breaking of the Si-N bond of the ring by the terminal surface hydroxyl groups. Additional information regarding this class of silane coupling agents can be found in the references below.

1. B. Arkles et al in "Silanes and Other Coupling Agents, Vol. 3," K. Mittal (Ed.) VSP-Brill, 2004, p 179.
2. M. Vedamuthu et al, J. Undergrad. Chem. Res., 1, 5, 2002
3. D. Brandhuber et al, J. Mater. Chem., 2005
4. Su, K. et al. U.S. Patent 2012 2672, 790, 2012




Figure 2: Extent of reaction of organosilanes with fumed silica.





Figure 3: DRIFT of untreated silica (blue) and SIB1932.4 treated silica (red) after 56 sec.


Thermal Stability of Silane Coupling Agents

The general order of thermal stability for silane coupling agents is depicted. Most commercial silane coupling agents have organic functionality separated from the silicon atom by three carbon atoms and are referred to as gamma-substituted silanes. The gamma-substituted silanes have sufficient thermal stability to withstand short-term process conditions of 350°C and long-term continuous exposure of 160°C. In some applications gamma-substituted silanes have insufficient thermal stability or other system requirements that can eliminate them from consideration. In this context, some comparative guidelines are provided for the thermal stability of silanes. Thermogravimetric Analysis (TGA) data for hydrolysates may be used for bench-marking. The specific substitution also plays a significant role in thermal stability. Electron withdrawing substitution reduces thermal stability, while electropositive groups enhance thermal stability.



Flexible multi-layer circuit boards for cell phones utilize polyimide films coupled w/chloromethylaromatic silanes.

Relative Thermal Stability of Silanes



Thermal Stability of Silanes

SIA0025.0	$\text{CH}_3\overset{\text{O}}{\parallel}\text{COCH}_2\text{CH}_2\text{Si}(\text{OC}_2\text{H}_5)_3$	220°
SIC2271.0	$\text{ClCH}_2\text{CH}_2\text{CH}_2\text{Si}(\text{OCH}_3)_3$	360°
SIM6487.4	$\text{H}_2\text{C}=\overset{\text{O}}{\parallel}\text{C}(\text{CH}_3)\text{COCH}_2\text{CH}_2\text{CH}_2\text{Si}(\text{OCH}_3)_3$	395°
SIA0591.0	$\text{H}_2\text{NCH}_2\text{CH}_2\overset{\text{H}}{\text{N}}\text{CH}_2\text{CH}_2\text{CH}_2\text{Si}(\text{OCH}_3)_3$	390°
SIA0588.0	$\text{H}_2\text{NCH}_2\text{CH}_2\text{NHCH}_2-\text{C}_6\text{H}_4-\text{CH}_2\text{CH}_2\text{Si}(\text{OCH}_3)_3$	435°
SIC2295.5	$\text{ClCH}_2-\text{C}_6\text{H}_4-\text{CH}_2\text{CH}_2\text{Si}(\text{OCH}_3)_3$	495°
SIA0599.1	$\text{H}_2\text{N}-\text{C}_6\text{H}_4-\text{Si}(\text{OC}_2\text{H}_5)_3$	485°
SIT8042.0	$\text{CH}_3-\text{C}_6\text{H}_4-\text{Si}(\text{OCH}_3)_3$	530°

25% weight loss of dried hydrolysates as determined by TGA

Aqueous Systems & Water-borne Silanes


Before most surface modification processes, alkoxy silanes are hydrolyzed forming silanol-containing species. The silanol-containing species are highly reactive intermediates which are responsible for bond formation with the substrate. In principal, if silanol species were stable, they would be preferred for surface treatments. Silanols condense with other silanols or with alkoxy silanes to form siloxanes. This can be observed when preparing aqueous treatment solutions. Initially, since most alkoxy silanes have poor solubility in water, two phases are observed. As the hydrolysis proceeds, a single clear phase containing reactive silanols forms. With aging, the silanols condense forming siloxanes and the solution becomes cloudy. Eventually, as molecular weight of the siloxanes increases, precipitation occurs.¹

Hydrolysis and condensation of alkoxy silanes is dependent on both pH and catalysts. The general objective in preparing aqueous solutions is to devise a system in which the rate of hydrolysis is substantially greater than the rate of condensation beyond the solubility limit of the siloxane oligomers. Other considerations are the work-time requirements for solutions and issues related to byproduct reactivity, toxicity or flammability.


Stable aqueous solutions of silanes are more readily prepared if byproduct or additional alcohol is present in the solution since they contribute to an equilibrium condition favoring monomeric species.

Water-borne coupling agent solutions are usually free of VOCs and flammable alcohol byproducts. Most water-borne silanes can be described as hydroxyl-rich silsesquioxane copolymers. Apart from coupling, silane monomers are included to control water-solubility and extent of polymerization. Water-borne silanes act as primers for metals, additives for acrylic latex sealants and as coupling agents for siliceous surfaces.


1. B. Arkles et al, "Factors contributing to the stability of alkoxy silanes in aqueous solutions", *J. Adhesion Science Technology*, **1992**, 6(1), 193.



Relative Hydrolysis Rates of Hydrolyzable Groups




Hydrolysis Profile of Phenylbis(2-methoxyethoxy)silanol



F. Osterholtz et al in *Silanes and Other Coupling Agents* ed K. Mittal, VSP, 1992, p119

Profile for Condensation of Silanols to Disiloxanes



E. Pohl et al in *Silanes Surfaces and Interfaces* ed., D. Leyden, Gordon and Breach, 1985, p481.

Water-borne Silsesquioxane Oligomers


Code	Functional Group	Mole %	Molecular Weight	Weight % in solution
WSA-7011	Aminopropyl	65-75	250-500	25-28
WSA-9911	Aminopropyl	100	270-550	22-25
WSA-7021	Aminoethylaminopropyl	65-75	370-650	25-28
WSAV-6511	Aminopropyl, Vinyl	60-65	250-500	25-28
WSAF-1511	Aminopropyl, Fluoroalkyl	15-20	—	15-20

Masked Silanes - Latent Functionality


Maximum bond strength in some adhesion and bonding systems requires that the organic functionality of a silane coupling agent becomes available during a discrete time period of substrate - matrix contact. Examples are epoxy adhesives in which reaction of the silane with the resin increases viscosity of an adhesive to the extent that substrate wet-out is inhibited and pretreated fillers for composites which can react prematurely with moisture before melt compounding or vulcanization. A general approach is to mask the organic functionality of the silane which converts it to a storage-stable form and then to trigger the demasking with moisture, or heat concomitant with bonding or composite formation.

Masked Silanes - Moisture Triggered

Single-component liquid-cure epoxy adhesives and coatings employ dimethylbutylidene blocked amino silanes. These materials show excellent storage stability in resin systems, but are activated by moisture provided by water adsorbed on substrate surfaces or from humidity. Deblocking begins in minutes and is generally complete within two hours in sections with a diffusional thickness of less than 1mm.




Storage Stability of Epoxy Coating Solutions with blocked and unblocked aminosilanes




Epoxy Resin Solution: 50 parts bisphenol A epoxide, 5 parts SID4068.0 or SIA0610.0, 50 parts toluene.

Hydrolysis of Blocked Aminosilane

(SID4068.0/H₂O/THF = 1/10/20wt%)




An alternative is to use the moisture adsorbed onto fillers to liberate alcohol which, in turn, demasks the organic functionality.



Masked Silanes - Heat Triggered

Isocyanate functionality is frequently delivered to resin systems during elevated temperature bonding or melt processing steps. Demasking temperatures are typically 160-200°C.




Modification of Metal Substrates

The optimum performance of silanes is associated with siliceous substrates. While the use of silanes has been extended to metal substrates, both the effectiveness and strategies for bonding to these less-reactive substrates vary. Four approaches of bonding to metals have been used with differing degrees of success. In all cases, selecting a dipodal or polymeric silane is preferable to a conventional trialkoxy silane.


Metals that form hydrolytically stable surface oxides, e.g. aluminum, tin, titanium. These oxidized surfaces tend to have sufficient hydroxyl functionality to allow coupling under the same conditions applied to the siliceous substrates discussed earlier.

Metals that form hydrolytically or mechanically unstable surface oxides, e.g. iron, copper, zinc. These oxidized surfaces tend to dissolve in water leading to progressive corrosion of the substrate or form a passivating oxide layer without mechanical strength. The successful strategies for coupling to these substrates typically involves two or more silanes. One silane is a chelating agent such as a diamine, polyamine or polycarboxylic acid. A second silane is selected which has a reactivity with the organic component and reacts with the first silane by co-condensation. If a functional dipodal or polymeric silane is not selected, 10-20% of a non-functional dipodal silane typically improves bond strength.

Metals that do not readily form oxides, e.g. nickel, gold and other precious metals. Bonding to these substrates requires coordinative bonding, typically a phosphine, sulfur (mercapto), or amine functional silane. A second silane is selected which has a reactivity with the organic component. If a functional dipodal or polymeric silane is not selected, 10-20% of a non-functional dipodal silane typically improves bond strength.




SIP6926.2



Octylsilane adsorbed on titanium

figure courtesy of
M. Banaszak-Holl

Metals that form stable hydrides, e.g. titanium, zirconium, nickel. In a significant departure from traditional silane coupling agent chemistry, the ability of certain metals to form so-called amorphous alloys with hydrogen is exploited in an analogous chemistry in which hydride functional silanes adsorb and then react with the surface of the metal.¹ Most silanes of this class possess only simple hydrocarbon substitution such as octylsilane. However they do offer organic compatibility and serve to markedly change wet-out of the substrate. Both hydride functional silanes and treated metal substrates will liberate hydrogen in the presence of base or with certain precious metals such as platinum and associated precautions must be taken.



Coupling Agents for Metals*			
Metal	Class	Screening	Candidates
Copper	Amine	SSP-060	SIT8398.0
Gold	Sulfur	SIT7908.0	SIP6926.2
	Phosphorus	SID4558.0	SIB1091.0
Iron	Amine	SIB1834.0	WSA-7011
	Sulfur	SIB1824.6	SIM6476.0
Tin	Amine	SIB1835.5	
Titanium	Epoxy	SIG5840.0	SIE6668.0
	Hydride	SIU9048.0	
Zinc	Amine	SSP-060	SIT8398.0
	Carboxylate	SIT8402.0	SIT8192.6

*These coupling agents are almost always used in conjunction with a second silane with organic reactivity or a dipodal silane.

1. B. Arkles et al J. *Adhesion Science Technol*, **2012**, 26, 41.


Difficult Substrates

Silane coupling agents are generally recommended for applications in which an inorganic surface has hydroxyl groups and the hydroxyl groups can be converted to stable oxane bonds by reaction with the silane. Substrates such as calcium carbonate, copper and ferrous alloys, and high phosphate and sodium glasses are not recommended substrates for silane coupling agents. In cases where a more appropriate technology is not available a number of strategies have been devised which exploit the organic functionality, film-forming and crosslinking properties of silane coupling agents as the primary mechanism for substrate bonding in place of bonding through the silicon atom. These approaches frequently involve two or more coupling agents.

Calcium carbonate fillers and marble substrates do not form stable bonds with silane coupling agents. Applications of mixed silane systems containing a dipodal silane or tetraethoxysilane in combination with an organofunctional silane frequently increases adhesion. The adhesive mechanism is thought to be due to the low molecular weight and low surface energy of the silanes which allows them initially to spread to thin films and penetrate porous structures followed by the crosslinking which results in the formation of a silica-rich encapsulating network. The silica-rich encapsulating network is then susceptible to coupling chemistry comparable to siliceous substrates. Marble and calciferous substrates can also benefit from the inclusion of anhydride-functional silanes which, under reaction conditions, form dicarboxylates that can form salts with calcium ions.

Metals and many metal oxides can strongly adsorb silanes if a chelating functionality such as diamine or dicarboxylate is present. A second organofunctional silane with reactivity appropriate to the organic component must be present. Precious metals such as gold and rhodium form weak coordination bonds with phosphine and mercaptan functional silanes.

High phosphate and sodium content glasses are frequently the most frustrating substrates. The primary inorganic constituent is silica and would be expected to react readily with silane coupling agents. However alkali metals and phosphates not only do not form hydrolytically stable bonds with silicon, but, even worse, catalyze the rupture and redistribution of silicon-oxygen bonds. The first step in coupling with these substrates is the removal of ions from the surface by extraction with deionized water. Hydrophobic dipodal or multipodal silanes are usually used in combination with organofunctional silanes. In some cases polymeric silanes with multiple sites for interaction with the substrate are used. Some of these, such as the polyethylenimine functional silanes can couple to high sodium glasses in an aqueous environment.



Substrates with low concentrations of non-hydrogen bonded hydroxyl groups, high concentrations of calcium, alkali metals or phosphates pose challenges for silane coupling agents.

Removing Surface Impurities

Eliminating non-bonding metal ions such as sodium, potassium and calcium from the surface of substrates can be critical for stable bonds. Substrate selection can be essential. Colloidal silicas derived from tetraethoxysilane or ammonia sols perform far better than those derived from sodium sols. Bulk glass tends to concentrate impurities on the surface during fabrication. Although sodium concentrations derived from bulk analysis may seem acceptable, the surface concentration is frequently orders of magnitude higher. Surface impurities may be reduced by immersion in 5% hydrochloric acid for 4 hours, followed by a deionized water rinse, and then immersion in deionized water overnight followed by drying.

Oxides with high isoelectric points can adsorb carbon dioxide, forming carbonates. These can usually be removed by a high temperature vacuum bake.

Increasing Hydroxyl Concentration

Hydroxyl functionalization of bulk silica and glass may be increased by immersion in a 1:1 mixture of 50% aqueous sulfuric acid : 30% hydrogen peroxide for 30 minutes followed by rinses in D.I. water and methanol and then air drying. Alternately, if sodium ion contamination is not critical, boiling with 5% aqueous sodium peroxodisulfate followed by acetone rinse is recommended¹.

1. K. Shirai et al, J. Biomed. Mater. Res. 53, 204, 2000.

Catalyzing Reactions in Water-Free Environments

Hydroxyl groups without hydrogen bonding react slowly with methoxy silanes at room temperature. Ethoxy silanes are essentially non-reactive. The methods for enhancing reactivity include transesterification catalysts and agents which increase the acidity of hydroxyl groups on the substrate by hydrogen bonding. Transesterification catalysts include tin compounds such as dibutyl diacetoxystyrene and titanates such as titanium isopropoxide. Incorporation of transesterification catalysts at 2-3 weight % of the silane effectively promotes reaction and deposition in many instances. Alternatively, amines can be premixed with solvents at 0.01-0.5 weight % based on substrate prior or concurrent to silane addition. Volatile primary amines such as butylamine can be used, but are not as effective as tertiary amines such as benzyl dimethylamine or diamines such as ethylenediamine. The more effective amines, however, are more difficult to remove after reaction¹.

1. S. Kanan et al, Langmuir, 18, 6623, 2002.

Hydroxylation by Water Plasma & Steam Oxidation

Various metals and metal oxides including silicon and silicon dioxide can achieve high surface concentrations of hydroxyl groups after exposure to H₂O/O₂ in high energy environments including steam at 1050°C and water plasma¹.

1. N. Alcanter et al, in "Fundamental & Applied Aspects of Chemically Modified Surfaces" ed. J. Blitz et al, 1999, Roy. Soc. Chem., p212.

Applying Silanes

Deposition from aqueous alcohol solutions is the most facile method for preparing silylated surfaces. A 95% ethanol-5% water solution is adjusted to pH 4.5-5.5 with acetic acid. Silane is added with stirring to yield a 2% final concentration. Five minutes should be allowed for hydrolysis and silanol formation. Large objects, e.g. glass plates, are dipped into the solution, agitated gently, and removed after 1-2 minutes. They are rinsed free of excess materials by dipping briefly in ethanol. Particles, e.g. fillers and supports, are silylated by stirring them in solution for 2-3 minutes and then decanting the solution. The particles are usually rinsed twice briefly with ethanol. Cure of the silane layer is for 5-10 mins at 110°C or 24 hours at room temperature (<60% relative humidity).

Figure 1: Vacuum tumble dryers can be used for slurry treatment of powders.



Deposition from aqueous solution is employed for most commercial fiberglass systems. The alkoxy silane is dissolved at 0.5-2.0% concentration in water. For less soluble silanes, 0.1% of a non-ionic surfactant is added prior to the silane and an emulsion rather than a solution is prepared. The solution is adjusted to pH 5.5 with acetic acid. The solution is either sprayed onto the substrate or employed as a dip bath. Cure is at 110-120°C for 20-30 minutes. Stability of aqueous silane solutions varies from 2-12 hours for the simple alkyl silanes. Poor solubility parameters limit the use of long chain alkyl and aromatic silanes by this method. Distilled water is not necessary, but water containing fluoride ions must be avoided.

Bulk deposition onto powders, e.g. filler treatment, is usually accomplished by a spray-on method. It assumes that the total amount of silane necessary is known and that sufficient adsorbed moisture is present on the filler to cause hydrolysis of the silane. The silane is prepared as a 25% solution in alcohol. The powder is placed in a high intensity solid mixer, e.g. twin cone mixer with intensifier. The methods are most effective. If the filler is dried in trays, care must be taken to avoid wicking or skinning of the top layer of treated material by adjusting heat and air flow.

Figure 2: Lab-scale Reactor for slurry treatment of powders. Separate filtration and drying steps are required.



Figure 3: Twin-cone blenders with intensive mixing bars are used for bulk deposition of silanes onto powders.



Applying Silanes


Integral blend methods are used in composite formulations. In this method the silane is used as a simple additive. Composites can be prepared by the addition of alkoxysilanes to dry-blends of polymer and filler prior to compounding. Generally 0.2 to 1.0 weight percent of silane (of the total mix) is dispersed by spraying the silane in an alcohol carrier onto a preblend. The addition of the silane to non-dispersed filler is not desirable in this technique since it can lead to agglomeration. The mix is dry-blended briefly and then melt compounded. Vacuum devolatilization of byproducts of silane reaction during melt compounding is necessary to achieve optimum properties. Properties are sometimes enhanced by adding 0.5-1.0% of tetrabutyl titanate or benzyldimethylamine to the silane prior to dispersal.

Anhydrous liquid phase deposition of chlorosilanes, methoxysilanes, aminosilanes and cyclic azasilanes is preferred for small particles and nano-featured substrates. Toluene, tetrahydrofuran or hydrocarbon solutions are prepared containing 5% silane. The mixture is refluxed for 12-24 hours with the substrate to be treated. It is washed with the solvent. The solvent is then removed by air or explosion-proof oven drying. No further cure is necessary. This reaction involves a direct nucleophilic displacement of the silane chlorines by the surface silanol. If monolayer deposition is desired, substrates should be predried at 150°C for 4 hours. Bulk deposition results if adsorbed water is present on the substrate. This method is cumbersome for large scale preparations and rigorous controls must be established to ensure reproducible results. More reproducible coverage is obtained with monochlorosilanes.

Chlorosilanes can also be deposited from alcohol solution. Anhydrous alcohols, particularly ethanol or isopropanol are preferred. The chlorosilane is added to the alcohol to yield a 2-5% solution. The chlorosilane reacts with the alcohol producing an alkoxysilane and HCl. Progress of the reaction is observed by halt of HCl evolution. Mild warming of the solution (30-40°C) promotes completion of the reaction. Part of the HCl reacts with the alcohol to produce small quantities of alkyl halide and water. The water causes formation of silanols from alkoxysilanes. The silanols condense on the substrate. Treated substrates are cured for 5-10 mins. at 110°C or allowed to stand 24 hours at room temperature.

Vapor Phase Deposition Silanes can be applied to substrates under dry aprotic conditions by chemical vapor deposition methods. These methods favor monolayer deposition. Although under proper conditions almost all silanes can be applied to substrates in the vapor phase, those with vapor pressures >5 torr at 100°C have achieved the greatest number of commercial applications. In closed chamber designs, substrates are supported above or adjacent to a silane reservoir and the reservoir is heated to sufficient temperature to achieve 5mm vapor pressure. Alternatively, vacuum can be applied until silane evaporation is observed. In still another variation the silane can be prepared as a solution in toluene, and the toluene brought to reflux allowing sufficient silane to enter the vapor phase through partial pressure contribution. In general, substrate temperature should be maintained above 50° and below 120° to promote reaction. Cyclic azasilanes deposit the quickest- usually less than 5 minutes. Amine functional silanes usually deposit rapidly (within 30 minutes) without a catalyst. The reaction of other silanes requires extended reaction times, usually 4-24 hours. The reaction can be promoted by addition of catalytic amounts of amines.

Figure 4: Apparatus for vapor phase silylation.



Applying Silanes

Spin-On Spin-On applications can be made under hydrolytic conditions which favor maximum functionalization and polylayer deposition or dry conditions which favor monolayer deposition. For hydrolytic deposition 2-5% solutions are prepared (see deposition from aqueous alcohol). Spin speed is low, typically 500 rpm. Following spin-deposition a hold period of 3-15 minutes is required before rinse solvent. Dry deposition employs solvent solutions such as methoxypropanol or ethyleneglycol monoacetate (EGMA). Aprotic systems utilize toluene or THF. Silane solutions are applied at low speed under a nitrogen purge. If strict monolayer deposition is preferred, the substrate should be heated to 50°. In some protocols, limited polylayer formation is induced by spinning under an atmospheric ambient with 55% relative humidity.

Spray Application Formulations for spray applications vary widely depending on end-use. They involve alcohol solutions and continuously hydrolyzed aqueous solutions employed in architectural and masonry applications. The continuous hydrolysis is effected by feeding mixtures of silane containing an acid catalyst such as acetic acid into a water stream by means of a venturi (aspirator). Stable aqueous solutions (see water-borne silanes), mixtures of silanes with limited stability (4-8 hours) and emulsions are utilized in textile and fiberglass applications. Complex mixtures with polyvinyl acetates or polyesters enter into the latter applications as sizing formulations.


Figure 5: Spin-coater for deposition on wafers.



Figure 6: Spray application of silanes on large structures.
















Figure 7: Spray & contact roller application of silanes on fiberglass.





Acrylate-silanes
in dental restorative
composites.

Silane Coupling Agents for Thermosets Selection Chart

	Coupling Agent Class	Suggestions for Primary Screening
Acrylate, UV cure 	Acrylate Vinyl/Olefin	SIA0200.0 SIM6487.4 SIS6964.0
Diallylphthalate 	Amine Vinyl/Olefin Dipodal	SIA0591.0 SIA0610.0 SIS6964.0 SIB1824.5
Epoxy 	Amine Anhydride Epoxy Dipodal	SIA0591.0 SIT8398.0 SIT8192.6 SIG5840.0 SIB1834.0
Epoxy, UV Cure 	Amine Epoxy	SIA0591.0 SIT8398.0 SIE4668.0 SIE4670.0
Polyimide 	Amine Halogen Dipodal	SIA0599.2 SIA0591.0 SIC2295.5 SIC2296.2 SIB1833.0
Furan 	Amine Epoxy	SIA0611.0 SIA0599.0 SIG5840.0
Melamine 	Amine Hydroxyl Dipodal	SIA0611.0 SIA0599.0 SIB1140.0 SIB1833.0 SIT8717.0
Parylene 	Halogen Vinyl/Olefinic Dipodal	SIC2295.5 SIS6990.0 SIM6487.4 SIB1832.0 VMM-010
Phenol-formaldehyde 	Amine Epoxy	SIA0611.0 SIT8187.5 SIE4670.0 SIG5840.0
Methylmethacrylate, cast 	Acrylate Amine	SIM6487.4 SIA0200.0 SIB1828.0
Polyester, unsaturated 	Acrylate Vinyl/Olefin	SIM6487.4 SIS6994.0 SIV9112.0
Urea-formaldehyde 	Amine Hydroxyl	SIA0610.0 SIU9055.0 SIB1140.0
Urethane 	Amine Isocyanate Sulfur	SIA0610.0 SIM6500.0 SII6455.0 SIM6476.0

Diamine-silanes
couple
polycarbonate
in CDs



Silane Coupling Agents for Thermoplastics Selection Chart


		Coupling Agent Class	Suggestions for Primary Screening	
Polyacetal	$\left[\text{CH}_2\text{O} \right]_n$	Vinyl/Olefin	SIS6994.0	
Polyacrylate	$\left[\begin{array}{c} \text{O} \\ \parallel \\ \text{C}-\text{OCH}_3 \\ \\ \text{CH}_2-\text{C} \\ \\ \text{H} \end{array} \right]_n$	Amine	SIU9058.0	SIA0610.0
Polyamide	$\left[\text{NH}(\text{CH}_2)_m \text{C} \begin{array}{c} \text{O} \\ \parallel \end{array} \right]_n$	Amine Dipodal Water-borne	SIA0610.0 SIB1834.1 WSA-7011	SIA0614.0 SSP-060
Polyamide-imide	$\left[\begin{array}{c} \text{O} \\ \parallel \\ \text{N} \\ \\ \text{C} \\ \\ \text{O} \end{array} \text{C}_6\text{H}_4 \begin{array}{c} \text{O} \\ \parallel \\ \text{N}-\text{R} \\ \\ \text{H} \end{array} \right]_n$	Amine Halogen	SIA0610.0 SIC2295.5	
Polybutylene terephthalate	$\left[\text{C} \begin{array}{c} \text{O} \\ \parallel \end{array} \text{C}_6\text{H}_4 \text{CO}(\text{CH}_2)_m\text{O} \right]_n$	Amine Isocyanate	SIA0610.0 SII6455.0	
Polycarbonate	$\left[\text{O}-\text{C}_6\text{H}_4-\text{C} \begin{array}{c} \text{CH}_3 \\ \\ \text{C} \\ \\ \text{CH}_3 \end{array} \text{C}_6\text{H}_4-\text{O}-\text{C} \begin{array}{c} \text{O} \\ \parallel \end{array} \right]_n$	Amine	SIA0591.0	SIA0610.0
Polyether ketone	$\left[\text{C}_6\text{H}_4-\text{O}-\text{C}_6\text{H}_4-\text{C} \begin{array}{c} \text{O} \\ \parallel \end{array} \right]_n$	Amine Dipodal	SIA0591.0 SIT8717.0	
Polyethylene	$\left[\text{CH}_2\text{CH}_2 \right]_n$	Amine Vinyl/Olefin	SIA0591.0 SSP-055	SIT8398.0 SIV9112.0
Polyphenylene sulfide	$\left[\text{C}_6\text{H}_4-\text{S} \right]_n$	Amine Halogen Sulfur	SIA0605.0 SIC2295.5 SIM6476.0	
Polypropylene	$\left[\text{CH}_2\text{CH} \begin{array}{c} \text{CH}_3 \\ \end{array} \right]_n$	Acrylate Azide Vinyl/Olefin	SIM6487.4 SIA0780.0 VEE-005	SSP-055
Polystyrene	$\left[\text{CH}_2\text{CH} \begin{array}{c} \text{C}_6\text{H}_5 \\ \end{array} \right]_n$	Acrylate Dipodal	SIM6487.4 SIB1831.0	
Polysulfone	$\left[\text{C}_6\text{H}_4-\text{C} \begin{array}{c} \text{CH}_3 \\ \\ \text{C} \\ \\ \text{CH}_3 \end{array} \text{C}_6\text{H}_4-\text{O}-\text{C}_6\text{H}_4-\text{C} \begin{array}{c} \text{O} \\ \parallel \\ \text{O} \end{array} \right]_n$	Amine	SIA0591.0	SIU9055.0
Polyvinyl butyral	$\left[\text{CH}_2-\text{C} \begin{array}{c} \text{O} \\ \\ \text{O} \\ \\ \text{CH}_2\text{CH}_2\text{CH}_3 \end{array} \right]_n$	Amine	SIA0611.0	SIU9058.0
Polyvinyl chloride	$\left[\text{CH}_2\text{CH} \begin{array}{c} \text{Cl} \\ \end{array} \right]_n$	Amine Sulfur	SIA0605.0 SIM6474.0	SIB1825.0

Water-borne aminosilanes increase bonding of acrylic latex sealants



Silane Coupling Agents for Sealants & Elastomers

Selection Chart

		Coupling Agent Class	Suggestions for Primary Screening	
Acrylic latex		Acrylate Vinyl/Olefin Water-borne	SIM6487.4 SIV9210.0 WSA-7021	SIV9218.0 WSA-6511
Butyl	$\left[\text{CH}_2\text{CH}=\text{CHCH}_2 \right]_n$	Acrylate Sulfur Vinyl/Olefin	SIM6487.4 SIB1825.0 SSP-055	SIM6476.0 VEE-005
Epichlorohydrin	$\left[\text{OCH}_2\text{CH} \left(\text{CH}_2\text{Cl} \right) \right]_n$	Amine Sulfur	SIA0605.0 SIM6474.0	
Fluorocarbon	$\text{---}(\text{CF}_2\text{CF}_2)_m(\text{CH}_2\text{CF}_2)_n\text{---}$	Amine Dipodal	SIB1834.1 SIT8717.0	
Isoprene	$\left[\text{CH}_2\text{C} \left(\text{CH}_3 \right) = \text{CHCH}_2 \right]_n$	Sulfur Vinyl/Olefin	SIM6474.0 SSP-055	SIM6476.0 VEE-005
Neoprene	$\left[\text{CH}_2\text{C} \left(\text{Cl} \right) = \text{CHCH}_2 \right]_n$	Sulfur Vinyl/Olefin	SIM6474.0 SSP-055	SIM6476.0 VEE-005
Nitrile	$\left[\text{CH}_2\text{CH} \left(\text{CN} \right) - \text{CH}_2 - \text{CH} = \text{CH} \right]_n$	Epoxy Sulfur	SIG5840.0 SIB1825.0	
Polysulfide	$\left[\text{---CH}_2\text{CH}_2\text{S---} \right]_n$	Epoxy Sulfur	SIG5840.0 SIB1825.0	SIM6476.0
SBR	$\left[\text{---CH}_2\text{CH} \left(\text{C}_6\text{H}_5 \right) - \text{CH}_2 - \text{CH} = \text{CH} \right]_n$	Amine Sulfur	SIA0605.0 SIB1825.0	SIM6486.0
Silicone (hydroxyl terminated)	$\text{HO-Si} \left(\text{CH}_3 \right)_2 \text{---O---} \left(\text{---Si} \left(\text{CH}_3 \right)_2 \text{---O---} \right)_n \text{---Si} \left(\text{CH}_3 \right)_2 \text{---OH}$	Amine Vinyl/Olefin Dipodal	SIA0605.0 SIV9098.0 SIB1824.0	SIA0589.0 VMM-010
Silicone (vinyl terminated)	$\text{H}_2\text{C}=\text{CH-Si} \left(\text{CH}_3 \right)_2 \text{---O---} \left(\text{---Si} \left(\text{CH}_3 \right)_2 \text{---O---} \right)_n \text{---Si} \left(\text{CH}_3 \right)_2 \text{---CH}=\text{CH}_2$	Acrylate Vinyl/Olefin	SIM6487.4 SIA0540.0	VMM-010



aldehyde-,
amino-, and
hydroxyl-
silanes couple
DNA in array
technology

Silane Coupling Agents for Biomaterials Selection Chart


	Site/Type	Coupling Class	Co-reactant	Suggestions for Screening	
 <p>Oligonucleotides</p>		hydroxyl		SIB1140.0	
		diamine	cobalt ethylenediamine	SIA0591.0	
		semicarbazide		SIS6944.0	
 <p>DNA</p>	terminal favored	vinyl/olefin		SIO6708.0	
	pendant amine	aldehyde		SIT8194.0	
	pendant amine	diamine		SIA0594.0	SID3543.0
	pendant amine	epoxy		SIE4675.0	SIG5838.0
 <p>Protein</p>	lysine	aldehyde		SIT8194.0	
	lysine	amine	glutaraldehyde	SIA0611.0	SIA0595.0
	lysine	amine	thiophosgene	SIA0611.0	
	cysteine	sulfur	dithionite	SIM6476.0	
	tyrosine	nitrobenzamide	NaNO ₂ /HCl	SIT8191.0	SIA0599.0
	heparinated	amine/quat		SSP-060	SIT8415.0
	immunoglobulin	pyridyl-thio		SIP6926.4	
	antibody	cyano		SIC2456.0	
	<p>J. Grobe et al, J. Chem. Soc. Chem. Commun, 2323, 1995. H. Weetall, US Pat. 3,652,761. G. Royer, CHEMTECH, 4, 699, 1974. S. Bhatia et al, Anal. Biochem., 178, 408, 1989. J. Venter et al, Proc. Nat. Acad. Soc., 69(5), 1141, 1972. R. Merker et al, Proc. Artificial Heart Prog. Conf., June 9-13, 1969 HEWNIH, p29. S. Falipou, Fundamental & Applied Aspects of Chemically Modified Surfaces, p389, 1999.</p>				
	 <p>mitochondria on silica bead</p>	chloroplast	alkyl		SIO6645.0
mitochondria		alkyl		SIO6645.0	
<p>B. Arkles et al, in "Silylated Surfaces" D. Leyden ed., Gordon & Breach, 1978, p363. B. Arkles et al, J. Biol. Chem., 250, 8856, 1975.</p>					
 <p>erythrocytes on glass wall</p>	Whole Cell	erythrocytes	short alkyl	SIE4901.4	
	<p>B. Arkles et al, in "Silylated Surfaces" D. Leyden ed., Gordon & Breach, 1978, p363.</p>				
Whole Cell (causing lysis)	procaryotic	alkyl-quat		SIO6620.0	
<p>W. White et al in "Silanes, Surfaces & Interfaces" ed. D. Leyden, Gordon & Breach, 1986, p. 107.</p>					
Tissue	histological samples			SIA0611.0	SIA0610.0



Silane Coupling Agent Properties



Adhesion promoter for structural polysulfide glass sealants



Epoxy-silanes are essential for performance of epoxy resin encapsulants for microchips.

Acrylate and Methacrylate Functional Silanes	29
Aldehyde Functional Silanes.....	32
Amino Functional Silanes	33
Cyclic Azasilanes	44
Water-borne Aminoalkyl Silsesquioxane Oligomers.....	45
Anhydride Functional Silanes	46
Azide Functional Silanes	46
Carboxylate, Phosphonate and Sulfonate Functional Silanes.....	47
Epoxy Functional Silanes	48
Ester Functional Silanes.....	49
Halogen Functional Silanes.....	50
Hydroxyl Functional Silanes.....	53
Isocyanate and Masked Isocyanate Functional Silanes	54
Phosphine and Phosphate Functional Silanes	56
Sulfur Functional Silanes.....	57
Vinyl and Olefin Functional Silanes	59
Multi-Functional and Polymeric Silanes.....	65
UV Active and Fluorescent Silanes	65
Chiral Silanes	67
Biomolecular Probes.....	68
Trihydrosilanes	69
Dipodal Silyl Hydrides.....	69
Dipodal Silanes - Non-Functional	70
Organosilane-Modified Silica Nanoparticles.....	71

Commercial Status—produced on a regular basis for inventory











Developmental Status—available to support development and commercialization









New Products—available to support development and commercialization





Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
------	----	----------	---------	------------------------------	------------------------------

Acrylate and Methacrylate Functional Silanes





Acrylate and Methacrylate Functional Silanes - Trialkoxy

	SIA0146.0 3-ACRYLAMIDOPROPYLTRIMETHOXYSILANE, tech-95 C ₉ H ₁₉ NO ₄ Si Inhibited with MEHQ [57577-96-5]	233.34		1.062	1.465	NEW
	SIA0180.0 N-(3-ACRYLOXY-2-HYDROXYPROPYL)-3-AMINOPROPYLTRIETHOXYSILANE, 50% in ethanol C ₁₅ H ₃₁ NO ₅ Si Inhibited with MEHQ Used to stabilize Stober silica suspensions. ¹ 1. Park, B. et al. <i>J. Mater. Sci.</i> 1992 , 27, 5692. [123198-57-2]	349.50	Flashpoint: 15°C (59°F)	0.931	1.4084	
	SIA0182.0 ACRYLOXYMETHYLTRIMETHOXYSILANE C ₇ H ₁₄ O ₅ Si Coupling agent, comonomer for ormosils [21134-38-3] TSCA-L	206.27	42-6° / 0.15	1.060		NEW
	SIA0184.0 (ACRYLOXYMETHYL)PHENETHYLTRIMETHOXYSILANE, tech-95 C ₁₅ H ₂₂ O ₅ Si Inhibited with MEHQ; mixed m-, p-, α-, β- isomers Coupling agent for UV cure systems with less oxygen sensitivity [141813-19-6](m-)[141813-20-9](p-)	310.42				NEW
	SIA0200.0 (3-ACRYLOXYPROPYL)TRIMETHOXYSILANE, 96% C ₉ H ₁₆ O ₅ Si Inhibited with BHT Coupling agent for epoxies, UV cure coatings Employed in optical fiber coatings. ¹ In combination with dipodal silane, SIB1833.0, increases strength and hydrolytic stability of dental composites. ² 1. Yokoshima, M. et al. <i>Chem. Abstr.</i> 113, 15746d; Jap. Pat. 02133338, 1990. 2. Matinlinna, J. et al. <i>Acta Odontol. Scand.</i> 2012 , 70, 405. [4369-14-6] TSCA	234.32	68° / 0.4 (<-20°) Flashpoint: 123°C (253°F)	1.06	1.4155	COMMERCIAL
	SIM6480.73 (3-METHACRYLAMIDOPROPYL)TRIETHOXYSILANE, tech-95 C ₁₃ H ₂₇ NO ₅ Si Contains 3-methacrylamidopropyltriethoxymethoxysilane Hydrophilic coupling agent Intermediate for contact lens monomers [109213-85-6]	289.44	Flashpoint: >110°C (>230°F)	1.017	1.454	NEW
	SIM6480.8 O-(METHACRYLOXYETHYL)-N-(TRIETHOXYSILYLPROPYL)CARBAMATE, 90% C ₁₆ H ₃₁ NO ₅ Si Inhibited with MEHQ Hydrophilic monomer [115396-93-5]	377.51		1.051 ²⁵	1.446 ²⁵	
	SIM6481.1 N-(3-METHACRYLOXY-2-HYDROXYPROPYL)-3-AMINOPROPYLTRIETHOXYSILANE, 50% in ethanol C ₁₆ H ₃₃ NO ₅ Si Inhibited with MEHQ Employed in conservation/consolidation of stone. ¹ 1. Wheeler, G. In <i>Ninth Int'l Cong. On Deterioration and Conservation of Stone</i> ; Fassina, Ed.; Elsevier: 2000; Vol. 2, 541. [96132-98-8]	363.53	Flashpoint: 11°C (52°F)	0.910	1.4084	
	SIM6482.0 METHACRYLOXYMETHYLTRIETHOXYSILANE C ₁₁ H ₂₂ O ₅ Si Inhibited with MEHQ Treatment of fumed silica in acrylic casting compositions accelerates polymerization. ¹ 1. Morozova, E. et al. <i>Chem. Abstr.</i> 95,98753g; <i>Plast. Massy</i> 1981 , 7. [5577-72-0] TSCA	262.38	65-8° / 2 Flashpoint: 100°C (212°F)	1.000	1.4225	
	SIM6483.0 METHACRYLOXYMETHYLTRIMETHOXYSILANE C ₈ H ₁₆ O ₅ Si Inhibited with MEHQ Viscosity: 1.5 cSt Hydrolysis rate > 10 X SIM6487.4 Modification of novolac resin affords bilevel resists having attributes of trilevel resists. ¹ 1. Reichmanis, E.; Smolinsky, G. U.S. Patent 4,481,049, 1984. [54586-78-6] TSCA	220.30	48-50° / 2 (-44°) Flashpoint: 92°C (198°F) TOXICITY: oral rat, LD50: >2,000 mg/kg Autoignition temperature: 285°C	1.070	1.4271	



Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰	
 SIM6487.3 METHACRYLOXYPROPYLTRIETHOXSILANE C ₁₃ H ₂₆ O ₅ Si Inhibited with MEHQ [21142-29-0]	290.43	130° / 4		0.985	1.4277	
		Flashpoint: 128°C (262°F)		10g	50g	2kg
TSCA EC 244-239-0		HMIS: 3-1-1-X store <5°C				
 SIM6487.35 METHACRYLOXYPROPYLTRIISOPROPOXSILANE C ₁₆ H ₃₂ O ₅ Si Inhibited with MEHQ, HQ [80750-05-6]	332.51			0.942	1.4234	
TSCA EC 279-538-5		HMIS: 2-2-1-X store <5°C		25g	100g	
		TOXICITY: oral rat, LD50: >5,000 mg/kg				
 SIM6487.4 METHACRYLOXYPROPYLTRIMETHOXSILANE C ₁₀ H ₂₀ O ₅ Si Inhibited with MEHQ, HQ Viscosity: 2 cSt Copolymerization parameters-e,Q: 0.07, 2.7 Specific wetting surface: 314m ² /g Coupling agent for radical cure polymer systems See SIA0200.0 for acrylate-functional UV cureable analog Widely used coupling agent for unsaturated polyester-fiberglass composites. ¹ Copolymerized with styrene in formation of sol-gel composites. ² Employed in dental polymer composites. ³ 1. Arkles, B. <i>Chemtech</i> 1977 , 7, 713. 2. Wei, Y. et al. <i>J. Mater. Res.</i> 1993 , 8, 1143. 3. Matinlinna, J. et al. <i>Int. J. Prosthodontics</i> 2004 , 17, 157. [2530-85-0]	248.35	78-81° / 1	(-48°)	1.045	1.4310	
TSCA EC 219-785-8		HMIS: 3-1-1-X store <5°C		100g	2kg	18kg
		Flashpoint: 108°C (226°F)				
		TOXICITY: oral rat, LD50: 3,000 mg/kg				
		Primary irritation index: 1.19				
 SIM6487.4LI METHACRYLOXYPROPYLTRIMETHOXSILANE, low inhibitor grade C ₁₀ H ₂₀ O ₅ Si Contains <10 ppm BHT, hydroquinones Employed in dental restorative composites 	248.35	78-81° / 1		1.045	1.4310	
TSCA EC 219-785-8		HMIS: 3-1-1-X store <5°C		25g		
 SIM6487.5 METHACRYLOXYPROPYLTRIS(METHOXYETHOXY)SILANE, tech-80 C ₁₆ H ₃₂ O ₈ Si Inhibited with MEHQ Contains methacryloxypropylbis(methoxyethoxy)methoxysilane and methacryloxypropyl(methoxyethoxy)dimethoxysilane [57069-48-4]	380.51	128° / 10		1.06		
TSCA EC 260-537-3		HMIS: 3-1-1-X store <5°C		50g		
		Flashpoint: 195°C (383°F)				
Acrylate and Methacrylate Functional Silanes - Dialkoxy						
 SIA0197.0 (3-ACRYLOXYPROPYL)METHYLDIETHOXSILANE, 95% C ₁₁ H ₂₂ O ₄ Si Inhibited with MEHQ [146666-71-9]	246.38	40° / 1		0.947	1.4139	
		Flashpoint: >110°C (>230°F)		25g		
		HMIS: 3-2-1-X store <5°C				
 SIA0198.0 (3-ACRYLOXYPROPYL)METHYLDIMETHOXSILANE, 95% C ₉ H ₁₈ O ₃ Si Inhibited with MEHQ Employed in fabrication of photoimageable, low shrinkage multimode waveguides. ¹ 1. Xu, C. et al. <i>Chem. Mater.</i> 1996 , 8, 2701. [13732-00-8]	218.33	65° / 0.35		1.01	1.431	
TSCA-L		HMIS: 3-2-1-X store <5°C		50g	250g	

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
 <p>SIM6481.43 (METHACRYLOXYMETHYL)METHYLDIETHOXSILANE C₁₀H₂₀O₄Si Inhibited with MEHQ</p>	232.35	221°		0.977	
[121177-93-3] TSCA		Flashpoint: 88°C (190°F) TOXICITY: oral rat, LD50: >2,000 mg/kg	25g		
 <p>SIM6481.46 (METHACRYLOXYMETHYL)METHYLDIMETHOXSILANE C₈H₁₆O₄Si Inhibited with MEHQ</p>	204.30	205°		1.020	1.4274
[3978-58-3]		Flashpoint: 82°C (180°F) TOXICITY: oral rat, LD50: >2,000 mg/kg Autoignition temperature: 300°C	25g	100g	
 <p>SIM6486.8 METHACRYLOXYPROPYLMETHYLDIETHOXSILANE, 95% C₁₂H₂₄O₅Si Inhibited with MEHQ</p>	260.40	95° / 1		0.965	1.4330
[65100-04-1]		Flashpoint: 136°C (277°F)	10g	50g	
 <p>SIM6486.9 METHACRYLOXYPROPYLMETHYLDIMETHOXSILANE, 95% C₁₀H₂₀O₄Si Monomer for hybrid inorganic-organic composites.¹ 1. Taylor-Smith, R. <i>Polym. Mater. Sci. Eng., Preprints</i> 1997, 77, 503.</p>	235.69	65° / 0.35		1.000	1.4355
[14513-34-9] TSCA-L EC 238-518-6		Flashpoint: 115°C (239°F)	25g	100g	2kg

Acrylate and Methacrylate Functional Silanes - Monoalkoxy


 <p>SIA0190.0 (3-ACRYLOXYPROPYL)DIMETHYLMETHOXSILANE, 95% C₉H₁₈O₃Si Inhibited with MEHQ</p>	202.32	54-5° / 0.1			
[111918-90-2]		Flashpoint: 85°C (185°F)	5g	25g	
 <p>SIM6481.3 (METHACRYLOXYMETHYL)DIMETHYLETHOXSILANE C₉H₁₈O₃Si Inhibited with MEHQ</p>	202.32	62-3° / 0.3		0.9447	1.4282 ²⁵
[5577-70-8]			10g	50g	
 <p>SIM6486.4 METHACRYLOXYPROPYLDIMETHYLETHOXSILANE, 95% C₁₁H₂₂O₃Si Inhibited with MEHQ</p>	230.38	75-6° / 0.4		0.926	1.4371
[13731-98-1]			10g		
 <p>SIM6486.5 METHACRYLOXYPROPYLDIMETHYLMETHOXSILANE, 95% C₁₀H₂₀O₃Si Inhibited with MEHQ Component in positive tone 157 nm resist.¹ 1. Tegou, E. et al. <i>Chem. Mater.</i> 2004, 16, 2567.</p>	216.35	70-2° / 0.5		0.944	1.4381
[66753-64-8] EC 266-472-7			10g	50g	

Acrylate and Methacrylate Functional Silanes - Polymeric

 <p>SIA0200.2 (3-ACRYLOXYPROPYL)TRIMETHOXSILANE, oligomeric hydrolysate</p>				1.2	
		Viscosity: 8-12 cSt Employed in rapid prototyping			
			25g	100g	
 <p>SIM6487.42 METHACRYLOXYPROPYLTRIMETHOXSILANE, oligomeric hydrolysate 1200-2500</p>				1.20	1.46
		Viscosity: 1,000-3,000 cSt Soluble: acetone, ethanol, THF			
[160185-24-0] TSCA			25g	100g	

Aldehyde Functional Silanes

Aldehyde Functional Silanes - Trialkoxy


SIT8185.3

TRIETHOXSILYL BUTYRALDEHYDE, tech-90

 $C_{10}H_{22}O_4Si$

234.37

85-71° / 1

0.96

1.414

Contains 3-TRIETHOXSILYL-2-METHYLPROPANAL isomer and cyclic siloxy acetal, 2,2,6-TRIETHOXY-1-OXA-2-SILACYCLOHEXANE

 Coupling agent for chitosan to titanium.¹


 1. Martin, H. et al. *Appl. Surf. Sci.* 2008, **254**, 4599.

[917773-12-7][88276-92-0]

TSCA-L

HMIS: 3-3-1-X

10g


SIT8194.0

TRIETHOXSILYL UNDECANAL, tech-95

 $C_{17}H_{36}O_4Si$

332.56

150-5° / 0.5

0.930

1.4343

Treated surface contact angle, water: 70°

Long chain coupling agent for DNA


 Provides greater stability for coupled proteins than shorter alkyl homologs.¹

 1. Seitz, O. et al. *J. Mater. Chem.* 2011, **21**, 4384.

[116047-42-8]

HMIS: 2-2-1-X

5g


SIT8194.5

TRIETHOXSILYL UNDECANAL, ETHYLENE GLYCOL ACETAL

 $C_{19}H_{40}O_5Si$

366.61

160-5° / 0.25

[866935-66-2]

HMIS: 2-2-1-X


5g








Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
------	----	----------	---------	------------------------------	------------------------------

Amino Functional Silanes


Monoamine Functional Silanes - Trialkoxy

	<p>SIA0587.0 4-AMINOBYLTRIETHOXYSILANE, 95% C₁₀H₂₅NO₃Si</p> <p>Amino-functional coupling agent [3069-30-5]</p>	235.40	114-6° / 14 Flashpoint: 109°C (228°F) TOXICITY: oral rat, LD50: 1,620 mg/kg	0.941 ²⁵ 1.4270 ²⁵	
	<p>SIA0587.07 4-AMINO-3,3-DIMETHYLBUTYLTRIMETHOXYSILANE AMINONEOHEXYLTRIMETHOXYSILANE C₉H₂₃NO₃Si</p> <p>Sterically hindered primary amine coupling agent Non-yellowing aminosilane coupling agent for flexible adhesives and sealants [157923-74-5]</p>	221.37	230° Flashpoint: 97°C (207°F)	0.977 1.4302	COMMERCIAL
	<p>SIA0590.5 N-(2-AMINOETHYL)-3-AMINOPROPYLTRIETHOXYSILANE, 95% C₁₁H₂₈N₂O₃Si</p> <p>Coupling agent with slower hydrolysis than SIA0591.0, SAI0592.6 [5089-72-5]</p>	264.55	156° / 15 Flashpoint: 148°C (298°F)	0.994 1.4367 ²⁵	
	<p>SIA0598.0 3-(m-AMINOPHENOXY)PROPYLTRIMETHOXYSILANE, tech-95 C₁₂H₂₁NO₃Si</p> <p>Amber liquid High temperature coupling agent [55648-29-8]</p>	271.39	125-135° / 0.5	1.02 1.495	
	<p>SIA0599.0 m-AMINOPHENYLTRIMETHOXYSILANE, 95% C₉H₁₅NO₃Si</p> <p>Contains other isomers Coupling agent for polyimides used in electronics [70411-42-6]</p>	213.31	110-4° / 0.6 Flashpoint: 180°C (356°F)	1.19 1.5187	
	<p>SIA0599.1 p-AMINOPHENYLTRIMETHOXYSILANE, 95% C₉H₁₅NO₃Si</p> <p>Contains other isomers Coupler for silica-poly(phenyleneterephthalamide) composite films.¹ Together with phenyltrimethoxysilane, SIP6822.0, can be used to increase the dispersibility of mesoporous silica.² 1. Mark, J. et al. <i>J. Mater. Chem.</i> 1997, 7, 259. 2. Banet, P. et al. <i>Langmuir</i> 2008, 24, 9030. [33976-43-1]</p>	213.31	110-4° / 0.6 (60-62°) Flashpoint: 180°C (356°F)		
	<p>SIA0599.2 AMINOPHENYLTRIMETHOXYSILANE C₉H₁₅NO₃Si</p> <p>Mixed isomers, typically 60-70% para, 30-40% meta Coupling agent for polyimides used in electronics [33976-43-1]</p>	213.31	110-4° / 0.6 (60-62°) Flashpoint: 180°C (356°F)	1.19	NEW
	<p>SIA0610.0 3-AMINOPROPYLTRIETHOXYSILANE GAPS, AMEO, A-1100 C₉H₂₃NO₃Si</p> <p>Viscosity: 1.6 cSt ΔHvap: 11.8 kcal/mole Treated surface contact angle, water: 59° yc of treated surfaces: 37.5 mN/m Specific wetting surface: 353 m²/g</p> <p>Widely used coupling agent for phenolic, epoxy, polyamide and polycarbonate resins Effects immobilization of enzymes.¹ Used to bind Cu(salicylaldehyde) to silica.²</p>	221.37	122-3° / 30 Flashpoint: 104°C (219°F) TOXICITY: oral rat, LD50: 1,780 mg/kg Primary irritation index: 6.50 Vapor pressure, 100°: 10 mm	0.951 1.4225	
<p>Aminosilanes are widely used as coupling agents for fiberglass insulation and composites</p>	<p>1. <i>Enzymes</i> 1976, 84, 55915. 2. Murphy, E. F. et al. <i>Inorg. Chem.</i> 2003, 42, 2559. [919-30-2]</p>	25g	25g	2kg 16kg	COMMERCIAL





Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
 <p>SIA0610.1 3-AMINOPROPYLTRIETHOXSILANE 99+% C₉H₂₃NO₃Si</p> <p>Low fluorescence grade for high throughput screening [919-30-2] TSCA EC 213-048-4 HMIS: 3-1-1-X</p> <p>* in fluoropolymer bottle</p>	221.37	122-3° / 30 Flashpoint: 104°C (219°F)		0.951	1.4225
 <p>SIA0611.0 3-AMINOPROPYLTRIMETHOXSILANE C₆H₁₇NO₃Si</p> <p>Coupling agent with superior reactivity in vapor phase and non-aqueous surface treatments Hydrolysis rate vs SIA0610.0: 6:1 Used to immobilize Cu and Zn Schiff base precatalysts for formation of cyclic carbonates.¹ 1. Tasci, Z.; Ulusoy, M. <i>J. Organomet. Chem.</i> 2012, 713, 104. [13822-56-5] TSCA EC 237-511-5 HMIS: 3-2-1-X</p>	179.29	80° / 8 Flashpoint: 83°C (181°F) TOXICITY: oral rat, LD50: >2,000 mg/kg Autoignition temperature: 300°C Vapor pressure, 67°: 5 mm		1.027	1.4240
 <p>SIA0614.0 3-AMINOPROPYLTRIS(METHOXYETHOXYETHOXY)SILANE, 95% C₁₈H₄₁NO₉Si</p> <p>Coupling agent for melt compounding of polyamides and epoxides [87994-64-7] EC 289-365-7 HMIS: 3-2-1-X</p>	443.61	Flashpoint: 68°C (154°F)		1.066	1.448
 <p>SIA0630.0 11-AMINOUNDECYLTRIETHOXSILANE C₁₇H₃₉NO₃Si</p> <p>Contains ~ 5% isomers [116821-45-5] HMIS: 2-2-1-X</p>	333.59	128-32° / 1		0.895 ²⁵	1.4352 ²⁵
 <p>SIP6928.0 2-(4-PYRIDYLETHYL)TRIETHOXSILANE C₁₃H₂₃NO₃Si</p> <p>Amber liquid Forms self-assembled layers which can be "nano-shaved" by scanning AFM.¹ 1. Rosa, L. et al. <i>Mater. Lett.</i> 2009, 63, 961. [98299-74-2] HMIS: 3-2-1-X</p>	269.43	105° / 0.9		1.00	1.4624 ²⁴
 <p>SIP6930.0 2-(2-PYRIDYLETHYL)TRIMETHOXSILANE 2-(TRIMETHOXSILYLETHYL)PYRIDINE C₁₀H₁₇NO₃Si</p> <p>[27326-65-4] HMIS: 3-1-1-X</p>	227.33	105° / 0.3 Flashpoint: >110°C (>230°F)		1.06	1.4755
 <p>SIT8410.0 N-(3-TRIMETHOXSILYLPROPYL)PYRROLE C₁₀H₁₉NO₃Si</p> <p>For electrode modification, polypyrrole adhesion.¹ 1. Simon, R. et al. <i>J. Am. Chem. Soc.</i> 1982, 104, 2031. [80906-67-8] HMIS: 3-1-1-X</p>	229.35	105-7° / 1 Flashpoint: >110°C (>230°F)		1.017	1.463

Monoamine Functional Silanes - Water-borne


 <p>SIA0608.0 3-AMINOPROPYLSILANETRIOL, 22-25% in water C₃H₁₁NO₃Si</p> <p>Mainly oligomers; monomeric at concentrations <5% pH: 10.0-10.5 Water-borne, VOC-free coupling agent Internal hydrogen bonding stabilizes solution See also WSA-7011 for greater hydrolytic stability [58160-99-9]/[29159-37-3] TSCA EC 261-145-5 HMIS: 2-0-0-X</p>	137.21	Flashpoint: >110°C (>230°F)		1.06	
--	--------	-----------------------------	--	------	--


Monoamine Functional Silanes - Dialkoxy


 <p>SIA0587.05 4-AMINO-3,3-DIMETHYLBUTYLMETHYLDIMETHOXSILANE C₉H₂₃NO₂Si</p> <p>Comonomer for non-yellowing textile lubricants [156849-43-3] TSCA HMIS: 3-2-1-X</p>	205.37	Flashpoint: 100°C (212°F) TOXICITY: oral rat, LD50: 631 mg/kg		0.9248 ²⁵	1.4365
--	--------	--	--	----------------------	--------

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰	
 <p>SIA0605.0 3-AMINOPROPYLMETHYLDIETHOXSILANE, 95% C₈H₂₁N₂O₂Si</p>	191.34	85-8° / 8		0.916	1.4272	COMMERCIAL
Flashpoint: 85°C (185°F) TOXICITY: oral rat, LD50: 4,760 mg/kg Coupling agent for foundry resins, including phenolic novolaks and resols Vapor phase deposition >150° on silica yields high density amine functionality. ¹ 1. Ek, S. et al. <i>Langmuir</i> 2003 , <i>19</i> , 3461. [3179-76-8] TSCA EC 221-660-8 HMIS: 3-2-1-X 25g 2kg 15kg						


Monoamine Functional Silanes - Monoalkoxy


 <p>SIA0587.1 1-AMINO-2-(DIMETHYLETHOXSILYL)PROPANE, 85% C₇H₁₉NOSi</p>	161.32	101-7° / 100		0.85	1.431	NEW
Contains 3-aminopropyl dimethylethoxysilane HMIS: 3-2-1-X 5g						


 <p>SIA0602.0 3-AMINOPROPYLDIISOPROPYLETHOXSILANE C₁₁H₂₇NOSi</p>	217.43	78-80° / 0.4		0.872 ²⁵	1.4489	
Forms hydrolytically stable amino-functional bonded phases and monolayers [117559-36-1] HMIS: 3-2-0-X 5g 25g						

 <p>SIA0603.0 3-AMINOPROPYLDIMETHYLETHOXSILANE C₇H₁₉NOSi</p>	161.32	78-9° / 24		0.857 ²⁵	1.4276 ²⁵	
ΔHform: 147.6 kcal/mole Flashpoint: 73°C (163°F) Coupling agent for DNA array technology [18306-79-1] TSCA HMIS: 3-2-1-X 5g 25g						










Diamine Functional Silanes - Trialkoxy










 <p>SIA0588.0 (AMINOETHYLAMINOMETHYL)PHENETHYLTRIMETHOXSILANE, tech-90 C₁₄H₂₆N₂O₃Si</p>	298.46	126-30° / 0.2		1.02	1.5083	
Mixed <i>m</i> -, <i>p</i> - isomers Flashpoint: >110°C (>230°F) Coupling agent for polyimides. ¹ Photochemically sensitive (194nm). ² Forms self assembled monolayers. ³ Reagent for charge heterogeneity in micropatterning. ⁴ 1. Arkles, B. et al. <i>Modern Plastics</i> 1980 , <i>57</i> , 64. 2. Dressick, W. et al. <i>Thin Solid Films</i> 1996 , <i>284</i> , 568. 3. Harnett, C. et al. <i>Appl. Phys. Lett.</i> 2000 , <i>76</i> , 2466. 4. Chen, J. et al. <i>Nano Lett.</i> 2002 , <i>2</i> , 393. [74113-77-2] TSCA HMIS: 3-1-1-X 25g 100g						

 <p>SIA0591.0 N-(2-AMINOETHYL)-3-AMINOPROPYLTRIMETHOXSILANE, tech-95 N-[3-(TRIMETHOXSILYL)PROPYL]ETHYLENEDIAMINE; DAMO C₈H₂₂N₂O₃Si</p>	222.36	140° / 15		1.019 ²⁵	1.450 ²⁵	COMMERCIAL
Contains 2-6% N,N'-bis[(trimethoxysilyl)propyl]ethylene-diamine Flashpoint: 150°C (302°F) Viscosity: 6.5 cSt TOXICITY: oral rat, LD50: 7,460 mg/kg Surface tension: 36.5 mN/m Autoignition temperature: 300°C Specific wetting surface: 358 m ² /g Coefficient of thermal expansion: 0.8 x 10 ⁻³ Coupling agent for polyamides, polycarbonates (e.g. in CDs), polyesters and copper/brass adhesion Film-forming coupling agent/primer, fiberglass size component For cyclic version see SID3543.0, for pre-hydrolyzed version see SIA0590.0 Used in the immobilization of copper (II) catalyst on silica. ¹ Used together w/ SID3396.0 to anchor PdCl ₂ catalyst to silica for acceleration of the Tsuji-Trost reaction. ² 1. Wu, Q.; Wang, L. <i>Synthesis</i> 2008 , 2007. 2. Noda, H. et al. <i>Angew. Chem., Int. Ed. Engl.</i> 2012 , <i>51</i> , 8017. [1760-24-3] TSCA EC 217-164-6 HMIS: 3-1-1-X 25g 2kg 16kg						


 <p>SIA0591.1 N-(2-AMINOETHYL)-3-AMINOPROPYLTRIMETHOXSILANE, 98% N-[3-(TRIMETHOXSILYL)PROPYL]ETHYLENEDIAMINE; DAMO C₈H₂₂N₂O₃Si</p>	222.36	140° / 15		1.019 ²⁵	1.450 ²⁵	COMMERCIAL
Flashpoint: 150°C (302°F) [1760-24-3] TSCA EC 217-164-6 HMIS: 3-1-1-X 25g 100g 2kg						

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
 SIA0592.6 N-(6-AMINOHEXYL)AMINOMETHYLTRIETHOXSILANE, 95% C ₁₃ H ₃₂ N ₂ O ₃ Si 292.49 160° / 0.1 Flashpoint: >110°C (>230°F) TOXICITY: oral rat, LD50: 500 mg/kg Long-chain amino coupling agent [15129-36-9] HMIS: 3-2-1-X 25g 100g	0.928 ²⁵	1.4385 ²⁵			
 SIA0594.0 N-(6-AMINOHEXYL)AMINOPROPYLTRIMETHOXSILANE, 95% C ₁₂ H ₃₀ N ₂ O ₃ Si 278.47 160-5° / 0.15 Flashpoint: >110°C (>230°F) Employed in immobilization of DNA. ¹ Employed for immobilization of PCR primers on beads. ² 1. Kneuer, C. et al. <i>Int'l. J. Pharmaceutics</i> 2000 , <i>196</i> , 257. 2. Andreadis, J. et al. <i>Nuc. Acid Res.</i> 2000 , <i>28</i> , E-5. [51895-58-0] HMIS: 3-1-1-X 10g 50g	1.11	1.4501			
 SIA0595.0 N-(2-AMINOETHYL)-11-AMINOUNDECYLTRIMETHOXSILANE C ₁₆ H ₃₈ N ₂ O ₃ Si 334.57 155-9° / 0.4 Coupling agent with extended spacer-group for remote substrate binding [121772-92-7] HMIS: 3-1-1-X 5g	0.873 ²⁵	1.4515			
 SIA0599.4 N-3-[(AMINO(POLYPROPYLENOXY))AMINOPROPYLTRIMETHOXSILANE, 60 - 65% 337-435 Contains amine-terminated polypropylene oxide 3-4 propylenoxy units Coupling agent with film-forming capability. HMIS: 2-2-1-X 25g	0.984	1.4508			
 SIB0956.0 N-(2-N-BENZYLAMINOETHYL)-3-AMINOPROPYLTRIMETHOXSILANE, tech-90 C ₁₅ H ₂₈ N ₂ O ₃ Si 312.48 Contains aminoethylaminopropyltrimethoxysilane Flashpoint: 69°C (156°F) [209866-89-7] TSCA HMIS: 3-2-1-X 25g 100g	1.035	1.4902			NEW
Diamine Functional Silanes - Water-borne					
 SIA0590.0 N-(2-AMINOETHYL)-3-AMINOPROPYLSILANETRIOL, 25% in water, mainly oligomers C ₅ H ₁₆ N ₂ O ₃ Si 180.28 pH: 10.0-10.5 Flashpoint: >110°C (>230°F) Internal hydrogen bonding stabilizes solution Additive for CMP slurries Aqueous primer, adhesion promoter for resin-to-metal applications See also WSA-7021 for greater hydrolytic stability [68400-09-9] TSCA HMIS: 2-0-0-X 100g 2kg 18kg	1.00				COMMERCIAL
Diamine Functional Silanes - Polymeric					
 SIA0591.3 N-(2-AMINOETHYL)-3-AMINOPROPYLTRIMETHOXSILANE-PROPYLTRIMETHOXSILANE, oligomeric co-hydrolysate C ₉ H ₂₄ N ₂ O ₅ Si 1.09 1.442 Flashpoint: >110°C (>230°F) TSCA HMIS: 3-2-1-X 25g 100g					NEW
Diamine Functional Silanes - Dialkoxy					
 SIA0587.5 N-(2-AMINOETHYL)-3-AMINOISOBUTYLMETHYLDIMETHOXSILANE, 95% C ₉ H ₂₄ N ₂ O ₂ Si 220.39 131° / 15 Flashpoint: 96°C (205°F) Amino-functional coupling agent [23410-40-4] TSCA EC 245-642-4 HMIS: 3-2-1-X 25g	0.960	1.4518			
 SIA0588.8 N-(2-AMINOETHYL)-3-AMINOPROPYLMETHYLDIETHOXSILANE C ₁₀ H ₂₆ N ₂ O ₂ Si 234.41 108-110° / 1.5 Flashpoint: >110°C (>230°F) Adhesion promoter for silanol-functional silicones on metal substrates [70240-34-5] EC 274-494-3 HMIS: 3-2-1-X 25g	0.923	1.445			NEW

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
 <p>SIA0589.0 N-(2-AMINOETHYL)-3-AMINOPROPYLMETHYLDIMETHOXYSILANE, tech-95 C₈H₂₂N₂O₂Si Specific wetting surface: 380 m²/g</p> <p>Comonomer for silicones in textile softeners and hair care formulations Coupling agent for furan-quartz sand floor coating systems Adhesion promoter for urea-formaldehyde binders on flexible substrates [3069-29-2] TSCA EC 221-336-6 HMIS: 3-1-1-X</p>	206.36	265°		0.975 ²⁵	1.4447 ²⁵
COMMERCIAL					
Diamine Functional Silanes - Monoalkoxy					
 <p>SIA0587.2 N-(2-AMINOETHYL)-3-AMINOISOBUTYLDIMETHYLMETHOXYSILANE, 95% C₉H₂₄N₂O₂Si Amino-functional coupling agent [31024-49-4]</p>	204.39	85-9° / 2		0.900 ²⁵	1.4513 ²⁵
COMMERCIAL					
Triamine Functional Silanes					
 <p>SIT8398.0 (3-TRIMETHOXYSILYL)PROPYL)DIETHYLENTRIAMINE, tech-95 C₁₀H₂₇N₃O₃Si yc of treated surface: 37.5 mN/m</p> <p>Hardener, coupling agent for epoxies [35141-30-1] TSCA EC 252-390-9 HMIS: 3-1-1-X</p>	265.43	114-8° / 2		1.030	1.4590
COMMERCIAL					
Secondary Amine Functional Silanes					
 <p>SIA0400.0 3-(N-ALLYLAMINO)PROPYLTRIMETHOXYSILANE, 95% C₉H₂₁NO₃Si</p> <p>Coupling agent for polyesters Coupling agent for acrylic coatings for glass containers.¹ 1. Hashimoto. Y. et al. Eur. Pat. Appl. EP 289,325, 1988. [31024-46-1] TSCA EC 250-435-7 HMIS: 3-2-1-X</p>	219.36	106-9° / 25		0.989 ²⁵	1.4990 ²⁵
COMMERCIAL					
 <p>SIB1932.2 n-BUTYLAMINOPROPYLTRIMETHOXYSILANE C₁₀H₂₅NO₃Si</p> <p>Reacts with isocyanate resins (urethanes) to form moisture cureable systems [31024-56-3] TSCA EC 250-437-8 HMIS: 2-1-1-X</p>	235.40	102° / 3.5		0.947	1.4246 ²⁵
COMMERCIAL					
 <p>SIB1932.3 t-BUTYLAMINOPROPYLTRIMETHOXYSILANE C₁₀H₂₅NO₃Si [174219-86-4]</p> <p>HMIS: 2-2-1-X</p>	235.40	98-9° / 3		0.924	1.4208
NEW					
 <p>SIC2464.16 (N-CYCLOHEXYLAMINOMETHYL)METHYLDIETHOXYSILANE, 95% C₁₂H₂₇NO₂Si</p> <p>End-cap modifier for moisture-cure urethane systems (SPUR) [27445-54-1] HMIS: 2-2-1-X</p>	245.40	235°		0.93	
NEW					
 <p>SIC2464.2 (N-CYCLOHEXYLAMINOMETHYL)TRIETHOXYSILANE, 95% C₁₃H₂₉NO₃Si [26495-91-0]</p> <p>TSCA EC 247-744-4 HMIS: 2-1-1-X</p>	275.46	236°		0.950	1.4377
NEW					
 <p>SIC2464.4 (N-CYCLOHEXYLAMINOPROPYL)TRIMETHOXYSILANE C₁₂H₂₇NO₃Si Viscosity: 5-7 cSt [3068-78-8]</p> <p>TSCA EC 221-329-8 HMIS: 3-2-1-X</p>	261.43	114° / 3		0.99	1.486 ²⁵
NEW					


Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰	
 <p>SIC2058.2 3-CARBAZOLYLPROPYLTRIETHOXSILANE C₂₁H₂₉NO₃Si For non-linear optic materials Employed in OLED fabrication.¹ 1. DeMais, T. et al. <i>SPIE Proc.</i> 1998, 3476, 338 [221105-38-0]</p>	371.55	185-195° / 0.3		1.072	1.5527 ²⁵	NEW
 <p>SID3395.4 (N,N-DIETHYLAMINOMETHYL)TRIETHOXSILANE C₁₁H₂₇NO₃Si Catalyst for neutral cure 1-part RTVs [15180-47-9] TSCA-L</p>	249.43	74-6° / 3		0.9336 ²⁵	1.4142 ²⁵	COMMERCIAL
 <p>SID3395.6 (N,N-DIETHYLAMINOMETHYL)TRIMETHOXSILANE, 95% C₉H₂₁NO₃Si Charge control agent for toner particles Crosslinker for moisture-cure silicone RTVs [67475-66-5] TSCA-L</p>	207.40			0.95	1.415	NEW
 <p>SID3396.0 (N,N-DIETHYL-3-AMINOPROPYL)TRIMETHOXSILANE C₁₀H₂₅NO₃Si Provides silica-supported catalyst for 1,4-addition reactions.¹ Used together w/ SIA0591.0 to anchor PdCl₂ catalyst to silica for acceleration of the Tsuji-Trost reaction.² 1. Mutukura, K. et al. <i>Chem.-Eur. J.</i> 2009, 15, 10871. 2. Noda, H. et al. <i>Angew. Chem., Int. Ed. Engl.</i> 2012, 51, 8017. [41051-80-3] TSCA EC 255-192-0</p>	235.40	120° / 20	Flashpoint: 100°C (212°F)	0.934	1.4245	COMMERCIAL
 <p>SID3546.92 3-(N,N-DIMETHYLAMINOPROPYL)AMINOPROPYLMETHYLDIMETHOXSILANE C₁₁H₂₆N₂O₂Si Combines secondary and tertiary amine functionality Comonomer for silicone textile finishes [224638-27-1]</p>	284.44	92-4° / 0.7		0.915	1.442	NEW
 <p>SID3546.94 N,N-DIMETHYL-3-AMINOPROPYLMETHYLDIMETHOXSILANE C₈H₂₁NO₂Si [67353-42-8]</p>	191.36	92° / 25		0.894	1.4203	NEW
 <p>SID3547.0 (N,N-DIMETHYL-3-AMINOPROPYL)TRIMETHOXSILANE C₈H₂₁NO₃Si Derivatized silica catalyzes Michael reactions.¹ 1. Mode, J. et al. <i>Synlett</i> 1998, 625. [2530-86-1] TSCA EC 219-786-3</p>	207.34	106° / 30	Flashpoint: 99°C (210°F)	0.948 ²⁵	1.4150	
 <p>SIM6572.0 N-METHYL-N-(TRIMETHYLSILYL)-3-AMINOPROPYLTRIMETHOXSILANE, 95% C₁₀H₂₇NO₃Si₂ Contains N-methylaminopropyltrimethoxysilane</p>	265.50					NEW
 <p>SIT8716.2 TRIS(TRIETHOXSILYLMETHYL)AMINE, tech-90 C₂₁H₅₁NO₃Si₃ Contains ~5% bis(trimethylsilylmethyl)amine Forms immobilized quaternary salts [1250435-76-7]</p>	545.90	156° / 0.7		0.99	1.426	NEW

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
------	----	----------	---------	------------------------------	------------------------------




SIT8716.3 TRIS(TRIETHOXYSILYLPROPYL)AMINE, tech-95 C ₂₇ H ₆₃ NO ₉ Si ₃ Coupling agent/primer for metal substrates [18784-74-2]	630.06	200-5° / 1		1.4322 ²⁵	NEW
	HMIS: 2-2-1-X	25g			


Quaternary Amine Functional Silanes




SIB0957.0 N-(2-N-BENZYLAMINOETHYL)-3-AMINOPROPYLTRIMETHOXYSILANE hydrochloride, 50% in methanol C ₁₅ H ₂₈ N ₂ O ₃ Si·HCl Amber liquid [623938-90-9]	348.95		Flashpoint: 9°C (48°F)	0.942	1.4104
	TSCA	HMIS: 3-3-1-X	25g	100g	




SID3392.0 N,N-DIDECYL-N-METHYL-N-(3-TRIMETHOXYSILYLPROPYL)AMMONIUM CHLORIDE, 40-42% in methanol C ₂₇ H ₆₀ ClNO ₃ Si Contains 3-5% Cl(CH ₂) ₃ Si(OMe) ₃ In combination with TEOS forms high pore volume xerogels with adsorptive capacity. ¹ 1. Markovitz, M. et al. <i>Langmuir</i> 2001 , <i>17</i> , 7085. [68959-20-6]	510.32		Flashpoint: 11°C (52°F)	0.863	1.4085
	TSCA	EC 273-403-4	HMIS: 3-4-0-X	25g	




SIO6620.0 OCTADECYLDIMETHYL(3-TRIMETHOXYSILYLPROPYL)AMMONIUM CHLORIDE, 60% in methanol C ₂₆ H ₅₈ ClNO ₃ Si Contains 3-5% Cl(CH ₂) ₃ Si(OMe) ₃ Employed as a glass lubricant Orients liquid crystals Provides an antistatic surface coating Dispersion/coupling agent for high density magnetic recording media. ¹ Application as immobilizable antimicrobial reported. ² 1. Vincent, H. In <i>Chemically Modified Oxide Surfaces</i> ; D. Leyden, D., Ed.; Gordon & Breach: 1990; p.305. 2. White, W. et al. In <i>Silanes, Surfaces & Interfaces</i> ; Leyden, D., Ed.; Gordon & Breach: 1986; p.107. [27668-52-6]	496.29		Flashpoint: 15°C (59°F) Autoignition temperature: 230°C	0.89	
	TSCA	EC 248-595-8	HMIS: 3-4-0-X	25g	2kg 16kg




SIS6989.0 (STYRYLMETHYL)BIS(TRIETHOXYSILYLPROPYL)AMMONIUM CHLORIDE, 40% in ethanol C ₂₇ H ₅₂ ClNO ₃ Si ₂ Inhibited with BHT, mixed <i>m</i> -, <i>p</i> -isomers Dipodal quaternary coupling agent [49377-00-3]	578.34		Flashpoint: 15°C (59°F)	0.909	
	HMIS: 3-4-1-X	store <5°C	25g	100g	











SIS6994.0 3-(N-STYRYLMETHYL-2-AMINOETHYLAMINO)PROPYLTRIMETHOXYSILANE HYDROCHLORIDE, 40% in methanol C ₁₇ H ₃₁ ClN ₂ O ₃ Si Inhibited with BHT Viscosity: 2.3 cSt Specific wetting surface area: 208 m ² /g Coupling agent for phenolic and epoxy fiberglass laminates (printed circuit boards) [34937-00-3]	374.98		Flashpoint: 11°C (52°F)	0.91	1.395
	TSCA	EC 252-297-3	HMIS: 3-4-1-X	store <5°C	25g 2kg 15kg













SIT7090.0 TETRADECYLDIMETHYL(3-TRIMETHOXYSILYLPROPYL)AMMONIUM CHLORIDE, 50% in methanol C ₂₂ H ₅₀ ClNO ₃ Si Contains 3-5% Cl(CH ₂) ₃ Si(OMe) ₃ [41591-87-1]	440.18		Flashpoint: 11°C (52°F)	0.878	1.3971
	TSCA	EC 255-451-8	HMIS: 3-4-0-X	25g	

















SIT8395.0 4-(TRIMETHOXYSILYLETHYL)BENZYLTRIMETHYLAMMONIUM CHLORIDE, 60% in methanol C ₁₅ H ₂₈ ClNO ₃ Si Candidate for exchange resins and extraction phases [49377-00-3]	333.93		Flashpoint: 11°C (52°F)	0.966	
	HMIS: 3-3-1-X	25g			

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
 <p>SIT8405.0 S-(TRIMETHOXYSILYL)PROPYL)ISOTHIOURONIUM CHLORIDE, 50% in water TRIHYDROXYPROPYL CARBAMIDOTHIOIC ACID HYDROCHLORIDE C₇H₁₉ClN₂O₃SSi Essentially silanetriol pH: 6 Antimicrobial activity reported [84682-36-0] TSCA EC 283-599-3 HMIS: 2-0-0-X 25g</p>	274.84			1.190	1.441
 <p>SIT8415.0 N,N,N-TRIMETHOXYSILYL)PROPYL-N,N,N-TRIMETHYLAMMONIUM CHLORIDE, 50% in methanol N,N,N-TRIMETHYL-3-(TRIMETHOXYSILYL)-1-PROPANAMINIUM CHLORIDE C₉H₂₄ClNO₃Si Flashpoint: 11°C (52°F) Employed for bonded chromatographic phases Anti-static agent Used to treat glass substrates employed in electroblotting Prevents contact electrification.¹ 1. Thomas, S. et al. <i>J. Am. Chem. Soc.</i> 2009, <i>131</i>, 8746. [35141-36-7] TSCA EC 252-393-5 HMIS: 2-4-1-X 25g 2kg 15kg</p>	257.83			0.927	1.3966
COMMERCIAL					
Dipodal Amine Functional Silanes					
 <p>SIA0587.6 1-[3-(2-AMINOETHYL)-3-AMINOISOBUTYL]-1,1,3,3,3-PENTAETHOXY-1,3-DISILAPROPANE, 95% C₁₇H₄₂N₂O₅Si₂ Pendant dipodal silane HMIS: 3-1-1-X 10g</p>	410.70	130-140° / 0.5		0.990	
NEW					
 <p>SIB1620.0 BIS(METHYLDIETHOXYSILYL)PROPYL)AMINE, 95% C₁₆H₃₉NO₄Si₂ Dipodal coupling agent [31020-47-0] HMIS: 2-1-1-X 25g 100g</p>	365.66	155° / 0.6		0.937	1.4385
 <p>SIB1645.0 BIS(METHYLDIMETHOXYSILYL)PROPYL)-N-METHYLAMINE, 95% C₁₃H₃₂NO₄Si₂ Viscosity: 6 - 7 cSt HMIS: 3-2-1-X 25g</p>	323.58	140° / 2		0.951	1.4368
 <p>SIB1824.5 BIS(3-TRIETHOXYSILYL)PROPYL)AMINE, 95% C₁₈H₄₃NO₅Si₂ Viscosity: 5.5 cSt Flashpoint: 162°C (324°F) Coupling agent for polyamides with improved hydrolytic stability Adhesion promoter, crosslinking agent for hot melt adhesives Adhesion promoter for aluminum-polyester multilayer laminates Adhesion promoter, crosslinker for 2-part condensation cure silicones [13497-18-2] TSCA EC 236-818-1 HMIS: 3-1-1-X 25g 100g 2kg</p>	425.71	160° / 0.6		0.97	1.4265
COMMERCIAL					
 <p>SIB1828.0 N,N'-BIS[3-(TRIETHOXYSILYL)PROPYL]UREA, 60% in ethanol C₁₉H₄₄N₂O₇Si₂ [69465-84-5] HMIS: 2-1-1-X 25g 100g</p>	468.73	Flashpoint: 17°C (63°F) 25g		0.923	
 <p>SIB1832.8 1,11-BIS(TRIMETHOXYSILYL)-4-OXA-8-AZAUNDECAN-6-OL, 50% in methanol 3,3,15,15-TETRAMETHOXY-2,7,16-TRIOXA-11-AZA-3,15-DISILAHEPTADECAN-9-OL C₁₅H₃₇NO₈Si₂ Flashpoint: 11°C (52°F) Provides strong glass-to-glass bonds [126552-46-3] HMIS: 3-4-1-X 25g 100g</p>	415.63	Flashpoint: 11°C (52°F)		0.922	1.384 ²⁵
NEW					





Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰	
<p>SIB1833.0 BIS(3-TRIMETHOXY-SILYL-PROPYL)AMINE, 96% C₁₂H₃₁N₂O₆Si₂ 341.56 Flashpoint: 113°C (235°F)</p> <p>Dipodal coupling agent Secondary amine allows more control of reactivity with isocyanates Low level incorporation with acryloxypropyltrimethoxysilane, SIA0200, increases strength and hydrolytic stability of dental composites.¹ 1. Matinlinna, J. et al. <i>Acta Odontol. Scand.</i> 2012, <i>70</i>, 405. [82985-35-1] TSCA EC 280-084-5 HMIS: 3-1-1-X 25g 2kg 18kg</p>				1.040	1.4320	COMMERCIAL
<p>SIB1834.0 N,N'-BIS[(3-TRIMETHOXY-SILYL)PROPYL]ETHYLENEDIAMINE, 62% in methanol C₁₄H₃₆N₂O₆Si₂ 384.62 Flashpoint: 20°C (68°F)</p> <p>Contains N,N-isomer Coupling agent for polyamides with enhanced hydrolytic stability Provides improved solder resistance for printed circuit boards [68845-16-9] TSCA EC 272-453-4 HMIS: 3-4-1-X 25g 2kg 16kg</p>				0.89		COMMERCIAL
<p>SIB1834.1 N,N'-BIS[(3-TRIMETHOXY-SILYL)PROPYL]ETHYLENEDIAMINE, 95% C₁₄H₃₆N₂O₆Si₂ 384.62 Flashpoint: >110°C (>230°F)</p> <p>Contains N,N-isomer Coupling agent for polyamides with enhanced hydrolytic stability Forms thin film environments for metal ions.¹ 1. He, J. et al. <i>RIKEN Review</i> 2002, <i>45</i>, 27. [68845-16-9] TSCA EC 272-453-4 HMIS: 3-1-1-X 10g 50g</p>				1.050	1.4428	
<p>SIB1835.0 BIS(3-TRIMETHOXY-SILYL-PROPYL)-N-METHYLAMINE C₁₃H₃₃N₂O₆Si₂ 355.58 Flashpoint: 175° / 10 Flashpoint: 106°C (223°F)</p> <p>[31024-70-1] HMIS: 2-1-0-X 25g 100g</p>				1.023	1.430	NEW
<p>SIB1835.2 N,N'-BIS(3-TRIMETHOXY-SILYL-PROPYL)THIOUREA, tech-90 C₁₄H₃₆N₂O₆Si₂ 400.66 Flashpoint: 113°C (235°F)</p> <p>Viscosity: 160-200 cSt Forms mesoporous silicas that chelate metals HMIS: 3-2-1-X 10g</p>				1.13	1.485	NEW
<p>SIB1835.5 N,N'-BIS(3-TRIMETHOXY-SILYL-PROPYL)UREA, 95% C₁₃H₃₂N₂O₆Si₂ 384.58 Flashpoint: >110°C (>230°F)</p> <p>Amber liquid Viscosity: 100 - 250 cSt Adhesion promoter for 2-part condensation cure silicone RTVs [18418-53-6] TSCA HMIS: 3-2-1-X 25g 100g 2kg</p>				1.10	1.4488	COMMERCIAL
<p>SIS6989.0 (STYRYLMETHYL)BIS(TRIETHOXY-SILYL-PROPYL)AMMONIUM CHLORIDE, 40% in ethanol C₂₇H₅₂ClNO₆Si₂ 578.34 Flashpoint: 15°C (59°F)</p> <p>Inhibited with BHT, mixed <i>m</i>-, <i>p</i>-isomers Dipodal quaternary coupling agent HMIS: 3-4-1-X store <5°C 25g 100g</p>				0.909		
Specialty Amine Functional Silanes						
<p>SIB1835.2 N,N'-BIS(3-TRIMETHOXY-SILYL-PROPYL)THIOUREA, tech-90 C₁₃H₃₂N₂O₆SSi₂ 400.66 Flashpoint: 113°C (235°F)</p> <p>Viscosity: 160-200 cSt Forms mesoporous silicas that chelate metals HMIS: 3-2-1-X 10g</p>				1.13	1.485	NEW
<p>SID4068.0 3-(1,3-DIMETHYLBUTYLIDENE)AMINOPROPYLTRIETHOXY-SILANE, tech-95 C₁₅H₃₃NO₃Si 303.52 Flashpoint: 134° / 5 Flashpoint: 131°C (268°F)</p> <p>Contains oligomers Coupling agent for epoxy coatings; blocked amine - moisture deblocked Preferred adhesion promoter for low viscosity epoxy systems [116229-43-7] TSCA HMIS: 2-1-1-X 25g 100g 2kg</p>				0.93	1.437 ²⁵	COMMERCIAL

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰	
 <p>SID4068.1 3-(1,3-DIMETHYLBUTYLIDENE)AMINOPROPYLTRIETHOXSILANE, 98% C₁₅H₃₃NO₃Si 303.52</p> <p>Coupling agent for epoxy coatings; blocked amine - moisture deblocked Preferred adhesion promoter for low viscosity epoxy systems [116229-43-7] TSCA HMIS: 2-1-1-X 25g</p> <p>Flashpoint: 131°C (268°F)</p> <p>0.93 1.437²⁵ NEW</p>						
 <p>SID4465.0 N,N-DIOCTYL-N'-TRIETHOXSILYLPROPYLUREA C₂₆H₅₆N₂O₄Si 488.83</p> <p>Forms hydrophobic phases with embedded hydrophilicity Forms organic-inorganic vesicles (cerasomers).¹ 1. Hashizume, M. et al. <i>J. Thin Solid Films</i> 2003, 438, 20. [259727-10-1] HMIS: 2-2-1-X 25g</p> <p>0.924²⁵ 1.4521²⁵</p>						
 <p>SIG5840.12 3-(GUANIDINYL)PROPYLTRIMETHOXSILANE C₇H₁₉N₃O₃Si 221.35</p> <p>White powder Water-soluble coupling agent [128310-21-4] HMIS: 3-1-1-X 5g</p> <p>(>200°)</p> <p>NEW</p>						
 <p>SIP6828.4 3-(1-PIPERAZINYL)PROPYLMETHYLDIMETHOXSILANE, 95% C₁₀H₂₄N₂O₂Si 232.40</p> <p>Adhesion promoter for metal substrates Comonomer for silicones [128996-12-3] HMIS: 3-2-1-X 25g</p> <p>110-2° / 1 Flashpoint: 123°C (253°F)</p> <p>0.986 1.4628²⁵ NEW</p>						
 <p>SIP6926.2 3-(2-PYRIDYLETHYL)THIOPROPYLTRIMETHOXSILANE C₁₃H₂₃NO₃SSi 301.48</p> <p>[29098-72-4] HMIS: 3-2-1-X 10g</p> <p>156-7° / 0.25</p> <p>1.089 1.498</p>						
 <p>SIP6926.4 3-(4-PYRIDYLETHYL)THIOPROPYLTRIMETHOXSILANE, 95% C₁₃H₂₃NO₃SSi 301.48</p> <p>[198567-47-4] HMIS: 3-2-1-X 10g</p> <p>160-2° / 0.2</p> <p>1.09 1.5037</p> <p>Immobilizable ligand for immunoglobulin IgG separation using hydrophobic charge induction chromatography (HCIC)</p>						
 <p>SIS6944.0 3-(4-SEMICARBAZIDYL)PROPYLTRIETHOXSILANE, tech-95 C₁₀H₂₅N₃O₄Si 279.41</p> <p>Employed in immobilization of oligonucleotides.¹ 1. Podyminogin, M. et al. <i>Nucleic Acids Res.</i> 2001, 29, 5090. [106868-88-6] HMIS: 3-1-1-X 5g</p> <p>Flashpoint: >110°C (>230°F)</p> <p>1.08 1.4593 NEW</p>						
 <p>SIS6995.0 11-(SUCCINIMIDYLOXY)UNDECYLDIMETHYLETHOXSILANE, 95% C₁₉H₃₇NO₄Si 371.59</p> <p>Reagent for immobilization of proteins via primary amines HMIS: 3-2-1-X 1.0g</p> <p>195-200° / 0.6 (28°)</p>						
 <p>SIT8186.45 4-(TRIETHOXSILYLPROPOXY)-2,2,6,6-TETRAMETHYLPYPERIDINE N-OXIDE, tech-85 <i>TEMPO-SILANE</i> C₁₈H₃₈NO₃Si 376.58</p> <p>Employed in nitroxyl mediated grafting of vinylsilanes to polyolefins.¹ 1. Weaver, J. et al. <i>J. Polym. Sci., Part A: Polym. Chem.</i> 2008, 46, 4542. [913375-27-6] HMIS: 2-2-1-X 10g</p> <p>NEW</p>						
 <p>SIT8186.7 N-[3-(TRIETHOXSILYL)PROPYL]-2-CARBOMETHOXYAZIRIDINE, 95% C₁₃H₂₇NO₅Si 305.44</p> <p>Reacts with fullerenes by 1,3-dipolar addition of azomethine ylide to yield pyrrolidine adduct.¹ 1. Bianco, A. et al. <i>J. Am. Chem. Soc.</i> 1997, 119, 7550. [193417-26-4] HMIS: 3-2-1-X 2.5g</p> <p>1.052²⁵ 1.4197²⁵ NEW</p>						



Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰	
 <p>SIT8187.5 N-(3-TRIETHOXSILYLPROPYL)-4,5-DIHYDROIMIDAZOLE 3-(2-IMIDAZOLIN-1-YL)PROPYLTRIETHOXSILANE, IMEO C₁₂H₂₆N₂O₃Si</p>	274.43	134° / 2		1.005	1.452	COMMERCIAL
Viscosity: 5 cSt. Flashpoint: >110°C (>230°F) Coupling agent for elevated temperature-cure epoxies Utilized in HPLC of metal chelates. ¹ Forms proton vacancy conducting polymers with sulfonamides by sol-gel. ² Ligand for molecular imprinting of silica with chymotrypsin transition state analog. ³ 1. Suzuki, T. et al. <i>Chem. Lett.</i> 1994 , 881. 2. De Zea Bermudez, V. et al. <i>Sol-Gel Optics II, SPIE Proc.</i> 1992 , 1728, 180. 3. Markowitz, M. et al. <i>Langmuir</i> 2000 , 16, 1759.						
[58068-97-6]	TSCA	EC 261-093-3	HMIS: 2-1-1-X	25g	100g	2kg
 <p>SIT8394.0 N-[5-(TRIMETHOXSILYL)-2-AZA-1-OXOPENTYL]CAPROLACTAM, 95% N-TRIMETHOXSILYLPROPYLCARBAMOYL-CAPROLACTAM C₁₃H₂₆N₂O₅Si</p>	318.45		(-39°)	1.14	1.4739	
Flashpoint: 136°C (277°F) [106996-32-1] HMIS: 3-1-1-X 25g						
 <p>SIT8409.0 N-(3-TRIMETHOXSILYLPROPYL)PERFLUOROHEXANAMIDE C₁₂H₁₆F₁₁NO₄Si</p>	475.33			1.367	1.360	NEW
Surface tension: 20.6 mN/m Contact angle, water on treated silica surface: 105-110° [154380-34-4] HMIS: 3-1-1-X 10g						
 <p>SIU9055.0 UREIDOPROPYLTRIETHOXSILANE, 50% in methanol C₁₀H₂₄N₂O₄Si</p>	264.40		(-97°)	0.92	1.386	COMMERCIAL
Contains ureidopropyltrimethoxysilane and related transesterification products Flashpoint: 11°C (52°F) Coupling agent for polyamides and urea-formaldehyde resins [23779-32-0] TSCA EC 245-876-7 HMIS: 2-4-1-X 25g 2kg 16kg						
 <p>SIU9058.0 UREIDOPROPYLTRIMETHOXSILANE C₇H₁₈N₂O₄Si</p>	222.32	217-225°	(-5°)	1.150	1.386 ²⁵	COMMERCIAL
Flashpoint: 99°C (210°F) TOXICITY: oral rat, LD50: >5,000 mg/kg Autoignition temperature: 300°C Component in primers for tin alloys Adhesion promoter for foundry resins [23843-64-3] TSCA EC 245-904-8 HMIS: 2-2-1-X 25g 100g 2kg						
Cyclic Azasilanes						
 <p>SIA0415.0 N-ALLYL-AZA-2,2-DIMETHOXSILACYCLOPENTANE C₈H₁₇NO₂Si</p>	187.31	52-4° / 3				
Coupling agent for nanoparticles [618914-49-1] HMIS: 3-3-1-X 10g						
 <p>SIA0592.0 N-(2-AMINOETHYL)-2,2,4-TRIMETHYL-1-AZA-2-SILACYCLOPENTANE C₈H₂₀N₂Si</p>	172.35	54-6° / 2		0.905	1.4769	
Coupling agent for vapor phase modification of nanoparticles [18246-33-8] HMIS: 3-2-1-X 10g						
 <p>SIA0604.0 N-(3-AMINOPROPYLDIMETHYLSILYL)AZA-2,2-DIMETHYL-2-SILACYCLOPENTANE, 95% 3-[(2,2-DIMETHYL-1-AZA-2-SILACYCLOPENT-1-YL)DIMETHYLSILYL]-1-PROPANAMINE, 3-[(2,2-DIMETHYL-1,2-AZASILOLIDIN-1-YL)DIMETHYLSILYL]-1-PROPANAMINE C₁₀H₂₆N₂Si₂</p>	230.50	158°			1.4705	
Flashpoint: 72°C (162°F) Employed in vapor phase derivatization of porous silica. ¹ 1. Brandhuber, D. et al. <i>J. Mater. Chem.</i> 2005 , 15, 3896. [388606-32-4] HMIS: 3-2-1-X 10g						
 <p>SIB1932.4 N-n-BUTYL-AZA-2,2-DIMETHOXSILACYCLOPENTANE C₉H₂₁NO₂Si</p>	203.36	63-71° / 3		0.941	1.438	
Flashpoint: 85°C (185°F) Coupling agent for nanoparticles. ¹ Interlayer bonding agent for anti-reflective lenses. ² 1. Arkles, B. et al. In <i>Silanes and Other Coupling Agents</i> ; Mittal, K., Ed.; 2004; Vol. 3, p.179. 2. Su, K. et al. U.S. Patent Appl. 2012 2672,790, 2012. [618914-44-6] TSCA HMIS: 3-2-1-X 25g						

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
 SID3543.0 2,2-DIMETHOXY-1,6-DIAZA-2-SILACYCLOOCTANE C ₇ H ₁₈ N ₂ O ₂ Si Volatile coupling agent [182008-07-7]	190.32	71-3° / 2.5	(61-2°)		
 SID3543.1 2,2-DIMETHOXY-1,6-DIAZA-2-SILACYCLOOCTANE, 10% in cyclohexane C ₇ H ₁₈ N ₂ O ₂ Si [182008-07-7]	190.32			0.79	NEW
 SID3546.93 (N,N-DIMETHYLAMINOPROPYL)-AZA-2-METHYL-2-METHOXY-SILACYCLOPENTANE, 95% C ₁₀ H ₂₄ N ₂ O ₂ Si Cyclic azasilane - acts as coupling agent for metals [182008-07-7]	216.39	65-7° / 0.5		0.919 ²⁵	NEW
 SIE4891.0 1-ETHYL-2,2-DIMETHOXY-4-METHYL-1-AZA-2-SILACYCLOPENTANE C ₈ H ₁₉ NO ₂ Si Cyclic azasilane utilized for vapor-phase surface treatment of nano-featured substrates [182008-07-7]	189.33				NEW
 SIT8187.2 (1-(3-TRIETHOXY-SILYL)PROPYL)-2,2-DIETHOXY-1-AZA-2-SILACYCLOPENTANE, tech-90 C ₁₆ H ₃₇ NO ₉ Si ₂ Contains bis(triethoxysilylpropyl)amine Cyclic azasilane that acts as coupling agent for metal-oxide films [1184179-50-7]	379.64	136-8° / 1		0.974	1.4322 ²⁵ NEW


Water-borne Aminoalkyl Silsesquioxane Oligomers

 WSA-7011 AMINOPROPYLSILSESQUIOXANE IN AQUEOUS SOLUTION Viscosity: 5-15 cSt Mole % functional group: 65-75 pH: 10-10.5 [1411854-75-5]	250-500			1.10	COMMERCIAL
 WSA-9911 AMINOPROPYLSILSESQUIOXANE IN AQUEOUS SOLUTION Viscosity: 5-15 cSt Mole % functional group: 100 pH: 10-10.5 [29159-37-3]	270-550			1.06	COMMERCIAL
 WSA-7021 AMINOETHYLAMINOPROPYLSILSESQUIOXANE IN AQUEOUS SOLUTION Viscosity: 5-10 cSt Mole % functional group: 65-75 pH: 10-11 [1411854-75-5]	370-650			1.10	COMMERCIAL
 WSAV-6511 AMINOETHYLAMINOPROPYL/VINYL/SILSESQUIOXANE IN AQUEOUS SOLUTION Viscosity: 3-10 cSt Mole % functional group: 60-65 pH: 10-11 [207308-27-8]	250-500			1.11	COMMERCIAL

Polyamine


 SSP-060 TRIMETHOXSILYLPROPYL MODIFIED (POLYETHYLENIMINE), 50% in isopropanol 1,500-1,800 Flashpoint: 12°C (54°F) [136856-91-2] / [37251-86-8]				0.92	
 SSP-065 DIMETHOXSILYLMETHYLPROPYL MODIFIED (POLYETHYLENIMINE), 50% in isopropanol 1,500-1,800 Flashpoint: 12°C (54°F) [125441-88-5]				0.92	

Anhydride Functional Silanes




SIT8192.6 (3-TRIETHOXSILYL)PROPYLSUCCINIC ANHYDRIDE, 95% C ₁₃ H ₂₄ O ₆ Si Viscosity: 20 cSt	304.41	135° / 0.2		1.070	1.4405	COMMERCIAL
Coupling agent for dibasic surfaces Acetic acid-catalyzed hydrolysis yields succinic acid derivatives [93642-68-3]	TSCA EC 297-566-6 HMIS: 2-1-1-X	25g	100g	2kg		


Azide Functional Silanes




SIA0770.0 (AZIDOMETHYL)PHENETHYLTRIMETHOXSILANE, tech-90 C ₁₂ H ₁₉ N ₃ O ₃ Si Mixed isomers CAUTION: CAN FORM EXPLOSIVE COMPOUNDS IN CONTACT WITH COPPER AND SILVER COMPOUNDS Reagent for "click" chemistry [1245946-78-4]	281.39			1.11	1.497 ²⁵	NEW
Used with click chemistry to introduce and immobilize discrete complexes onto the SBA-15 surface. ¹ Used in the preparation of poly-L-lysine to silica nanoparticles. ² 1. Nakazawa, J. et al. <i>J. Am. Chem. Soc.</i> 2012 , <i>134</i> , 2750. 2. Kar, M. et al. <i>Langmuir</i> , 2010 , <i>26</i> , 5772. [83315-69-9]	HMIS: 4-2-1-X	5g				




SIA0774.0 p-AZIDOMETHYLPHENYLTRIMETHOXSILANE, 90% C ₁₀ H ₁₅ N ₃ O ₃ Si AVOID CONTACT WITH METALS [83315-74-6]	253.33			1.14	1.497 ²⁵	NEW
	HMIS: 4-2-1-X	2.5g				




SIA0777.0 3-AZIDOPROPYLTRIETHOXSILANE C ₉ H ₂₁ N ₃ O ₃ Si AVOID CONTACT WITH METALS Used with click chemistry to introduce and immobilize discrete complexes onto the SBA-15 surface. ¹ Used in the preparation of poly-L-lysine to silica nanoparticles. ² 1. Nakazawa, J. et al. <i>J. Am. Chem. Soc.</i> 2012 , <i>134</i> , 2750. 2. Kar, M. et al. <i>Langmuir</i> , 2010 , <i>26</i> , 5772. [83315-69-9]	247.37	72-5° / 0.5		0.981	1.4257	NEW
	HMIS: 4-1-1-X	1.0g				












SIA0780.0 6-AZIDOSULFONYLHEXYLTRIETHOXSILANE, tech-95 1-TRIETHOXSILYL-6-SULFONAZIDE-n-HEXANE C ₁₂ H ₂₇ N ₃ O ₅ SSi AVOID CONTACT WITH METALS Amber-brown liquid Flashpoint: 114°C (237°F) Inserts nitrenes into aliphatics and aromatics at temperatures >110°C [96550-26-4]	353.51			1.147	1.4634	NEW
	HMIS: 3-1-1-X	25g				





SIA0790.0 4-(AZIDOSULFONYL)PHENETHYLTRIMETHOXSILANE, 22-25% in methanol/toluene 4-[2-(TRIMETHOXSILYL)ETHYL]-1-BENZENESULFONYL AZIDE C ₁₁ H ₁₇ N ₃ O ₅ SSi AVOID CONTACT WITH METALS Contains hydrolysis oligomers Extremely reactive coupling agent Flashpoint: 29°C (4°F) Inserts nitrenes into aliphatics and aromatics at temperatures >110°C [68479-60-7]	331.42			0.90	1.550	NEW
	TSCA EC 270-862-2 HMIS: 3-4-1-X	25g				











SIA0795.0 11-AZIDOUNDECYLTRIMETHOXSILANE, 95% C ₁₄ H ₃₁ N ₃ O ₃ Si AVOID CONTACT WITH METALS Forms "click" functionalized surfaces [334521-23-2]	317.50	111° / 0.05				NEW
	HMIS: 3-1-1-X	2.5g				

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
 <p>SIB1824.56 BIS(3-TRIETHOXSILYLPROPYL)CARBONATE C₁₉H₄₂O₉Si₂ Forms thermally labile bridged silsesquioxanes.¹ 1. Loy, D.A. et al. <i>Chem. Mater.</i> 1999, <i>11</i>, 3333. [88321-11-3]</p>	470.71			1.020	
 <p>SIB1834.5 BIS(3-TRIMETHOXSILYLPROPYL) FUMARATE, 96% C₁₆H₃₂O₁₀Si₂ Dipodal silane Adhesion promoter for Pt-cure silicone RTVs [3371-62-8] TSCA</p>	440.59		Flashpoint: >110°C (>230°F)	1.118	1.443
 <p>SIC2263.0 CARBOXYETHYLSILANETRIOL, DISODIUM SALT, 25% in water C₃H₆Na₂O₅Si pH: 12 - 12.5 In combination with aminofunctional silanes forms amphoteric silicas.¹ 1. Han, L. et al. <i>Chem. Mater.</i> 2007, <i>19</i>, 2860. [18191-40-7]</p>	196.14			1.170 ²⁵	
 <p>SIC2417.0 2-(4-CHLOROSULFONYLPHENYL)ETHYLTRIMETHOXSILANE, 50% in methylene chloride C₁₁H₁₇ClO₅SSi Amber color Contains free sulfonic acid Treated silica acts as etherification catalyst.¹ Reagent for surface initiated ATRP.² Employed in mesostructured fuel-cell membranes.³ 1. Sow, B. et al. <i>Microporous and Mesoporous Mat'ls.</i> 2005, <i>79</i>, 129. 2. Fukuda, J. et al. <i>Macromolecules</i> 2000, <i>33</i>, 2870. 3. Pereira, F. et al. <i>Chem. Mater.</i> 2008, <i>20</i>, 1710. [126519-89-9]</p>	324.85			1.30 ²⁵	
 <p>SIC2417.4 2-(4-CHLOROSULFONYLPHENYL)ETHYLTRIMETHOXSILANE, 50% in toluene C₁₁H₁₇ClO₅SSi Flashpoint: 4°C (39°F) 25g [126519-89-9]</p>	324.85				
 <p>SIT8189.8 TRIETHOXSILYLPROPYLMALEAMIC ACID, tech-90 C₁₃H₂₅NO₆Si Contains condensation products and 10% ethanol Viscosity: 600-900 cSt Reagent for immobilization of anti-bodies on silicon nitride sensor chips.¹ 1. Kurihara, Y. et al. <i>Langmuir</i> 2012, <i>28</i>, 13609. [33525-68-7] TSCA EC 251-564-1</p>	319.43			1.11	1.472
 <p>SIT8192.2 TRIETHOXSILYLPROPYL(POLYETHYLENEOXY)PROPYLPOTASSIUM SULFATE, 45-50% in ethanol C₂₅H₅₃KO₁₃SSi - C₃₁H₆₅KO₁₆SSi Flashpoint: 15°C (59°F) 2.5g HMIS: 2-3-0-X</p>	661 - 793			0.95	
 <p>SIT8378.1 TRIHYDROXSILYLETHYL PHENYLSULPHONIC ACID, 25% in water C₈H₁₂O₆SSi Mixed isomers Forms hybrid organic-inorganic ionic membranes [143282-00-2][876338-08-8]</p>	264.33				
 <p>SIT8378.3 3-(TRIHYDROXSILYL)-1-PROPANESULFONIC ACID, 30-35% in water C₃H₁₀O₆SSi pH: <1 Employed in preparation of nanoscale ionic silicas.¹ Forms sulfonated fuel cell electrodes.² 1. Giannelis, E. et al. <i>Appl. Organomet. Chem.</i> 2010, <i>24</i>, 581. 2. Eastcott, J. et al. <i>J. Power Sources</i> 2012, <i>197</i>, 102. [70942-24-4] TSCA</p>	202.26			1.12	

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
 <p>SIT8401.0 N-(TRIMETHOXYSILYLPROPYL)ETHYLENEDIAMINETRIACETATE, TRIPOTASSIUM SALT, 30% in water C₁₄H₂₅K₃N₂O₉Si Essentially silanetriol, contains KCl Chelates metal ions [1309595-29-6]</p>	510.75			1.22	NEW
 <p>SIT8402.0 N-(TRIMETHOXYSILYLPROPYL)ETHYLENEDIAMINETRIACETATE, TRISODIUM SALT, 35% in water C₁₄H₂₅Na₃O₉Si Essentially silanetriol, contains NaCl Chelates metal ions [128850-89-5][1138444-28-6]</p>	462.42			1.26	




Epoxy Functional Silanes

Epoxy Functional Silanes - Trialkoxy


	<p>SIE4668.0 2-(3,4-EPOXYCYCLOHEXYL)ETHYLTRIETHOXYSILANE C₁₄H₂₈O₄Si Adhesion promoter for water-borne coatings on alkaline substrates [10217-34-2]</p>	288.46	114-7° / 0.4 Flashpoint: 104°C (219°F)	1.015	1.4455	COMMERCIAL
	<p>SIE4670.0 2-(3,4-EPOXYCYCLOHEXYL)ETHYLTRIMETHOXYSILANE C₁₁H₂₂O₄Si Viscosity: 5.2 cSt γc of treated surfaces: 39.5 mN/m Coefficient of thermal expansion: 0.8 x 10⁻³</p>	246.38	95-7° / 0.25 Flashpoint: 146°C (295°F) TOXICITY: oral rat, LD50: 12,300 mg/kg Vapor pressure, 152°: 10 mm Specific wetting surface: 317 m ² /g	1.065	1.4490	COMMERCIAL
 <p>Epoxy-silanes are essential for performance of epoxy resin encapsulants for microchips</p>	<p>Ring epoxide more reactive than glycidoxypropyl systems UV initiated polymerization of epoxy group with weak acid donors Forms UV-curable coating resins by controlled hydrolysis.¹ Used to make epoxy-organosilica particles w/ high positive Zeta potential.² 1. Crivello, J.; Mao. <i>Z. Chem. Mater.</i> 1997, 9, 1554. 2. Nakamura, M.; Ishimura, K. <i>Langmuir</i> 2008, 24, 12228.</p>	[3388-04-3]	TSCA EC 222-217-1 HMIS: 3-1-1-X	100g	2kg	18kg
	<p>SIE4675.0 5,6-EPOXYHEXYLTRIETHOXYSILANE C₁₂H₂₆O₄Si [86138-01-4]</p>	262.42	115-9° / 1.5 Flashpoint: 99°C (210°F)	0.960 ²⁵	1.4254 ²⁵	
	<p>SIG5839.0 (3-GLYCIDOXYPROPYL)TRIETHOXYSILANE C₁₂H₂₆O₅Si Viscosity: 3 cSt</p>	278.42	124° / 3 Flashpoint: 144°C (291°F) Autoignition temperature: 225°C	1.00	1.425	COMMERCIAL
	<p>SIG5840.0 (3-GLYCIDOXYPROPYL)TRIMETHOXYSILANE 3-(2,3-EPOXYPROPOXY)PROPYLTRIMETHOXYSILANE C₉H₂₀O₅Si Viscosity: 3.2 cSt Surface tension: 38.5 mN/m Specific wetting surface area: 331 m²/g</p>	236.34	120° / 2 (<-70°) Flashpoint: 135°C (276°F) TOXICITY: oral rat, LD50: 8,400 mg/kg Autoignition temperature: 231°C	1.070	1.4290	COMMERCIAL
	<p>Component in aluminum metal bonding adhesives Coupling agent for epoxy composites employed in electronic "chip" encapsulation Component in abrasion resistant coatings for plastic optics Used to prepare epoxy-containing hybrid organic-inorganic materials.¹ 1. Innocenzi, P. et al. <i>Chem. Mater.</i> 1999, 11, 1672.</p>	[2530-83-8]	TSCA EC 219-784-2 HMIS: 3-1-1-X	100g	2kg	18kg
	<p>SIG5840.1 (3-GLYCIDOXYPROPYL)TRIMETHOXYSILANE, 99+% C₉H₂₀O₅Si</p>	236.34	120° / 2 (<-70°) Flashpoint: 135°C (276°F)	1.070	1.4290	
	<p>Low fluorescence grade for high-throughput screening [2530-83-8]</p>	TSCA EC 219-784-2 HMIS: 3-1-1-X	25g			* in fluoropolymer bottle

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
------	----	----------	---------	------------------------------	------------------------------







Epoxy Functional Silanes - Dialkoxy




	SIE4666.0 2-(3,4-EPOXYCYCLOHEXYL)ETHYLMETHYLDIETHOXSILANE C ₁₃ H ₂₆ O ₃ Si	258.43	114-7° / 1	0.976 ²⁵	1.4248 ²⁵	NEW
	UV polymerizeable monomer [14857-35-3]	HMIS: 2-2-1-X	25g			
	SIG5832.0 (3-GLYCIDOXYPROPYL)METHYLDIETHOXSILANE C ₁₁ H ₂₄ O ₄ Si	248.39	122-6° / 5 Flashpoint: 122°C (252°F) TOXICITY: oral rat, LD50: >2,000 mg/kg	0.978 ²⁵	1.431	COMMERCIAL
	Employed in scratch resistant coatings for eye glasses Coupling agent for latex systems with reduced tendency to gel compared to SIG5840.0 [2897-60-1]	TSCA EC 220-780-8 HMIS: 2-1-1-X	25g	100g	2kg	
	SIG5836.0 (3-GLYCIDOXYPROPYL)METHYLDIMETHOXSILANE C ₉ H ₂₀ O ₄ Si	220.34	100° / 4 Flashpoint: 105°C (221°F)	1.02	1.431 ²⁵	
	Relative hydrolysis rate vs. SIG5840.0: 7.5:1 [65799-47-5]	TSCA EC 265-929-8 HMIS: 3-1-1-X	25g	100g		

Epoxy Functional Silanes - Monoalkoxy

	SIG5825.0 (3-GLYCIDOXYPROPYL)DIMETHYLETHOXSILANE C ₁₀ H ₂₂ O ₃ Si	218.37	100° / 3 Flashpoint: 87°C (189°F)	0.950	1.4337 ²⁵	
	[17963-04-1]	TSCA EC 241-889-7 HMIS: 3-2-1-X	10g	50g		








Ester Functional Silanes











	SIA0050.0 ACETOXYMETHYLTRIETHOXSILANE C ₉ H ₂₀ O ₅ Si	236.34	106° / 15	1.042 ²⁵	1.4092	
	Hydrolyzes to form stable silanol solutions in neutral water [5630-83-1]	TSCA-L HMIS: 2-2-1-X	25g	100g		
	SIA0055.0 ACETOXYMETHYLTRIMETHOXSILANE, 95% C ₆ H ₁₄ O ₅ Si	194.26	190-1° Flashpoint: 56°C (133°F)	1.085	1.4031	
	[65625-39-0]	TSCA-L HMIS: 3-2-1-X	10g	50g		
	SIA0078.0 2-[(ACETOXY(POLYETHYLENEOXY)PROPYL)TRIETHOXSILANE, 95% 500 - 700			1.071	1.4527	
	Viscosity: 30 cSt	HMIS: 2-1-1-X	25g			
	SIA0100.0 3-ACETOXYPROPYLTRIMETHOXSILANE C ₈ H ₁₆ O ₅ Si	222.31	92° / 2 Flashpoint: 93°C (199°F)	1.062	1.4146	COMMERCIAL
	yc of treated surfaces: 37.5 mN/m Forms sol-gel derived epoxidation catalyst in combination with titanium diisopropoxide bispentanedionate. ¹ 1. Müller, C. et al. <i>Catal. Lett.</i> 2000 , <i>64</i> (1), 9-14, DOI 10.1023/A:109074617565 [59004-18-1]	EC 261-552-8 HMIS: 3-1-1-X	25g	100g	2kg	
	SIB0959.0 BENZOYLOXYPROPYLTRIMETHOXSILANE C ₁₃ H ₂₀ O ₅ Si	284.38	145° / 0.2	1.104	1.4806	
	[76241-02-6]	TSCA HMIS: 3-2-1-X	25g			
	SIC2067.0 10-(CARBOMETHOXY)DECYLDIMETHYLMETHOXSILANE C ₁₅ H ₃₂ O ₃ Si	288.50	130° / 0.3	0.903	1.4399	
	Long chain organofunctional silane [1211488-83-3]	HMIS: 2-1-1-X	10g	50g		











Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
 <p>SIC2072.0 2-(CARBOMETHOXY)ETHYLTRIMETHOXYSILANE METHYL (3-TRIMETHOXYSILYL)PROPIONATE) C₇H₁₆O₅Si</p>	208.29	75° / 1.5		1.069	1.410
Contains ~ 20% 1-(carbomethoxy)ethyltrimethoxysilane isomer [76301-00-3] HMIS: 3-3-1-X 10g Flashpoint: >43°C (>110°F)					
 <p>SIT8186.3 TRIETHOXYSILYLPROPOXY (POLYETHYLENEOXY) DODECANOATE, tech-95 C₂₇H₅₆O₉Si</p>	536.82			0.977	1.4479 ²⁵
2-4 EO Units Contact angle (treated surface), water: 61-2° Contact angle (treated surface), 2-ethylhexyl palmitate: <15° Provides embedded hydrophilicity with oleophilic compatibility Surface treatments stabilize particle dispersions. ¹ 1. Arkles, B. et al. in <i>Silanes and Other Coupling Agents</i> ; Mittal, K., Ed.; VSP (Brill), 2009, Vol. 5, p. 51. [1041420-54-5] TSCA-L HMIS: 2-1-1-X 25g 100g					
					

Halogen Functional Silanes

Halogen Functional Silanes - Trialkoxy






 <p>SIB1879.7 4-BROMOBUTYLTRIMETHOXYSILANE C₇H₁₇BrO₃Si</p>	257.20	68-9° / 1		1.259	1.4437
[226558-82-3] HMIS: 3-2-1-X 10g					
 <p>SIB1886.0 7-BROMOHEPTYLTRIMETHOXYSILANE C₁₀H₂₃BrO₃Si</p>	299.28	83-5° / 1		1.1695 ²⁵	1.4493 ²⁵
HMIS: 3-2-1-X 10g					
 <p>SIB1894.2 5-BROMOPENTYLTRIMETHOXYSILANE C₈H₁₉BrO₃Si</p>	271.22	110° / 1		1.222 ²⁵	1.446
[773893-02-0] HMIS: 3-2-1-X 10g					
 <p>SIB1906.0 3-BROMOPROPYLTRIMETHOXYSILANE C₆H₁₅BrO₃Si</p>	243.17	130° / 45		1.293	1.4400
Flashpoint: 82°C (180°F) Forms self-assembled monolayers which can be modified with pyridine ligands. ¹ 1. Paulson, S. et al. <i>J. Chem. Soc., Chem. Commun.</i> 1992 , 21, 1615. [51826-90-5] EC 257-453-4 HMIS: 2-2-1-X 10g 2.5kg					
 <p>SIB1909.0 11-BROMOUNDECYLTRIMETHOXYSILANE, 95% C₁₄H₃₁BrO₃Si</p>	355.39	158° / 0.8		1.119	1.4559
Contains undecyltrimethoxysilane [17947-99-8] TSCA-L HMIS: 2-1-0-X 10g					
 <p>SIC2280.4 3-CHLOROISOBUTYLTRIMETHOXYSILANE C₇H₁₅ClO₃Si</p>	212.75	59-60° / 2		1.054 ²⁵	1.4277 ²⁵
[17256-27-8] HMIS: 3-2-1-X 25g					
 <p>SIC2282.0 2-(CHLOROMETHYL)ALLYLTRIMETHOXYSILANE C₇H₁₅ClO₃Si</p>	210.73	128° / 70		1.09	
Flashpoint: 89°C (192°F) Versatile coupling agent [39197-94-9] HMIS: 3-2-1-X 2.5g					

	Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰	
	SIC2295.5 ((CHLOROMETHYL)PHENYLETHYL)TRIMETHOXYSILANE C ₁₂ H ₁₉ ClO ₃ Si Mixed <i>m</i> -, <i>p</i> -isomers Adhesion promoter for polyphenylenesulfide and polyimide coatings Employed as a high temperature coupling agent. ¹ 1. Arkles, B. et al. <i>Modern Plastics</i> 1980 , 57, 64. [68128-25-6]	274.82	115° / 1.5	Flashpoint: 130°C (266°F)	1.09 ²⁵	1.4930 ²⁵	COMMERCIAL
	SIC2295.7 CHLOROMETHYLPHENETHYLTRIS(TRIMETHYLSILOXY)SILANE, mixed <i>m</i> , <i>p</i> ; α , β isomers C ₁₈ H ₃₇ ClO ₃ Si ₄ Mixed <i>m</i> -, <i>p</i> -isomers [167426-89-3]	449.28			1.007	1.469	NEW
	SIC2296.2 (<i>p</i> -CHLOROMETHYL)PHENYLTRIMETHOXYSILANE, 95% C ₁₀ H ₁₅ ClO ₃ Si Coupling agent for polyimides, e.g in printed circuit boards Modifies silica for high-throughput peptide synthesis. ¹ 1. Houghten, R. et al. <i>J. Am. Chem. Soc.</i> 2005 , 127, 8582. [24413-04-5]	246.77	134-43° / 10	Flashpoint: 183°C (361°F)	1.14	1.4965	COMMERCIAL
							
	SIC2298.4 CHLOROMETHYLTRIETHOXYSILANE C ₇ H ₁₇ ClO ₃ Si Grignard reacts with chlorosilanes or intermolecularly to form carbosilanes. ¹ 1. Brondani, D. et al. <i>Tetrahedron Lett.</i> 1993 , 34, 2111. [15267-95-5]	212.75	90-1° / 25	Flashpoint: 47°C (117°F) TOXICITY: oral rat, LD50: 2,400 mg/kg	1.048	1.4069 ²⁵	COMMERCIAL
	SIC2298.5 CHLOROMETHYLTRIISOPROPOXYSILANE C ₁₀ H ₂₃ ClO ₃ Si Forms Grignard reagents [18162-82-8]	254.82	195-8°		0.9836	1.4145	NEW
	SIC2298.6 CHLOROMETHYLTRIMETHOXYSILANE C ₄ H ₁₁ ClO ₃ Si [5926-26-1]	170.67	156°	Flashpoint: 26°C (79°F) 10g	1.125	1.4070	
	SIC2407.0 3-CHLOROPROPYLTRIETHOXYSILANE C ₉ H ₂₁ ClO ₃ Si Adhesion promoter for polyamine and epoxy resins [5089-70-3]	240.80	100-2° / 10	Flashpoint: 78°C (172°F)	1.009	1.420	COMMERCIAL
	SIC2410.0 3-CHLOROPROPYLTRIMETHOXYSILANE C ₆ H ₁₅ ClO ₃ Si Viscosity, 20°: 0.56 cSt γ c of treated surface: 40.5 mN/m Specific wetting surface: 394 m ² /g Adhesion promoter for SBR hot-melt adhesives Powder flow control additive for dry powder fire extinguishing media [2530-87-2]	198.72	100° / 40	Flashpoint: 78°C (172°F) TOXICITY: oral rat, LD50: 5,628 mg/kg Vapor pressure, 100°: 40 mm	1.077 ²⁵	1.4183 ²⁵	COMMERCIAL
	SIC2428.0 11-CHLOROUNDACYLTRIETHOXYSILANE C ₁₇ H ₃₇ ClO ₃ Si [120876-31-5]	353.01	112° / 0.01	10g	0.944	1.4390	NEW


Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰	
 <p>SIC2429.0 11-CHLOROUNDACYLTRIMETHOXYSIANE C₁₄H₃₁ClO₃Si [17948-05-9]</p>	310.93	156-8° / 2 Flashpoint: >110°C (>230°F) 10g		0.9855	1.4226	NEW
 <p>SII6452.0 3-IODOPROPYLTRIMETHOXYSIANE C₆H₁₅IO₃Si [14867-28-8]</p> <p>Couples zeolite monolayers to glass.¹ 1. Ha, K. et al. <i>Adv. Mater.</i> 2002, <i>12</i>(15), 1114.</p>	290.17	79-80° / 2 Flashpoint: 78°C (172°F)		1.475	1.4714	
 <p>SIT8397.0 (3-TRIMETHOXYSIYL)PROPYL 2-BROMO-2-METHYLPROPIONATE C₁₀H₂₁BrO₅Si [314021-97-1]</p> <p>For surface initiated ATRP polymerization.^{1,2} 1. Mulvihill, M. et al. <i>J. Am. Chem. Soc.</i> 2005, <i>127</i>, 16040. 2. Huck, J. et al. <i>J. Mater. Chem.</i> 2004, <i>14</i>, 730.</p>	329.27	90-5° / 0.5		1.243 ²⁵		
 <p>SIV9064.0 VINYL(CHLOROMETHYL)DIMETHOXYSIANE C₅H₁₁ClO₂Si [1314981-48-0]</p> <p>Multi-functional coupling agent</p>	166.68					NEW
Halogen Functional Silanes - Dialkoxy						
 <p>SIC2292.0 CHLOROMETHYLMETHYLDIETHOXYSIANE C₆H₁₅ClO₂Si [2212-10-4]</p>	182.72	160-1° Flashpoint: 38°C (100°F) TOXICITY: oral rat, LD50: 1,300 mg/kg Vapor pressure, 70°: 20 mm 25g 100g		1.000 ²⁵	1.407	
 <p>SIC2295.2 ((CHLOROMETHYL)PHENYLETHYL)METHYLDIMETHOXYSIANE C₁₂H₁₉ClO₂Si [160676-60-8]/[160676-58-4]</p> <p>Mixed <i>m</i>-, <i>p</i>-isomers Intermediate for silicone analog of Merrifield resins</p>	258.82	120-5° / 0.5				
 <p>SIC2352.0 3-CHLOROPROPYLMETHYLDIETHOXYSIANE C₈H₁₉ClO₂Si [13501-76-3]</p> <p>Intermediate for functional silicone polymers</p>	210.77	81-3° / 8 Flashpoint: 80°C (176°F) 100g		0.9744	1.4260	NEW
 <p>SIC2353.0 3-CHLOROPROPYLMETHYLDIISOPROPOXYSIANE C₁₀H₂₃ClO₂Si [18171-19-2]</p>	238.84	92-4° / 10 25g		0.93		NEW
 <p>SIC2355.0 3-CHLOROPROPYLMETHYLDIMETHOXYSIANE C₆H₁₅ClO₂Si [18171-19-2]</p> <p>Specific wetting surface: 428 m²/g</p>	182.72	70-2° / 11 Flashpoint: 80°C (176°F) 100g		1.0250	1.4253	COMMERCIAL
 <p>SII6451.2 (3-IODOPROPYL)METHYLDIISOPROPOXYSIANE C₁₀H₂₃IO₂Si [18171-19-2]</p>	330.27	50-3° / 0.3 10g		1.257	1.4623	NEW

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
------	----	----------	---------	------------------------------	------------------------------

Halogen Functional Silanes - Monoalkoxy




	SIC2278.0 3-CHLOROISOBUTYLDIMETHYLMETHOXY-SILANE C ₇ H ₁₇ ClOSi [18244-08-1]	180.75 TSCA HMIS: 3-3-1-X	182° 25g	0.950	1.4331 ²⁵
	SIC2286.0 CHLOROMETHYLDIMETHYLETHOXY-SILANE C ₉ H ₁₃ ClOSi Dipole moment: 2.14 debye [13508-53-7]	152.70 TSCA EC 236-835-4 HMIS: 3-3-1-X	132-3° Flashpoint: 26°C (79°F) TOXICITY: oral rat, LD50: 1,550 mg/kg 25g	0.944 ²⁵	1.412 ²⁵
	SIC2295.2 ((CHLOROMETHYL)PHENYLETHYL)-METHYLDIMETHOXY-SILANE C ₁₂ H ₁₉ ClO ₂ Si Mixed <i>m</i> -, <i>p</i> -isomers Intermediate for silicone analog of Merrifield resins [160676-60-8]/[160676-58-4]	258.82 HMIS: 2-1-1-X	120-5° / 0.5 25g		
	SIC2337.0 3-CHLOROPROPYLDIMETHYLETHOXY-SILANE C ₇ H ₁₇ ClOSi [13508-63-9]	180.75 EC 236-837-5 HMIS: 2-3-1-X	87° / 30 Flashpoint: 46°C (115°F) 25g	0.932 ²⁵	1.427 ²⁵
	SIC2338.0 3-CHLOROPROPYLDIMETHYLMETHOXY-SILANE, 95% C ₆ H ₁₅ ClOSi [18171-14-7]	166.73 EC 242-055-5 HMIS: 3-2-1-X	170-1° Flashpoint: 39°C (102°F) 10g	0.941	1.4278

Halogen Functional Silanes - Dipodal

	SIC2279.3 1-(3-CHLOROISOBUTYL)-1,1,3,3,3-PENTAETHOXY-1,3-DISILAPROPANE, 95% C ₁₅ H ₃₅ ClO ₅ Si ₂ Pendant dipolar silane [160676-60-8]/[160676-58-4]	387.06 HMIS: 3-1-1-X	115-7° / 0.5 10g	1.020	
--	--	-------------------------	---------------------	-------	--

Hydroxyl Functional Silanes

Hydroxyl Functional Silanes - Trialkoxy

	SIH6172.0 N-(HYDROXYETHYL)-N-METHYLAMINOPROPYLTRIMETHOXY-SILANE, 75% in methanol C ₉ H ₂₃ NO ₄ Si [330457-46-0]	237.37 HMIS: 3-4-1-X	Flashpoint: 11°C (52°F) 25g	100g	0.99	1.417
	SIH6175.0 HYDROXYMETHYLTRIETHOXY-SILANE, 50% in ethanol TRIETHOXY-SILYL METHANOL C ₇ H ₁₈ O ₄ Si [162781-70-6]	194.31 TSCA-L HMIS: 2-4-0-X	Flashpoint: 15°C (59°F) 25g		0.866	
	SIT8189.0 N-(3-TRIETHOXY-SILYL-PROPYL)-GLUCONAMIDE, 50% in ethanol GLUCONAMIDOPROPYLTRIETHOXY-SILANE C ₁₅ H ₃₃ NO ₉ Si [104275-58-3]	399.51 HMIS: 2-4-1-X	Flashpoint: 15°C (59°F) 25g	100g	0.951	2kg

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
------	----	----------	---------	------------------------------	------------------------------

	SIT8189.5				
	N-(3-TRIETHOXSILYLPROPYL)-4-HYDROXYBUTYRAMIDE C ₁₃ H ₂₉ NO ₅ Si 307.47			1.02	1.4533
Anchoring reagent for light directed synthesis of DNA on glass. ¹ 1. McGall, G. et al. <i>J. Am. Chem. Soc.</i> 1997 , <i>119</i> , 5081.					
[156214-80-1]	HMIS: 2-2-1-X	10g	50g		

	SIT8192.0				
	N-(TRIETHOXSILYLPROPYL)-O-POLYETHYLENE OXIDE URETHANE, 95% C ₁₀ H ₂₂ NO ₄ SiO(CH ₂ CH ₂ O) ₄ -H 400-500			1.09	1.4540 ²⁵
Contains some bis(urethane) analog Viscosity: 75-125 cSt Hydrophilic surface modifier Forms PEGylated glass surfaces suitable for capillary electrophoresis. ¹ 1. Razunguzwa, T. et al. <i>Anal. Chem.</i> 2006 , <i>78</i> , 4326.					
[74695-91-3]	TSCA HMIS: 2-1-1-X	25g	100g	2kg	COMMERCIAL

Dipodal Hydroxyl Functional Silanes

	SIH6171.5				
	N-(HYDROXYETHYL)-N,N-BIS(TRIMETHOXSILYLPROPYL)AMINE, 65% in methanol C ₁₄ H ₃₅ NO ₇ Si ₂ 385.61		Flashpoint: 15°C (59°F)	0.97	NEW
Dipodal silane with hydroxyl functionality					
[264128-94-1]	TSCA HMIS: 3-4-1-X	10g			

Masked Hydroxyl Functional Silanes

	SIT8572.8				
	11-(TRIMETHYLSILOXY)UNDECYLTRIETHOXSILANE C ₂₀ H ₄₆ O ₄ Si ₂ 406.75	145° / 0.3		0.887 ²⁵	1.4264 ²⁵
Masked hydroxyl - deprotected after deposition with acidic aqueous ethanol					
[75389-03-6]	HMIS: 2-1-1-X	5g			

Isocyanate and Masked Isocyanate Functional Silanes

Isocyanate Functional Silanes - Trialkoxy








	SI16455.0				
	3-ISOCYANATOPROPYLTRIETHOXSILANE, 95% C ₁₀ H ₂₁ NO ₄ Si 247.37	130° / 20		0.990	1.4190
Flashpoint: 80°C (176°F) TOXICITY: oral rat, LD50: 710 mg/kg					
Component in hybrid organic/inorganic urethanes. ¹ 1. Cuney, S. et al. <i>Better Ceramics Through Chemistry VII (MRS. Symp. Proc.)</i> 1996 , <i>435</i> , 143.					
[24801-88-5]	TSCA EC 246-467-6 HMIS: 3-2-1-X	25g	100g	2kg	COMMERCIAL

	SI16456.0				
	3-ISOCYANATOPROPYLTRIMETHOXSILANE, 95% C ₇ H ₁₅ NO ₄ Si 205.29	95-8° / 10		1.073	1.4219
Flashpoint: 108°C (226°F) TOXICITY: oral rat, LD50: 878 mg/kg Autoignition temperature: 265°C					
[15396-00-6]	TSCA EC 239-415-9 HMIS: 3-2-1-X	25g	100g	2kg	COMMERCIAL

Isocyanate Functional Silanes - Dialkoxy

	SI16454.45				
	3-ISOCYANATOPROPYLMETHYLDIETHOXSILANE, 95% C ₉ H ₁₉ NO ₃ Si 217.34	110-5° / 10			NEW
Reacts rapidly with amine and hydroxyl functional species that can hydrolyze to form siloxane polymers					
[33491-28-0]	HMIS: 3-2-1-X	10g			

	SI16454.5				
	3-ISOCYANATOPROPYLMETHYLDIMETHOXSILANE, tech-95 C ₇ H ₁₅ NO ₃ Si 189.29	61° / 1		1.03	NEW
Contains isomers					
[26115-72-0]	HMIS: 3-2-1-X	10g			

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰	
 <p>SIT7907.7 (THIOCYANATOMETHYL)PHENETHYLTRIMETHOXYSIANE, tech-95 C₁₃H₁₉NO₃SSi 297.44 Contains isomers On exposure to UV light of 254 nm undergoes conversion to isothiocyanate which reacts w/ amines, etc. HMIS: 3-2-1-X</p>		100-5° / 0.5				NEW
 <p>SIT7908.0 3-THIOCYANATOPROPYLTRIETHOXYSIANE, 96% C₁₀H₂₁NO₃SSi 263.43 Complexing agent for Ag, Au, Pd, Pt.¹ Potential adhesion promoter for gold.² 1. Schilling, T. et al. <i>Mikrochemica Acta</i> 1996, 124, 235. 2. Ciszek, J. W. et al. <i>J. Am. Chem. Soc.</i> 2004, 126, 13172. [34708-08-2] TSCA EC 252-161-3 HMIS: 3-1-1-X</p>		95° / 0.1		1.03	1.4460	COMMERCIAL
 <p>SIT8186.5 N-(3-TRIETHOXYSIYLPROPYL)-O-t-BUTYL CARBAMATE C₁₄H₃₁NO₅Si 321.49 Masked isocyanate [137376-38-6] HMIS: 2-1-1-X</p>		110-5° / 0.2		0.990	1.4334	
 <p>SIT8188.0 TRIETHOXYSIYLPROPYL ETHYL CARBAMATE C₁₂H₂₇NO₅Si 293.44 Masked isocyanate [17945-05-0] TSCA EC 241-872-4 HMIS: 2-1-1-X</p>		124-6° / 0.5		1.015	1.4321	COMMERCIAL
 <p>SIT8407.0 N-TRIMETHOXYSIYLPROPYLMETHYL CARBAMATE METHYL [3-(TRIMETHOXYSIYL)PROPYL]CARBAMATE C₈H₁₉NO₅Si 237.32 Viscosity: 12 cSt [23432-62-4] HMIS: 3-2-1-X</p>		102° / 0.75		1.1087		NEW
 <p>SIT8411.0 2-(3-TRIMETHOXYSIYLPROPYLTHIO)THIOPHENE C₁₀H₁₈O₃S₂Si 278.46 Contact angle, water on treated silica surface: 76° [1364140-50-0] HMIS: 3-2-1-X</p>		125-7° / 0.4		1.136 ²⁵	1.5123 ²⁵	NEW
 <p>SIT8717.0 TRIS(3-TRIMETHOXYSIYLPROPYL)ISOCYANURATE, tech-95 C₂₁H₄₅N₃O₁₂Si₃ 615.86 Viscosity: 325-350 cSt. Coupling agent for polyimides to silicon metal Adhesion promoter for hotmelt adhesives Forms periodic mesoporous silicas.¹ 1. Zhang, W. et al. <i>Chem. Mater.</i> 2007, 19, 2663. [26115-70-8] TSCA EC 247-465-8 HMIS: 2-1-1-X</p>				1.170	1.4610	COMMERCIAL


Phosphine and Phosphate Functional Silanes

	<p>SIB1091.0 BIS(2-DIPHENYLPHOSPHINOETHYL)METHYLSILYLETHYLTRIETHOXYSIANE, mixed isomers C₃₇H₅₀O₃P₂Si₂ 660.92</p>	1.07 1.5746
<p>Analogous structures form ruthenium(II) complexes with high selectivity for hydrogenation and non-leachable binding to solid supports.¹ 1. Wu, D. et al. <i>Chem. Mater.</i> 2005, <i>17</i>, 3951.</p>		
HMIS: 2-2-1-X 1.0g		
	<p>SID3385.0 (2-DICYCLOHEXYLPHOSPHINOETHYL)TRIETHOXYSIANE C₂₀H₄₁O₃PSi 388.60 140° / 0.03</p>	0.979 ²⁵ 1.4811 ²⁵
<p>Ligand for immobilization of precious metal catalytic complexes [55289-47-9] HMIS: 3-1-0-X 5g</p>		
	<p>SID3411.0 (2-DIETHYLPHOSPHATOETHYL)METHYLDIETHOXYSIANE, tech-95 C₁₁H₂₇O₅PSi 298.39 124° / 2</p>	1.020 1.4270 ²⁵
<p>Comonomer for hydrophilic coatings [18048-06-1] HMIS: 3-2-1-X 10g</p>		
	<p>SID3412.0 (2-DIETHYLPHOSPHATOETHYL)TRIETHOXYSIANE, tech-95 DIETHYLPHOSPHONATOETHYLTRIETHOXYSIANE C₁₂H₂₉O₆PSi 328.41 141° / 2</p>	1.031 ²⁵ 1.4216
<p>Water-soluble silane; anti-pilling agent for textiles. Hydrolysis product catalytically hydrates olefins, forming alcohols.¹ Forms corrosion resistant films for magnesium alloys.² 1. Young, F. et al. U.S. Patent 3,816,550, 1974. 2. Kramov, A. et al. <i>Thin Solid Films</i> 2006, <i>174</i>, 514.</p>		
[757-44-8] TSCA EC 212-056-5 HMIS: 3-2-1-X 25g 100g		
	<p>SID4557.5 (2-DIPHENYLPHOSPHINO)ETHYLDIMETHYLETHOXYSIANE C₁₈H₂₅O₃PSi 316.46 160° / 1</p>	1.004 1.5630
[359859-29-3] HMIS: 2-2-1-X 10g		
	<p>SID4558.0 2-(DIPHENYLPHOSPHINO)ETHYLTRIETHOXYSIANE C₂₀H₂₉O₃PSi 376.50 182° / 1.3</p>	1.05 1.5384
<p>Immobilizing ligand for precious metals Adhesion promoter for gold substrates in microelectronic applications.¹ Forms stable bonds to silica and basic alumina suitable for catalyst immobilization.² Forms luminescent gels on hydrolysis with (EtO)₄Si and Eu(NO₃)₃.³ Used to immobilize an iridium catalyst for the enantioselective hydrogenation of aryl ketones.⁴ Used in the preparation of solid-phase Pd catalyst for Suzuki-Miyaura cross-coupling.⁵ 1. Helbert, J. U.S. Patent 4,497,890, 1985. 2. Merchle, C. H. et al. <i>Chem. Mater.</i> 2001, <i>13</i>, 3617. 3. Corriu, R. et al. <i>J. Chem. Soc., Chem. Commun.</i> 2001, 1116. 4. Liu, G. et al. <i>Adv. Synth. Catal.</i> 2008, <i>350</i>, 1464. 5. Zhang, X. et al. <i>Synthesis</i>, 2011, 2975.</p>		
[18586-39-5] TSCA EC 242-427-7 HMIS: 3-1-1-X 5g 25g		
	<p>SID4558.2 3-(DIPHENYLPHOSPHINO)PROPYLTRIETHOXYSIANE C₂₁H₃₁O₃PSi 390.53 190° / 1</p>	1.05 1.5384
[52090-23-0] HMIS: 3-1-1-X 1.0g		
	<p>SIT8378.5 3-(TRIHYDROXYSILYL)PROPYL METHYLPHOSPHONATE, MONOSODIUM SALT, 42% in water C₄H₁₂NaO₆PSi 238.18 1.25</p>	1.25
<p>Contains 4-5% methanol, sodium methylphosphonate Flashpoint: 79°C (174°F) Forms functionalized silica nanoparticles employed in amperometric glucose sensor.¹ 1. Zhao, W. et al. <i>Electrochim. Acta</i> 2013, <i>89</i>, 278.</p>		
[84962-98-1] TSCA EC 284-799-3 HMIS: 1-2-0-X 100g 500g 2.5kg		

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
------	----	----------	---------	------------------------------	------------------------------

Sulfur Functional Silanes

Sulfur Functional Silanes - Trialkoxy


SID3545.0

2,2-DIMETHOXY-1-THIA-2-SILACYCLOPENTANE

 $C_5H_{12}O_2SSi$

164.29

57-8° / 7


1.094

 Reagent for modification of silver and gold surfaces
 Coupling agent for rubber

[26903-85-5]

HMIS: 3-3-1-X

25g


SIM6475.0

3-MERCAPTOPROPYLTRIETHOXYLSILANE, 95%

 $C_9H_{22}O_3SSi$

238.42

210°

0.9325

1.4331

Flashpoint: 88°C (190°F)

TOXICITY: oral rat, LD50: >2,000 mg/kg

For blocked version see SIO6704.0

Used to make thiol-organosilica nanoparticles.¹1. Nakamura, M.; Ishimura, K. *Langmuir* **2008**, *24*, 5099.


[14814-09-6]

TSCA EC 238-883-1 HMIS: 2-2-1-X

25g

100g

2kg


SIM6476.0

3-MERCAPTOPROPYLTRIMETHOXYLSILANE

 $C_6H_{16}O_3SSi$

196.34

93° / 40

1.051²⁵1.4502²⁵

Viscosity: 2 cSt

yc of treated surfaces: 41 mN/m

Specific wetting surface: 348 m²/g

Flashpoint: 96°C (205°F)

TOXICITY: oral rat, LD50: 2,380 mg/kg

Primary irritation index: 0.19

Coupling agent for EPDM and mechanical rubber applications

Adhesion promoter for polysulfide adhesives

For enzyme immobilization.¹Treatment of mesoporous silica yields highly efficient heavy metal scavenger.²Couples fluorescent biological tags to semiconductor CdS nanoparticles.³Modified mesoporous silica supports Pd in coupling reactions.⁴Used to make thiol-organosilica nanoparticles.⁵Forms modified glass and silica surfaces suitable for SILAR fabrication of CdS thin films.⁶1. Stjernlöf, P. et al. *Tetrahedron Lett.* **1990**, *31*, 5773.2. Liu, J. et al. *Science* **1997**, *276*, 923.3. Bruchez, M. et al. *Science* **1998**, *281*, 2013.4. Crudden, C. et al. *J. Am. Chem. Soc.* **2005**, *127*, 10045.5. Nakamura, M.; Ishimura, K. *Langmuir* **2008**, *24*, 5099.6. Sun, H. et al. *J. Dispersion Sci. Technol.* **2005**, *26*, 719.

[4420-74-0]

TSCA EC 224-588-5 HMIS: 3-2-1-X


100g

2kg

18kg



Adhesion promoter for structural polysulfide glass sealants


SIM6476.1

3-MERCAPTOPROPYLTRIMETHOXYLSILANE, 99+%

 $C_6H_{16}O_3SSi$

196.34

93° / 40

1.051²⁵1.4502²⁵

Flashpoint: 96°C (205°F)


Low fluorescence grade for high-throughput screening

[4420-74-0]

TSCA EC 224-588-5 HMIS: 3-2-1-X

25g

* in fluoropolymer bottle


SIM6480.0

11-MERCAPTOUNDECYLTRIMETHOXYLSILANE, 95%

 $C_{14}H_{32}O_3SSi$

308.55

150° / 0.5


0.955

Stabilizes ionic liquid drop micro-reactors.¹1. Zhang, X. et al. *J. Nanotechnol.* **2012**, *3*, 33.

[877593-17-4]

HMIS: 3-2-1-X

2.5g


SIO6704.0

S-(OCTANOYL)MERCAPTOPROPYLTRIETHOXYLSILANE

 $C_{17}H_{36}O_4SSi$

364.62

Flashpoint: 176°C (349°F)

0.9686

1.4515

TOXICITY: oral rat, LD50: >2,000 mg/kg

Masked mercaptan - deblocked with alcohols

Latent coupling agent for butadiene rubber

[220727-26-4]


TSCA

HMIS: 2-1-1-X

25g

100g

18kg


SIP6926.2

3-(2-PYRIDYLETHYL)THIOPROPYLTRIMETHOXYLSILANE

 $C_{13}H_{23}NO_3SSi$

301.48

156-7° / 0.25


1.089

1.498

[29098-72-4]

HMIS: 3-2-1-X

10g


SIP6926.4

3-(4-PYRIDYLETHYL)THIOPROPYLTRIMETHOXYLSILANE, 95%

 $C_{13}H_{23}NO_3SSi$

301.48

160-2° / 0.2

1.09

1.5037










pKa: 4.8

Immobilizable ligand for immunoglobulin IgG separation using hydrophobic charge induction chromatography (HCIC)

[198567-47-4]

HMIS: 3-2-1-X

10g

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰	
 <p>SIT7908.0 3-THIOCYANATOPROPYLTRIETHOXSILANE, 96% C₁₀H₂₁NO₃SSi</p>	263.43	95° / 0.1		1.03	1.4460	COMMERCIAL
<p>Complexing agent for Ag, Au, Pd, Pt.¹ Potential adhesion promoter for gold.² 1. Schilling, T. et al. <i>Mikrochemica Acta</i> 1996, 124, 235. 2. Ciszek, J. W. et al. <i>J. Am. Chem. Soc.</i> 2004, 126, 13172. [34708-08-2] TSCA EC 252-161-3 HMIS: 3-1-1-X</p>		50g	250g	2kg		
 <p>SIT8411.0 2-(3-TRIMETHOXSILYLPROPYLTHIO)THIOPHENE C₁₀H₁₈O₃S₂Si</p>	278.46	125-7° / 0.4		1.136 ²⁵	1.5123 ²⁵	NEW
<p>Contact angle, water on treated silica surface: 76° [1364140-50-0] HMIS: 3-2-1-X</p>		10g				
Sulfur Functional Silanes - Dialkoxy						
 <p>SIM6473.0 (MERCAPTOMETHYL)METHYLDIETHOXSILANE, 95% C₆H₁₆O₂SSi</p>	180.34	60° / 10		0.975	1.4446	
[55161-63-2] HMIS: 3-2-1-X		10g				
 <p>SIM6474.0 3-MERCAPTOPROPYLMETHYLDIMETHOXSILANE, 96% C₆H₁₆O₂SSi</p>	180.34	96° / 30		1.000	1.4502	COMMERCIAL
<p>Intermediate for silicones in thiol-ene UV-cure systems Adhesion promoter for polysulfide sealants Used to make thiol-organosilica nanoparticles.¹ 1. Nakamura, M.; Ishimura, K. <i>Langmuir</i> 2008, 24, 5099. [31001-77-1] TSCA EC 250-426-8 HMIS: 3-2-1-X</p>		100g	2kg	18kg		
Sulfur Functional Silanes - Dipodal						
 <p>SIB1820.5 BIS[m-(2-TRIETHOXSILYLETHYL)TOLYL]POLYSULFIDE, tech-90 C₃₀H₅₀O₆(S₂₋₄)₂Si₂</p>	627-691			1.10	1.533	
<p>Dark, viscous liquid Coupling agent for SBR rubber [198087-81-9]/[85912-75-0]/[67873-85-2] TSCA HMIS: 2-2-1-X</p>		25g	2kg			
 <p>SIB1824.6 BIS[3-(TRIETHOXSILYL)PROPYL]DISULFIDE, 90% BIS(3-TRIETHOXSILYL)-4,5-DITHIOOCTANE C₁₈H₄₂O₆S₂Si₂</p>	474.82			1.025	1.457	COMMERCIAL
<p>Contains sulfide and tetrasulfide Dipodal coupling agent/vulcanizing agent for rubbers Intermediate for mesoporous silicas with acidic pores.¹ 1. Alauzun, J. et al. <i>J. Am. Chem. Soc.</i> 2006, 128, 8718. [56706-10-6] TSCA EC 260-350-7 HMIS: 2-2-1-X</p>		25g	100g	2kg		
 <p>SIB1825.0 BIS[3-(TRIETHOXSILYL)PROPYL]TETRASULFIDE, tech-95 TESPT C₁₈H₄₂O₆S₄Si₂</p>	538.94	250° dec		1.095	1.49	COMMERCIAL
<p>Contains distribution of S₂ - S₁₀ species; average 3.8 Viscosity: 11 cSt Adhesion promoter for precious metals Coupling agent/vulcanizing agent for "green" tires Adhesion promoter for PVD copper on parylene.¹ 1. Pimanpang, S. et al. <i>J. Vac. Sci. Technol. A</i> 2006, 24, 1884. [40372-72-3] TSCA EC 254-896-5 HMIS: 2-2-1-X</p>		100g	2kg	18kg		
						
 <p>SIB1827.0 N,N'-BIS[3-(TRIETHOXSILYL)PROPYL]THIOUREA, 90% C₁₈H₄₄N₂O₆SSi₂</p>	484.73			1.047	1.4696	
<p>Forms films on electrodes for determination of mercury.¹ 1. Guo, Y. et al. <i>J. Pharm. Biol. Anal.</i> 1999, 19 175. [69952-89-2] HMIS: 2-1-1-X</p>		25g				










Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
------	----	----------	---------	------------------------------	------------------------------

Vinyl and Olefin Functional Silanes

Vinyl and Olefin Functional Silanes - Trialkoxy








	<p>SIA0482.0 11-ALLYLOXYUNDECYLTRIMETHOXYSILANE C₁₇H₃₆O₃Si ω-olefin for functional self-assembled monolayers (SAMs) [1196453-35-6]</p>	<p>332.56 HMIS: 2-1-0-X</p>	<p>140° / 0.5 5g</p>	<p>0.914</p>	<p>1.4415</p>
	<p>SIA0489.0 m-ALLYLPHENYLPROPYLTRIETHOXYSILANE C₁₈H₃₀O₃Si Coupling agent for amine functional aromatic optical coatings HMIS: 2-2-1-X</p>	<p>322.52</p>	<p>5g</p>		<p>NEW</p>
	<p>SIA0525.0 ALLYLTRIETHOXYSILANE 3-(TRIETHOXYSILYL)-1-PROPENE C₉H₂₀O₃Si Dipole moment: 1.79 debye Extensive review on the use in silicon-based cross-coupling reactions.¹ 1. Denmark, S. E. et al. <i>Organic Reactions</i>, Vol. 75, Denmark, S. E. Ed., John Wiley and Sons, 233, 2011. [2550-04-1] TSCA EC 219-843-2 HMIS: 2-2-1-X</p>	<p>204.34</p>	<p>176° Flashpoint: 47°C (117°F) Vapor pressure, 100°: 50 mm</p>	<p>0.9030</p>	<p>1.4074</p>
	<p>SIA0540.0 ALLYLTRIMETHOXYSILANE C₆H₁₄O₃Si Adhesion promoter for vinyl-addition silicones Allylation of ketones, aldehydes and imines with dual activation of a Lewis Acid and fluoride ion.¹ Used in the regioselective generation of the thermodynamically more stable enol trimethoxysilyl ethers, which in turn are used in the asymmetric generation of quaternary carbon centers.² Converts arylselenyl bromides to arylallylselenides.³ Allylates aryl iodides.⁴ 1. Yamasaki, S. et al. <i>J. Am. Chem. Soc.</i> 2002, <i>124</i>, 6536. 2. Ichibakase, T. et al. <i>Tetrahedron Lett.</i> 2008, <i>49</i>, 4427. 3. Bhadra, S. et al. <i>J. Org. Chem.</i> 2010, <i>75</i>, 4864. 4. Mowery, M. E.; DeShong, P. <i>J. Org. Chem.</i> 1999, <i>64</i>, 1684. F&F: Vol 18, p 14; Vol 19, p 360; Vol 20, p 85; Vol 21, p 3, Vol 12, p 395 [2551-83-9] TSCA EC 219-855-8 HMIS: 3-2-1-X</p>	<p>162.26</p>	<p>146-8° Flashpoint: 46°C (115°F)</p>	<p>0.963²⁵</p>	<p>1.4036²⁵ COMMERCIAL</p>
	<p>SIB0988.0 [(5-BICYCLO[2.2.1]HEPT-2-ENYL)ETHYL]TRIMETHOXYSILANE, tech-95, endo/exo isomers C₁₂H₂₂O₃Si [68323-30-8] HMIS: 2-1-1-X</p>	<p>242.39</p>	<p>65° / 10 25g</p>	<p>1.02</p>	<p>1.458</p>
	<p>SIB0990.0 (5-BICYCLO[2.2.1]HEPT-2-ENYL)METHYLDICHLOROSILANE, 95% C₉H₁₂Cl₂Si [18245-94-8] TSCA EC 242-122-9 HMIS: 3-2-1-X</p>	<p>207.17</p>	<p>74-5° / 10 Flashpoint: 86°C (187°F) 10g</p>	<p>1.151</p>	<p>1.4938 NEW</p>
	<p>SIB0992.0 (5-BICYCLO[2.2.1]HEPT-2-ENYL)TRIETHOXYSILANE NORBORNENYLTRIETHOXYSILANE C₁₃H₂₄O₃Si Coupling agent for norbornadiene rubbers Component in low dielectric constant films Undergoes ring-opening metathetic polymerization (ROMP) with RuCl₂(P(C₆H₅)₃)₃.¹ 1. Finkelstein, E. <i>10th Int'l Organosilicon Symp. Proc.</i> 1993, P-120. [18401-43-9] TSCA EC 242-278-8 HMIS: 2-2-1-X</p>	<p>256.42</p>	<p>106-8° / 8 Flashpoint: 98°C (208°F)</p>	<p>0.960</p>	<p>1.4486</p>
	<p>SIB1928.0 3-BUTENYLTRIETHOXYSILANE, 95% C₁₀H₂₂O₃Si Mixed isomers (mainly 3-butenyl) [57813-67-9] HMIS: 2-2-1-X</p>	<p>218.37</p>	<p>64° / 6 (-80°) Flashpoint: 73°C (163°F) TOXICITY: oral rat, LD50: >5,000 mg/kg 25g</p>	<p>0.90</p>	









	Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
	SIC2282.0 2-(CHLOROMETHYL)ALLYLTRIMETHOXYSIANE C ₇ H ₁₅ ClO ₃ Si Versatile coupling agent [39197-94-9]	210.73	128° / 70 Flashpoint: 89°C (192°F)		1.09	NEW
	SIC2459.5 [2-(3-CYCLOHEXENYL)ETHYL]TRIETHOXYSIANE C ₁₄ H ₂₈ O ₃ Si Contains isomers [77756-79-7]	272.46	Flashpoint: 120°C (248°F) 10g	50g	0.948	1.444
	SIC2460.0 [2-(3-CYCLOHEXENYL)ETHYL]TRIMETHOXYSIANE C ₁₁ H ₂₂ O ₃ Si Contains isomers Orients liquid crystals in display devices. ¹ Coupling agent for aramid fiber reinforced epoxy. ² 1. Sharp, <i>Chem. Abstr.</i> 101,81758g; Jap. Patent JP 58122517, 1983. 2. Lechner, U. <i>Chem. Abstr.</i> 112, 218118x; Germ. Offen. DE 3820971, 1989. [67592-36-3]	230.38	109° / 6 Flashpoint: 80°C (176°F)	10g	1.02	1.4476
	SIC2464.1 3-CYCLOHEXENYLTRIMETHOXYSIANE C ₉ H ₁₈ O ₃ Si [21619-76-1]	202.32	78-9° / 6 Flashpoint: 60°C (140°F) 5g		1.039	NEW
	SIC2520.0 (3-CYCLOPENTADIENYLPROPYL)TRIETHOXYSIANE C ₁₄ H ₂₆ O ₃ Si Dimer; may be cracked to monomer at ~ 190° at 100mm Employed in silica-supported purification of fullerenes. ¹ 1. Nie, B. et al. <i>J. Org. Chem.</i> 1996 , 61, 1870. [102056-64-4]	270.44	115° / 0.5 Flashpoint: 100°C (212°F)	10g	0.99	1.4513
	SID4610.3 2-(DIVINYLMETHYLSILYL)ETHYLTRIETHOXYSIANE C ₁₃ H ₂₈ O ₃ Si ₂	288.54	79-81° / 0.15 5g		0.895	
	SID4618.0 DOCOSENYLTRIETHOXYSIANE, 95% C ₂₈ H ₅₈ O ₃ Si Contains isomers Forms self-assembled monolayers that can be modified to hydroxyls. ¹ 1. Penansky, J. et al. <i>Langmuir</i> 1995 , 11, 953. [330457-44-8]	470.88	187-195° / 0.05	1.0g		
	SIH5919.0 HEXADECAFLUORODODEC-11-EN-1-YLTRIMETHOXYSIANE C ₁₅ H ₁₆ F ₁₆ O ₃ Si Forms self-assembled monolayers; reagent for immobilization of DNA HMIS: 3-1-1-X	576.35	90° / 0.5	1.0g	1.4145 ²⁵	1.3526 ²⁵
	SIH6164.2 5-HEXENYLTRIETHOXYSIANE, 95% C ₁₂ H ₂₆ O ₃ Si Primarily α-olefin [52034-14-7]	246.43	97° / 1 Flashpoint: 86°C (187°F) 10g		0.883	1.4185
	SIH6164.3 5-HEXENYLTRIMETHOXYSIANE, 95% C ₉ H ₂₀ O ₃ Si Adhesion promoter for Pt-cure silicones [58751-56-7]	204.34	193-4°	10g	0.927	NEW

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
SIO6709.0 7-OCTENYLTRIMETHOXYSILANE, tech-95  C ₁₁ H ₂₄ O ₃ Si Contains 10-15% internal olefin isomers Coupling agent for "in situ" polymerization of acrylamide for capillary electrophoresis. ¹ Employed in stretched DNA fibers for FISH (fluorescent in situ hybridization) mapping. ² Surface treatment for FISH and replication mapping on DNA fibers. ³ 1. Cifuentes, A. et al. <i>J. Chromatogr., A</i> 1999 , 830(2), 423. 2. Labit, H. et al. <i>BioTechniques</i> 2008 , 45, 649. 3. Labit, H. et al. <i>Biotechniques Protocol Guide</i> 2010 (48) DOI 10.2144/000113255. [52217-57-9] TSCA HMIS: 3-1-1-X 5g 25g	232.39	48-9° / 0.1		0.94	1.4305
SIP6902.6 O-(PROPARGYL)-N-(TRIETHOXYSILYL)PROPYL CARBAMATE, 90%  C ₁₃ H ₂₅ NO ₅ Si Inhibited with MEHQ Surface derivatization reagent enabling "click" chemistry of nanoparticles. ¹ 1. Achatz, D. et al. <i>Sensors and Actuators B</i> 2010 , 150, 211. [870987-68-1] HMIS: 2-2-1-X 25g	303.43	110-20° / 0.2		0.990	1.4461 ²⁵
SIS6990.0 STYRYLETHYLTRIMETHOXYSILANE, tech-90  C ₁₃ H ₂₀ O ₃ Si Inhibited with t-butyl catechol Copolymerization parameter, e _Q : -0.880, 1.500 Mixed m-, p-isomers and α-, β-isomers Contains ethylphenethyltrimethoxysilane [119181-19-0]/[52783-38-7] TSCA-E HMIS: 2-1-1-X store <5°C 10g 50g	252.38	98° / 0.1		1.02	1.505
SIS6993.0 3-(N-STYRYLMETHYL-2-AMINOETHYLAMINO)PROPYLTRIMETHOXYSILANE, 40% in methanol  C ₁₇ H ₃₀ N ₂ O ₃ Si Inhibited with BHT Coupling agent for unsaturated polyesters [34937-00-3] TSCA EC 252-297-3 HMIS: 3-4-1-X store <5°C 25g 100g 2kg	338.52		Flashpoint: 11°C (52°F)	0.871	1.3900
SIU9049.0 10-UNDECENYLTRIMETHOXYSILANE  C ₁₄ H ₃₀ O ₃ Si When treated on glass provides a contact angle of 100° [872575-06-9] HMIS: 2-1-1-X 5g	274.48	102-5° / 1		0.908	1.4334
SIV9088.4 O-(VINILOXYBUTYL)-N-TRIETHOXYSILYLPROPYL CARBAMATE, tech-95  C ₁₆ H ₃₃ NO ₅ Si Inhibited with MEHQ UV reactive coupling agent [159856-61-8] HMIS: 3-2-1-X 10g	363.53			1.015	1.4454
SIV9098.0 VINYLTRIACTOXYLSILANE  C ₈ H ₁₂ O ₆ Si Crosslinker for moisture-cure silicone RTVs with greater liquid range for formulation, faster moisture-cure rate and better substrate adhesion properties than methyltriactoxysilane [4130-08-9] TSCA EC 223-943-1 HMIS: 3-2-1-X store <5°C 100g 2kg 18kg	232.26	112-3° / 1	(10-13°) Flashpoint: 88°C (190°F)	1.167	1.423
SIV9099.0 VINYLTRI-t-BUTOXYSILANE  C ₁₄ H ₃₀ O ₃ Si [5356-88-7] EC 226-343-8 HMIS: 2-2-1-X 10g 50g	274.47	54° / 2	Flashpoint: 79°C (174°F)	0.869	
SIV9112.0 VINYLTRIETHOXYSILANE  C ₈ H ₁₈ O ₃ Si ΔH _{form} : -463.5 kcal/mole ΔH _{vap} : 6.8 kcal/mole Dipole moment: 1.69 debye Specific wetting surface area: 412 m ² /g Copolymerization parameters- e _Q : -0.42, 0.028 yc of treated glass surface: 25 mN/m Relative hydrolysis rate versus SIV9220.0, vinyltrimethoxysilane; 0.05 Forms copolymers with ethylene Couples fillers or fiberglass to resins See VEE-005 for polymeric version Reacts with enamines to give (E)-β-silylenamines, which cross-couple w/ aryl iodides to give β-aryl enamines. ¹ Extensive review on the use in silicon-based cross-coupling reactions. ² 1. Marciniak, B. et al. <i>J. Org. Chem.</i> 2005 , 70, 8550. 2. Denmark, S. E. et al. <i>Organic Reactions</i> , Vol. 75, Denmark, S. E. ed., John Wiley and Sons, 233, 2011 . [78-08-0] TSCA EC 201-081-7 HMIS: 1-2-1-X 25g 2kg 16kg	190.31	160-1°	Flashpoint: 44°C (111°F) TOXICITY: oral rat, LD50: 22,500 mg/kg Autoignition temperature: 268°C Vapor pressure, 20°: 5 mm Specific heat: 0.25 cal/g/°	0.903	1.3960

COMMERCIAL

COMMERCIAL

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰	
 SIV9209.0 VINYLTRIISOPROPENOXYSIANE, tech-95 C ₁₁ H ₁₈ O ₃ Si Employed as a cross-linker and in vapor phase derivatization; byproduct is acetone [15332-99-7] TSCA EC 239-362-1 HMIS: 1-3-1-X	226.35	73-5° / 12		0.926	1.4373	COMMERCIAL
 SIV9210.0 VINYLTRIISOPROPOXYSIANE C ₁₁ H ₂₄ O ₃ Si Copolymerization parameters- e,Q: -0.36, 0.031 Used for free-radical cure water-borne resin systems Adhesion promotor for vinyl acetate/ethylene latex Relative hydrolysis rate versus vinyltrimethoxysilane; 0.0015 [18023-33-1] TSCA EC 241-931-4 HMIS: 1-2-1-X	232.39	179-81°		0.8659	1.3961 ²⁵	COMMERCIAL
 SIV9220.0 VINYLTRIMETHOXYSIANE C ₅ H ₁₂ O ₃ Si Viscosity: 0.6 cSt Copolymerization parameters- e,Q: -0.38, 0.031 Specific wetting surface area: 528 m ² /g Employed in two-stage ¹ and one-stage ² graft polymerization/crosslinking for PE Copolymerizes with ethylene to form moisture crosslinkable polymers. ³ Converts arylselenyl bromides to arylvinylselenides. ⁴ Reacts with anhydrides to transfer both vinyl and methoxy and thus form the mixed diester. ⁵ Cross-couples w/ α-bromo esters to give α-vinyl esters in high ee. ⁶ 1. Scott, H. U.S. Patent 3,646,155, 1972. 2. Swarbrick, P. et al. U.S. Patent 4,117,195, 1978. 3. Isaka, T. et al. U.S. Patent 4,413,066, 1983. 4. Bhadra, S. et al. <i>J. Org. Chem.</i> 2010 , <i>75</i> , 4864. 5. Luo, F. et al. <i>J. Org. Chem.</i> 2010 , <i>75</i> , 5379. 6. Strotman, N. A.; Sommer, S.; Fu, G. C. <i>Angew. Chem., Int. Ed. Engl.</i> 2007 , <i>46</i> , 3556. [2768-02-7] TSCA EC 220-449-8 HMIS: 3-3-1-X	148.23	123°		0.970	1.3930	COMMERCIAL
 Vinylsilanes are used in PE and EPDM insulated wire and cable						
 SIV9275.0 VINYLTRIS(2-METHOXYETHOXY)SIANE C ₁₁ H ₂₄ O ₆ Si Employed in peroxide graft-moisture crosslinking of polyethylene Relative hydrolysis rate versus SIV9220.0, vinyltrimethoxysilane; 0.50 Coupling agent for kaolin in EPDM/PE cable formulations. ¹ 1. Arkles, B. et al. <i>Modern Plastics</i> 1987 , <i>64</i> , 138. [1067-53-4] TSCA EC 213-934-0 HMIS: 3-1-1-X	280.39	284-6°		1.0336 ²⁵	1.4271 ²⁵	COMMERCIAL
 SIV9277.0 VINYLTRIS(1-METHOXY-2-PROPOXY)SIANE C ₁₄ H ₃₀ O ₆ Si [303746-21-6] HMIS: 2-1-1-X	322.47	Flashpoint: 122°C (252°F)		0.981	1.424	
 SIV9280.0 VINYLTRIS(METHYLETHYLKETOXIMINO)SIANE, tech-95 C ₁₄ H ₂₇ N ₃ O ₃ Si Neutral cross-linker/coupling agent for condensation cure silicones Byproduct: methylethylketoxime [2224-33-1] TSCA EC 218-747-8 HMIS: 3-3-1-X	313.47	113° / 0.1	(-22°)	0.982 ²⁵	1.465	COMMERCIAL

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰	
Vinyl and Olefin Functional Silanes - Dialkoxy						
	SIA0415.0 N-ALLYL-AZA-2,2-DIMETHOXY-SILACYCLOPENTANE C ₈ H ₁₇ NO ₂ Si Coupling agent for nanoparticles [618914-49-1]	187.31	52-4° / 3			NEW
	SIA0485.0 ALLYLMETHYLDIMETHOXY-SILANE, tech-95 C ₆ H ₁₄ O ₂ Si Contains β-methylvinyl isomer [67681-66-7]	146.26	127-8°	0.871 ²⁵	1.4055 ²⁵	NEW
	SIB0990.3 (5-BICYCLO[2.2.1]HEPT-2-ENYL)METHYLDIETHOXY-SILANE C ₁₂ H ₂₂ O ₂ Si Comonomer for ROMP reactive resins [113276-73-6]	226.39	65-7° / 1	0.959	1.4525 ²⁵	NEW
	SIV9085.0 VINYL METHYLDIETHOXY-SILANE C ₇ H ₁₆ O ₂ Si Dipole moment: 1.27 debye Copolymerization parameters- e,Q; -0.86, 0.020 Chain extender, crosslinker for silicone RTVs and hydroxy-functional resins [5507-44-8]	160.29	133-4° Flashpoint: 16°C (61°F)	0.858	1.3998	COMMERCIAL
	SIV9086.0 VINYL METHYLDIMETHOXY-SILANE C ₈ H ₁₂ O ₂ Si Viscosity: 0.7 cSt Additive to moisture-cure silane modified polyurethanes as a water scavenger to prevent premature cure [16753-62-1]	132.23	103° Flashpoint: 15°C (59°F) Autoignition temperature: 245°C Vapor pressure, 20°: 38 mm	0.889	1.395	COMMERCIAL
Vinyl and Olefin Functional Silanes - Monoalkoxy						
	SIB0981.0 (5-BICYCLO[2.2.1]HEPT-2-ENYL)DIMETHYLETHOXY-SILANE C ₁₁ H ₂₀ O ₂ Si Undergoes ROMP reactions [899823-76-8]	196.36	67° / 3.5	0.926 ²⁵	1.4604 ²⁵	NEW
	SIT8732.0 TRIVINYLMETHOXY-SILANE, 95% C ₇ H ₁₂ O ₂ Si [193828-96-5]	140.25	131-3° Flashpoint: 26°C (79°F)		1.4400	
	SIV9072.0 VINYL DIMETHYLETHOXY-SILANE C ₆ H ₁₄ O ₂ Si Dipole moment: 1.23 debye Vinylates aryl halides. ¹ 1. Denmark, S. E.; Butler, C. R. <i>J. Am. Chem. Soc.</i> 2008 , <i>130</i> , 3690. [5356-83-2]	130.26	99-100° Flashpoint: 4°C (39°F)	0.790	1.3983	

Vinyl and Olefin Functional Silanes - Dipodal

	SIB1618.0 1,2-BIS(METHYLDIETHOXYSILYL)ETHYLENE C ₁₂ H ₂₈ O ₄ Si ₂ Mixed cis/trans isomers - primarily trans [124279-15-8]	292.52	80° / 2	0.918	1.4214	NEW
	HMIS: 2-3-1-X	25g				
	SIB1818.0 BIS(TRIETHOXYSILYLETHYL)VINYLMETHYLSILANE C ₁₉ H ₄₄ O ₆ Si ₃	452.82	141° / 0.15	0.943		
	HMIS: 2-1-1-X	5g				
	SIB1820.0 1,2-BIS(TRIETHOXYSILYL)ETHYLENE, 95% 4,4,7,7-TETRAETHOXY-3,8-DIOXA-4,7-DISILADEC-5-ENE C ₁₄ H ₃₂ O ₆ Si ₂ ~80% trans isomer Forms ethylene-bridged mesoporous silicas. ¹ 1. Vercaemst, C. et al. <i>Chem. Mater.</i> 2009 , <i>21</i> , 5792. [87061-56-1]	352.57	122-5° / 4	0.958	1.4168	
	HMIS: 2-2-1-X	5g		25g		
	SIB1824.9 1,3-BIS(3-TRIETHOXYSILYLPROPYL)POLYETHYLENOXY-2-METHYLENEPROPANE C ₅₀ H ₁₀₄ O ₂₀ Si ₂ (av) 1113.5 Vinyl functional hydrophilic dipodal coupling agent for protein immobilization HMIS: 2-2-1-X					1.0g
	SIB1832.5 1,1-BIS(TRIMETHOXYSILYLMETHYL)ETHYLENE, tech-95 C ₁₀ H ₂₄ O ₆ Si ₂ Contains ~ 10% cyclic siloxane condensation products Dipodal coupling agent [143727-20-2]	296.47	230-3°	1.05	1.429 ²⁵	
	HMIS: 3-3-1-X	5g				
	SIB1834.5 BIS(3-TRIMETHOXYSILYLPROPYL) FUMARATE, 96% C ₁₆ H ₃₂ O ₁₀ Si ₂ Dipodal silane Adhesion promoter for Pt-cure silicone RTVs [3371-62-8]	440.59	Flashpoint: >110°C (>230°F)	1.118	1.443	NEW
	TSCA	25g				





Vinyl and Olefin Functional Silanes - Polymeric

	SIV9112.2 VINYLTRIETHOXYSILANE, oligomeric hydrolysate			1.02		NEW
	Viscosity: 4-7 cSt [29434-25-1]	TSCA	HMIS: 2-2-1-X	100g	1kg	
	SIV9112.3 VINYLTRIETHOXYSILANE-PROPYLTRIETHOXYSILANE, oligomeric co-hydrolysate	600-1,000		1.02		NEW
	Viscosity: 3-7 cSt 9-11 wgt% vinyl [201615-10-3]	TSCA	HMIS: 2-1-1-X	100g	1kg	
	SIV9220.2 VINYLTRIMETHOXYSILANE, oligomeric hydrolysate			1.10	1.428	NEW
	Viscosity: 8-10 cSt Graft-coupling agent for grafted polyethylene composites. ¹ 1. Arkles, B. et al, <i>Modern Plastics</i> 1987 , <i>64</i> , 38. [131298-48-1]	TSCA	HMIS: 2-2-1-X	100g	2kg	

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
------	----	----------	---------	------------------------------	------------------------------

Multi-Functional and Polymeric Silanes

Polybutadiene

	SSP-055 TRIETHOXSILYL MODIFIED POLY-1,2-BUTADIENE, 50% in toluene	3,500-4,500		Flashpoint: 4°C (39°F)	0.90
	[72905-90-9]		TSCA	HMIS: 2-4-1-X store <5°C 100g	2kg
	SSP-056 TRIETHOXSILYL MODIFIED POLY-1,2-BUTADIENE, 50% in volatile silicone	3,500-4,500			0.93
	[72905-90-9]		TSCA	HMIS: 2-3-1-X store <5°C 100g	2kg
	SSP-058 DIETHOXYMETHYLSILYL MODIFIED POLY-1,2-BUTADIENE, 50% in toluene	3,500-4,500		Flashpoint: 4°C (39°F)	0.90
				HMIS: 2-4-1-X store <5°C 100g	3kg
	SSP-255 (30-35% TRIETHOXSILYLETHYL)ETHYLENE-(35-40% 1,4-BUTADIENE)-(25-30% STYRENE) terpolymer, 50% in toluene	4,500-5,500		Flashpoint: 4°C (39°F)	
				HMIS: 2-3-1-X 100g	



UV Active and Fluorescent Silanes

SIB1824.8

1,7-BIS(4-TRIETHOXSILYLPROPOXY-3-METHOXYPHENYL)-1,6-HEPTADIENE-3,5-DIONE, tech-90
 $C_{39}H_{60}O_{12}Si_2$ 777.07

UV max: 220, 232(vs), 354(broad)

Metal chelating chromophore



HMIS: 2-1-1-X 0.5g

[947329-82-0]

SIC2058.2

3-CARBAZOLYLPROPYLTRIETHOXSILANE

$C_{21}H_{29}NO_3Si$

371.55 185-195° / 0.3

1.072 1.5527²⁵


For non-linear optic materials

Employed in OLED fabrication.¹

1. DeMais, T. et al. *SPIE Proc.* **1998**, 3476, 338

[221105-38-0]

HMIS: 2-2-1-X 2.5g



SID4352.0

3-(2,4-DINITROPHENYLAMINO)PROPYLTRIETHOXSILANE, 95%

N-[3-(TRIETHOXSILYL)PROPYL]-2,4-DINITROPHENYLAMINE

$C_{15}H_{25}N_3O_7Si$

387.46

(27-30°)



1.5665

Viscous liquid or solid

UV max: 222, 258, 350(s), 410

Flashpoint: >110°C (>230°F)

Forms χ^2 non-linear optical sol-gel materials by corona poling.^{1,2}




1. Toussaere, E. et al. *Non-Linear Optics* **1992**, 2, 37.

2. Lebeau, J. et al. *J. Mater. Chem.* **1994**, 4, 1855.


[71783-41-0]

HMIS: 2-1-0-X 25g


Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
------	----	----------	---------	------------------------------	------------------------------




SIH6198.0
 2-HYDROXY-4-(3-METHYLDIETHOXSILYLPROPOXY)DIPHENYLKETONE, tech-90
 C₂₁H₂₈O₅Si 388.54
 Amber liquid
 Viscosity: 100-125 cSt
 HMIS: 2-1-1-X 25g
 D₄²⁰ 1.116²⁵ n_D²⁰ 1.5601²⁵




SIH6200.0
 2-HYDROXY-4-(3-TRIETHOXSILYLPROPOXY)DIPHENYLKETONE, tech-90
 C₂₂H₃₀O₆Si 418.56
 Amber liquid
 Viscosity, 25°: 125-150 cSt
 UV max: 230, 248, 296(s), 336
 Strong UV blocking agent for optically clear coatings, absorbs from 210-420 nm
 UV blocking agent.¹




1. Anthony, B. U.S. Patent 4,495,360, 1985.
 [79876-59-8] TSCA HMIS: 2-1-1-X 25g 100g 2kg




SIM6502.0
 O-4-METHYLCOUMARINYL-N-[3-(TRIETHOXSILYL)PROPYL]CARBAMATE
 C₂₀H₂₉NO₇Si 423.54 (88-90°)
 UV max: 223, 281, 319.5 (vs)
 Soluble: THF
 Immobilizeable fluorescent compound.¹




1. Arkles, B. U.S. Patent 4,918,200, 1990.
 [129119-78-4] HMIS: 2-2-1-X 10g




SIN6597.25
 NITROVERATRYLOXYCARBONYLAMIDOPROPYLTRIETHOXSILANE, 10% in tetrahydrofuran
 N-TRIETHOXSILYLPROPYL-O-4,5-DIMETHOXY-2-NITROBENZYL CARBAMATE
 C₁₉H₃₂N₂O₉Si 460.56
 UV max: 365 nm Flashpoint: -14°C (7°F)
 Photosensitive silane for lithography.¹
 1. del Campo, A. et al. *Angew. Chem.* **2005**, *44*, 4707
 [188541-09-5] HMIS: 3-4-1-X 1.0g



SIT8186.2
 7-TRIETHOXSILYLPROPOXY-5-HYDROXYFLAVONE, 50% in xylene
 C₂₄H₃₀O₇Si 458.58
 Contains non-reactive dyestuffs Flashpoint: 30°C (86°F)
 UV max: 350 nm
 [945761-08-0] HMIS: 2-3-1-X 1.0g 5g



SIT8187.0
 N-(TRIETHOXSILYLPROPYL)DANSYLAMIDE
 5-DIMETHYLAMINO-N-(3-TRIETHOXSILYLPROPYL)NAPHTHALENE-1-SULFONAMIDE
 C₂₁H₃₄N₂O₅SSi 454.66 115-9° / 0.1 1.12 1.5421
 Viscous liquid
 UV max: 222(s), 256, 354
 Fluorescent - employed as a tracer in UV cure composites
 Soluble in toluene, tetrahydrofuran
 Fluorescence probe for crosslinking in silicones.¹
 Employed in a chemically modified logic gate.²



1. Leezenberg, P. et al. *Chem. Mater.* **1995**, *7*, 1784.
 2. Mu, L. et al. *Angew. Chem., Int. Ed. Engl.* **2009**, *48*, 3469.
 [70880-05-6] TSCA EC 274-980-5 HMIS: 2-1-1-X 1.0g

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
------	----	----------	---------	------------------------------	------------------------------

SIT8191.0

3-(TRIETHOXY-SILYL-PROPYL)-p-NITROBENZAMIDE



C₁₆H₂₆N₂O₆Si

370.48

(54-5°)

1.5127

UV max: 224, 260, 292(s)

Used to prepare diazotizable supports for enzyme immobilization.¹

1. Weetall, H. U.S. Patent 3,652,761, 1972.

[60871-86-5]

TSCA EC 262-508-0 HMIS: 2-1-1-X

25g

SIT8192.4

(R)-N-TRIETHOXY-SILYL-PROPYL-O-QUININEURETHANE, 95%

C₃₀H₄₅N₃O₆Si



571.79

(82-4°)

UV max: 236(s), 274, 324, 334

Fluorescent, optically active silane

Soluble: warm toluene



[200946-85-6]

HMIS: 2-1-1-X

5g

Chiral Silanes**SIP6731.5**

(R)-N-1-PHENYLETHYL-N'-TRIETHOXY-SILYL-PROPYLUREA

C₁₈H₃₂N₂O₄Si

368.55

1.05²⁵

Viscous liquid or solid


Flashpoint: >110°C (>230°F)

Optically active silane; treated surfaces resolve enantiomers

[131206-15-0]

HMIS: 2-1-0-X

25g

**SIP6731.6**

(S)-N-1-PHENYLETHYL-N'-TRIETHOXY-SILYL-PROPYLUREA

C₁₈H₃₂N₂O₄Si

368.55

1.05²⁵

Viscous liquid or solid


Flashpoint: >110°C (>230°F)

Optically active silane; treated surfaces resolve enantiomers

[67240-22-2]

HMIS: 2-1-0-X

25g

**SIT8190.0**

(S)-N-TRIETHOXY-SILYL-PROPYL-O-MENTHOCARBAMATE

C₂₀H₄₁NO₅Si

406.63

Flashpoint: >110°C (>230°F)

0.985²⁵


1.4526

Optically active

[68479-61-8]

TSCA EC 270-863-8 HMIS: 2-1-1-X

10g

**SIT8192.4**

(R)-N-TRIETHOXY-SILYL-PROPYL-O-QUININEURETHANE, 95%

C₃₀H₄₅N₃O₆Si



571.79

(82-4°)

UV max: 236(s), 274, 324, 334

Fluorescent, optically active silane

Soluble: warm toluene








[200946-85-6]

HMIS: 2-1-1-X

5g

Biomolecular Probes







	<p>SIA0120.2 N-(ACETYLGLYCYL)-3-AMINOPROPYLTRIMETHOXSILANE, 5% in methanol C₁₀H₂₂N₂O₅Si 278.38 (171-3°) 0.80 Flashpoint: 11°C (52°F) Amino acid-tipped silane [1274903-53-5] HMIS: 3-4-1-X 25g</p>	NEW
	<p>SIA0126.0 3-(N-ACETYL-4-HYDROXYPROLYLOXY)PROPYLTRIETHOXSILANE, 25% in ethanol C₁₆H₃₁NO₇Si 377.51 0.872 Flashpoint: 15°C (59°F) Amino acid-tipped silane Hydrophilic reagent for biomimetic surface modification [1300591-79-0] HMIS: 2-3-0-X 5g</p>	NEW
	<p>SIA0127.0 N-(N-ACETYLLEUCYL)-3-AMINOPROPYLTRIETHOXSILANE, 12-15% in ethanol C₁₇H₃₆N₂O₅Si 376.58 0.816 Flashpoint: 15°C (59°F) Hydrophobic amino acid-tipped silane [1367348-25-1] HMIS: 2-3-1-X 2.5g</p>	NEW
	<p>SIT7909.7 3-(3-THYMINYLPROPIONOXY)PROPYLTRIMETHOXSILANE C₁₄H₂₄N₂O₇Si 360.44 Derivatized surfaces bind adenine modified polymers.¹ 1. Viswanathan, K. et al. <i>Polymer Preprints</i> 2005, 4602, 1133. [879197-67-8] HMIS: 2-2-1-X 1.0g</p>	
	<p>SIT8012.0 O-DL-a-TOCOPHEROLYLPROPYLTRIETHOXSILANE, tech-90 C₃₈H₇₀O₅Si 635.04 0.956 1.485 HMIS: 2-2-1-X 10g</p>	

Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
------	----	----------	---------	------------------------------	------------------------------



Trihydrosilanes

Silyl Hydrides are a distinct class of silanes that behave and react very differently than conventional silane coupling agents. They react with the liberation of byproduct hydrogen. Silyl hydrides can react with hydroxylic surfaces under both non-catalyzed and catalyzed conditions by a dehydrogenative coupling mechanism^{1,2}. Trihydrosilanes react with a variety of pure metal surfaces including gold, titanium, zirconium and amorphous silicon, by a dissociative adsorption mechanism.³ The reactions generally take place at room temperature and can be conducted in the vapor phase or with the pure silane or solutions of the silane in aprotic solvents. Deposition should not be conducted in water, alcohol or protic solvents.


1. Fadeev, A. et al. *J. Am. Chem. Soc.* **1999**, *121*, 12184.
2. N. Morita, N. et al. *J. Am. Chem. Soc.* 2014, *136*, 11370.
3. B. Arkles, B. et al. *J. Adhesion Sci. Technol.* **2012**, *26*, 141.

	SIB1907.8 11-BROMOUNDECYLSILANE C ₁₁ H ₂₅ BrSi	265.31	100-2° / 0.5	1.02	NEW
	Forms SAMs on titanium, gold and silicon surfaces. ¹ 1. Arkles, B. et al. <i>J. Adhes. Sci. Technol.</i> 2012 , <i>26</i> , 41. [469904-33-4]		HMS: 2-3-1-X	5g	
	SIC2268.5 2-CHLOROETHYLSILANE C ₂ H ₇ ClSi	94.61	69-71°	0.904	1.4232
	Intermediate for H ₃ SiCl; employed in CVD of SiN. ¹ 1. Arkles, B. et al. U.S. Patent 5,968,611, 1999. [18165-19-0]		HMS: 3-4-1-X store <5°C	10g	
	SID4629.6 DODECYLSILANE C ₁₂ H ₂₆ Si	200.44	80° / 7	0.7753	1.438 ²⁵
	Forms SAMs on gold and titanium surfaces [872-19-5]		HMS: 2-2-1-X	10g	
	SIO6635.0 n-OCTADECYLSILANE C ₁₈ H ₄₀ Si	284.60	195° / 15 (29°) Flashpoint: >110°C (>230°F)	0.794	
	Contains 4-6% C ₁₈ isomers Forms self-assembled monolayers on titanium. ¹ Reacts onto a gold surface to form monolayers of long alkyl chains. ² Forms SAMs on titanium, gold and silicon surfaces. ³ 1. Fadeev, A. et al. <i>J. Am. Chem. Soc.</i> 1989 , <i>121</i> , 12184. 2. Owens, T. M. et al. <i>J. Am. Chem. Soc.</i> 2002 , <i>124</i> , 6800. 3. Arkles, B. et al. <i>J. Adhes. Sci. Technol.</i> 2012 , <i>26</i> , 41. [18623-11-5]		TSCA EC 242-453-9	HMS: 2-1-1-X	25g
	SIT8173.0 (TRIDECAFLUORO-1,1,2-TETRAHYDROOCTYL)SILANE C ₈ H ₇ F ₁₃ Si	378.21	75° / 25	1.446	1.3184
	Provides vapor-phase hydrophobic surfaces on titanium, gold, silicon [469904-32-3]		HMS: 3-3-1-X	10g	
	SIU9048.0 10-UNDECENYLSILANE C ₁₁ H ₂₄ Si	184.40		0.768 ²⁵	1.4415 ²⁵
	Forms self-assembled monolayers on gold		HMS: 2-3-1-X	2.5g	

Dipodal Silyl Hydrides


	SIB1770.0 1,2-BIS(TETRAMETHYLDISILOXY)ETHANE, 95% C ₁₀ H ₃₀ O ₂ Si ₄	294.69	90-5° / 20	0.845	1.41
	[229621-70-9]		HMS: 2-2-1-X	25g	
	SID4593.5 1,10-DISILADECANE C ₈ H ₂₂ Si ₂	174.44	35° / 0.3	0.772	1.4461
	[4364-10-7]		HMS: 2-3-1-X	10g	

Dipodal Silanes - Non-Functional




SIB1660.0 BIS[(3-METHYLDIMETHOXYSILYL)PROPYL]POLYPROPYLENE OXIDE 600-800				1.00	1.452 ²⁵
Viscosity: 6,000-10,000 cSt. Hydrophilic dipodal silane With tin catalyst forms moisture-cross-linkable resins Flashpoint: >110°C (>230°F)					
[75009-88-0]	TSCA	HMIS: 3-1-1-X	100g	2kg	18kg

COMMERCIAL




SIB1817.0 1,2-BIS(TRIETHOXYSILYL)ETHANE HEXAETHOXYDISILETHYLENE, BSE C ₁₄ H ₃₄ O ₆ Si ₂ ΔHvap: 101.5 kJ/mole	354.59	96° / 0.3	(-33°)	0.957	1.4052	
Additive to silane coupling agents formulations that enhances hydrolytic stability Employed in corrosion resistant coatings/primers for steel and aluminum. ^{1,2} Sol-gels of α,ω-bis(trimethoxysilyl)alkanes reported. ³ Component in evaporation-induced self-assembly of mesoporous structures. ⁴ Forms mesoporous, derivatizable molecular sieves. ^{5,6} Hydrolysis kinetics studied. ⁷ 1. Van Ooij, W. et al. <i>J. Adhes. Sci. Technol.</i> 1997 , <i>11</i> , 29. 2. Van Ooij, W. et al. <i>Chemtech</i> 1999 , <i>28</i> , 3302. 3. Loy, D. A. et al. <i>J. Am. Chem. Soc.</i> 1999 , <i>121</i> , 5413. 4. Lu, Y. et al. <i>J. Am. Chem. Soc.</i> 2000 , <i>122</i> , 5258. 5. Molde, B. et al. <i>Chem. Mater.</i> 1999 , <i>11</i> , 3302. 6. Cho, E. et al. <i>Chem. Mater.</i> 2004 , <i>16</i> , 270. 7. Diaz-Benito, B. <i>Colloids and Surfaces A: Physicochemical Aspects</i> 2010 , <i>369</i> , 53.						
[16068-37-4]	TSCA	EC 240-212-2	HMIS: 3-1-1-X	25g	100g	2kg


COMMERCIAL




SIB1821.0 BIS(TRIETHOXYSILYL)METHANE 4,4,6,6-TETRAETHOXY-3,7-DIOXA-4,6-DISILANONANE C ₁₃ H ₃₂ O ₆ Si ₂	340.56	114-5° / 3.5		0.9741	1.4098	
Intermediate for sol-gel coatings, hybrid inorganic-organic polymers Forms methylene-bridged mesoporous structures. ¹ Forms modified silica membranes that separate propylene/propane mixtures. ² 1. Zhang, W. et al. <i>Chem. Mater.</i> 2005 , <i>17</i> , 6407. 2. Kanezashi, M. et al. <i>J. Membr. Sci.</i> 2012 , <i>415-416</i> , 478.						
[18418-72-9]	TSCA-L		HMIS: 3-2-1-X	5g	25g	100g



SIB1824.0 1,8-BIS(TRIETHOXYSILYL)OCTANE C ₂₀ H ₄₆ O ₆ Si ₂	438.76	172-5° / 0.75		0.926	1.4240	
Employed in sol-gel synthesis of mesoporous structures Crosslinker for moisture-cure silicone RTVs with improved environmental resistance Sol-gels of α,ω-bis(trimethoxysilyl)alkanes reported. ¹ 1. Loy, D.A. et al. <i>J. Am. Chem. Soc.</i> 1999 , <i>121</i> , 5413.						
[52217-60-4]	TSCA		HMIS: 2-1-1-X	25g	100g	2kg




SIB1829.0 1,2-BIS(TRIMETHOXYSILYL)DECANE C ₁₆ H ₃₈ O ₆ Si ₂	382.65	130-2° / 0.4		0.984	1.4303	
Pendant dipodal silane; employed in high pH HPLC Employed in the fabrication of luminescent molecular thermometers. ¹ 1. Brites, C. et al. <i>New J. Chem.</i> 2011 , <i>35</i> , 1173.						
[832079-33-1]	TSCA-L		HMIS: 3-2-1-X	25g	100g	





SIB1830.0 1,2-BIS(TRIMETHOXYSILYL)ETHANE C ₈ H ₂₂ O ₆ Si ₂ CAUTION: INHALATION HAZARD AIR TRANSPORT FORBIDDEN	270.43	103-4° / 5		1.068	1.4091	
Employed in fabrication of multilayer printed circuit boards. ¹ 1. Palladino, J. U.S. Patent 5,073,456, 1991.						
[18406-41-2]	TSCA	EC 242-285-6	HMIS: 4-2-1-X	25g	100g	2kg

Flashpoint: 65°C (149°F)
TOXICITY: ihl rat, LC50: 2.4 ppm
Vapor pressure, 20°: 0.08 mm




SIB1831.0 BIS(TRIMETHOXYSILYLETHYL)BENZENE, tech-95 C ₁₆ H ₃₀ O ₆ Si ₂ Mixed isomers	374.58	148-50° / 0.1		1.08	1.4734	
Forms high refractive index coatings Forms resins that absorb organics from aqueous media. ¹ 1. Edmiston, P. et al. <i>Sep. Purif. Technol.</i> 2009 , <i>66</i> , 532.						
[266317-71-9]	TSCA		HMIS: 2-1-0-X	10g	50g	2kg


Name	Mw	bp °C/mm	(mp °C)	D ₄ ²⁰	n _D ²⁰
 <p>SIB1832.0 1,6-BIS(TRIMETHOXSILYL)HEXANE C₁₂H₃₀O₆Si₂</p> <p>Sol-Gels of α,ω-bis(trimethoxysilyl)alkanes reported.¹ 1. Loy, D.A. et al. <i>J. Am. Chem. Soc.</i> 1999, 121, 5413. [87135-01-1]</p>	326.54	161° / 2		1.014	1.4213
			Flashpoint: 95°C (203°F)		
		HMS: 3-2-1-X	10g	50g	2kg
 <p>SIT8185.8 1-(TRIETHOXSILYL)-2-(DIETHOXYMETHYLSILYL)ETHANE C₁₅H₃₂O₅Si</p> <p>Dipodal silane</p> <p>Lower toxicity, easier to handle than bis(triethoxysilyl)ethane Improves hydrolytic stability of silane adhesion promotion systems [18418-54-7]</p>	324.56	100° / 0.5		0.946	1.4112
			Flashpoint: 102°C (216°F) TOXICITY: oral rat, LD50: >500 mg/kg		
	TSCA	HMS: 2-1-1-X	25g	100g	2kg

Organosilane-Modified Silica Nanoparticles


A range of silica structures from 20 nm to 1 micron have been modified with silanes to reduce hydroxyl content allowing improved dispersion. Other versions have monolayers with isolated secondary amine functionality, providing controlled interactions with resins. Systems that maintain low levels of hydroxyls have improved electrical properties. Introduction of low levels of secondary amines impart improved mechanical properties particularly in high humidity environments.




<p>SIS6960.0 SILICON DIOXIDE, amorphous <i>Fumed silica</i> SiO₂</p> <p>Ultimate particle size: 12 - 20 nm Surface area, 200 m²/g Isoelectric point: 2.2</p> <p>[112945-52-5]</p>	60.09		(>1,600°)	2.2	1.46
			TOXICITY: oral rat, LD50: 8,160 mg/kg yc: 44 mN/m Bulk density: ~50 g/l pH, (4% aqueous slurry): 3.5-4.5		
	TSCA	HMS: 2-0-0-X	500g	2kg	



<p>SIS6962.0 SILICON DIOXIDE, amorphous, HEXAMETHYLDISILAZANE TREATED <i>Fumed silica, HMDZ treated</i> SiO₂</p> <p>Surface area, 150-200 m²/g</p> <p>[68909-20-6]/[7631-86-9]</p>	60.09		(>1,600°)	2.2	1.45
			Carbon content: 3% Calculated ratio: (CH ₃) ₃ Si/HO-Si: 2/1		
	TSCA	EC 272-697-1	HMS: 2-0-0-X	500g	2kg



<p>SIS6962.1M30 SILICON DIOXIDE, amorphous, HEXAMETHYLDISILAZANE TREATED <i>Fumed silica, HMDZ treated</i> SiO₂</p> <p>Surface area, 150-200 m²/g</p> <p>[68909-20-6]/[7631-86-9]</p>	60.09		(>1,600°)	2.2	1.45
			Carbon content: 2-3% Calculated ratio: (CH ₃) ₃ Si/HO-Si: 1/1		
	TSCA	EC 272-697-1	HMS: 2-0-0-X	500g	2kg



<p>SIS6962.1N30 SILICON DIOXIDE, amorphous, CYCLIC AZASILANE/HEXAMETHYLDISILAZANE TREATED SiO₂</p> <p>Surface area, 150-200 m²/g</p> <p>[68909-20-6]/[7631-86-9]</p>	60.09		(>1,600°)	2.2	1.45
			Carbon content: 4-7% Calculated ratio: CH ₃ NHCH ₂ CH ₂ CH ₂ Si/(CH ₃) ₃ Si/HO-Si: 1/2/1		
	TSCA	HMS: 2-0-0-X	500g		

Gelest provides custom surface treatment services. We can handle a wide range of materials with special process considerations including: inert atmospheres, highly flammable and corrosive treatments, as well as thermal and vacuum drying.

Surface Modification with Silanes: What's not covered in "Silane Coupling Agents"?

Chlorosilane, silazane and dialkylaminosilane coupling agents are not discussed in this brochure. These materials can be found in the Gelest catalog entitled "Silicon Compounds: Silanes & Silicones." The use of these materials is limited commercially due to the difficulty in handling the corrosive, flammable or toxic byproducts associated with hydrolysis.

Hydrophobic, Hydrophilic and Polar silanes, although important in surface modification, do not have reactive organic functionality and are not discussed with coupling agents. Please see the Gelest brochure entitled "Hydrophobicity, Hydrophilicity and Silane Surface Modification" includes these materials.

Further Reading

Silane Coupling Agents - General References and Proceedings

1. B. Arkles, Tailoring Surfaces with Silanes, CHEMTECH, 7, 766-778, 1977
2. E. Plueddemann, "Silane Coupling Agents," Plenum, 1982.
3. K. Mittal, "Silanes and Other Coupling Agents," VSP, 1992.
4. D. Leyden and W. Collins, "Silylated Surfaces," Gordon & Breach, 1980.
5. D. E. Leyden, "Silanes, Surfaces and Interfaces," Gordon & Breach 1985.
6. J. Steinmetz and H. Mottola, "Chemically Modified Surfaces," Elsevier, 1992.
7. J. Blitz and C. Little, "Fundamental & Applied Aspects of Chemically Modified Surfaces," Royal Society of Chemistry, 1999.
8. B. Arkles, Y. Pan, G. Larson, M. Singh, Chemistry - A European Journal, 20, 9442-9450, 2014.

Substrate Chemistry - General References and Proceedings

9. R. Iler, "The Chemistry of Silica," Wiley, 1979.
10. S. Pantelides, G. Lucovsky, "SiO₂ and Its Interfaces," MRS Proc. 105, 1988.

Product Information		Molecular Weight		Refractive Index	
Product Code	SIA0540.0	Boiling point/mm (Melting Point)	162.26	146-8	Specific Gravity
Product Name	ALLYLTRIMETHOXSILANE	Flashpoint: 46°C (115°F)	0.963 ²⁵		1.4036 ²⁵
Empirical Formula	C ₆ H ₁₄ O ₃ Si	Other Physical Properties			
<chem>H2C=CHCH2Si(OCH3)3</chem>	Adhesion promoter for vinyl-addition silicones	References			
	Allylation of ketones, aldehydes and imines w/ dual activation of a Lewis Acid and fluoride ion. ¹				
	1. Yamasaki, S.; et al. <i>J. Am. Chem. Soc.</i> 2002 , <i>124</i> , 6536.				
	F&F: Vol 18, p 14; Vol 19, p 360; Vol 20, p 85; Vol 21, p 3, Vol 12, p 395				
	HYDROLYTIC SENSITIVITY: 7 reacts slowly with moisture/water				
CAS #	[2551-83-9]	10g	\$28.00	50g	\$112.00
European Registration #	TSCA EC 219-855-8	2kg	\$840.00	Hazardous Material Information	
	Indicates Product Listed in TSCA Inventory (E= Exempt - Naturally Occurring Substance) (L= Low Volume Exemption) (S= Significant New Use Restriction)	System Ratings (Health-Flammability-Reactivity) See Below			
	HYDROLYTIC SENSITIVITY: 10 most sensitive to water; 0 least sensitive (see p.13 for details)				

Commercial Status—produced on a regular basis for inventory

Developmental Status—available to support development and commercialization

New Products—available to support development and commercialization

Gelest Product Lines




Silicon Compounds: Silanes & Silicones

Detailed chemical properties and reference articles for over 3,000 compounds. The 608 page handbook of silane and silicone chemistry includes scholarly reviews as well as detailed application information.



Metal-Organics for Material & Polymer Technology

A reference manual for optical and electronic and nanotechnology applications. The literature provides information on metallization, electroceramic, and dielectric applications of silicon, germanium, aluminum, gallium, copper and other metal chemistries. Deposition techniques include ALD, CVD, spin coating and self-assembled monolayers (SAMs). Presents chemistry and physics in the context of device applications ranging from ULSI semiconductors to DNA array devices to flat-panel displays.




Reactive Silicones – Forging New Polymer Links

The 64 page brochure describes reactive silicones that can be formulated into coatings, membranes, cured rubbers and adhesives for mechanical, optical, electronic and ceramic applications. Information on reactions and cures of silicones as well as physical property shortens product development time for chemists and engineers.



Silicone Fluids – Stable, Inert Media

Design and Engineering properties for conventional silicone fluids as well as thermal, fluorosilicone, hydrophilic and low temperature grades are presented in a 30 page selection guide. The brochure provides data on thermal, rheological, electrical, mechanical and optical properties for silicones. Silicone fluids are available in viscosities ranging from 0.65 to 2,500,000 cSt.



Hydrophobicity, Hydrophilicity and Silane Surface Modification

A description of non-functional silanes that are used to prepare hydrophobic and water repellent surfaces, as well as polar and hydroxylic silanes used to prepare wettable surfaces.



Enabling Your Technology

11 East Steel Road
Morrisville, PA 19067
USA
Tel: (215) 547-1015
Toll-Free: (888) 734-8344
Fax: (215) 547-2484
Internet: www.gelest.com

46 Pickering Street
Maidstone
Kent ME15 9RR
United Kingdom
Tel: +44-(0)-1622-741115
Fax: +44-(0)-8701-308421
Email: europe@gelest.com

Stroofstrasse 27 Geb.2901
65933 Frankfurt am Main
Germany
Tel: +49-(0)-69-3800-2150
Fax: +49-(0)-69-3800-2300
Email: info@gelestde.com
Internet: www.gelestde.com